



Technical Catalogue
SCALA laboratory furniture system

SCALA



WALDNER



Technical Catalogue

The design of our **SCALA** range of laboratory furniture will set the trend for future laboratory design.

But only if design and functionality work together effectively, real values will result that can contribute to enrich the laboratory environment.

We have redesigned our range of laboratory furniture based on innovative ideas, sophisticated detailed solutions and high-quality materials, thus meeting the requirements of our users with respect to ergonomics and profitability more than ever.

Our **SCALA** laboratory furniture system with its flexible application units can easily be adapted to new room situations. In this way we can provide a large number of different design and furniture variants for every functional area of the laboratory.

With our **SCALA** laboratory furniture we offer innovative, mature technology, maximum operational safety, ergonomic design and perfect service. Discover all details of our furniture on the following pages.

Not without good reason have customers from all over the world relied on us and our service for more than 70 years.

With this technical catalogue, we are providing you with the basis for your future laboratory.

Contact us. Our specialists will always be pleased to talk to you.





Technical Catalogue Table of contents

1

Fume hoods and extraction devices

| | |
|--|----|
| Fume hoods and extraction devices | 6 |
| Bench-mounted fume hoods with rear wall installation | 12 |
| Bench-mounted fume hood..... | 12 |
| Low ceiling bench-mounted fume hood..... | 15 |
| Secuflow bench-mounted fume hood..... | 18 |
| Secuflow low ceiling bench-mounted fume hood..... | 21 |
| Bench-mounted fume hoods with side installation | 24 |
| Secuflow bench-mounted fume hood with side installation..... | 24 |
| Secuflow low ceiling bench-mounted fume hood with side installation..... | 27 |
| Secuflow bench-mounted fume hood with side installation for work performed while seated..... | 30 |
| Fume hood with side installation, made of steel SI 3 steel..... | 33 |
| Walk-in fume hoods | 36 |
| Walk-in fume hood with side installation..... | 36 |
| Low level fume hoods | 39 |
| Low level fume hood with side installation..... | 39 |

Fume hoods and extraction devices

| | |
|--|----|
| Special fume hoods | 42 |
| Secuflow EN7 fume hood for high thermal loads..... | 42 |
| EN7 fume hood for high thermal loads in connection with acid digestions (special application fume hood)..... | 45 |
| Fume hood for handling of perchloric acid..... | 48 |
| Fume hood for handling of hydrofluoric acid..... | 51 |
| Radio-isotope fume hood..... | 54 |
| Filter fume hood..... | 58 |
| Safety weighing cabinet | 62 |
| AKKURAT..... | 62 |
| Mobile fume hoods | 65 |
| AeroEm..... | 65 |
| MobilAir..... | 68 |
| Housings | 71 |
| Permanent enclosure..... | 71 |
| Local extraction devices | 73 |
| Underbench unit extraction system..... | 73 |
| Extraction arm..... | 74 |
| Snorkel hood..... | 75 |
| Extractor hood..... | 76 |

2

Service modules

| | |
|------------------------------------|-----------|
| Service modules | 78 |
| Service duct element..... | 84 |
| Service spine..... | 86 |
| Service wing..... | 90 |
| Suspended service boom..... | 93 |
| Service column..... | 95 |
| Service distribution terminal..... | 98 |
| Service wall duct..... | 99 |
| Bench-mounted service duct..... | 100 |

3

Laboratory benches and sinks

| | |
|---|------------|
| Laboratory benches and sinks | 104 |
| Laboratory benches | 110 |
| Combinations of worktop and bench frame materials..... | 110 |
| Worktop material..... | 112 |
| Bench with H-frame..... | 117 |
| Bench with C-frame..... | 118 |
| Bench with self-supporting underbench units..... | 119 |
| Mobile tables | 120 |
| Laboratory sinks | 121 |
| Laboratory sink..... | 121 |
| End sink..... | 122 |
| Drip cup on service spine..... | 123 |
| Drip cup in worktop..... | 124 |
| Mobile sink..... | 125 |
| AquaEl..... | 126 |
| Special tables | 127 |
| Add-on table for low level fume hoods..... | 127 |
| Balance table..... | 128 |
| Rack..... | 129 |
| Heavy duty rack..... | 130 |
| Heavy duty mobile table..... | 131 |
| Height-adjustable table..... | 132 |

4

Storage cupboards

| | |
|--|------------|
| Storage cupboards | 134 |
| Underbench units | 138 |
| Underbench unit on plinth..... | 138 |
| Underbench unit on castors..... | 140 |
| Suspended underbench unit..... | 142 |
| Push-in underbench unit for fume hoods..... | 144 |
| Underbench unit for sinks..... | 146 |
| Overbench cabinets | 149 |
| Laboratory cabinets | 151 |
| Laboratory cabinet..... | 151 |
| Emergency cabinet..... | 155 |
| Top-mounted cabinets | 156 |
| Pull-out cabinets | 157 |

4

Storage cupboards

| | |
|---|------------|
| Special cabinets | 159 |
| Laboratory cabinet for storing acids and alkalis..... | 159 |
| Underbench safety unit for storing acids and alkalis..... | 161 |
| FWF 90 underbench safety unit for fume hoods for storing flammable liquids..... | 163 |
| FWF 90 safety cabinet for storing flammable liquids..... | 165 |
| G 90 gas cylinder cabinet..... | 167 |

5

Supply and disposal

| | |
|---|------------|
| Supply and disposal | 170 |
| Supply system for flammable liquids..... | 172 |
| Waste disposal system for acids and alkalis..... | 175 |
| Waste disposal system for flammable liquids..... | 178 |
| Waste disposal system for solid matter and domestic waste..... | 181 |
| Waste disposal system for radio-isotope residual material..... | 183 |

6

Education

| | |
|---|------------|
| Waldner educational institutions | 184 |
|---|------------|

7

Services

| | |
|-----------------------------|------------|
| Services | 186 |
| Control and Monitoring..... | 190 |

8

Dimensions

| | |
|-------------------------|------------|
| Dimensions | 200 |
|-------------------------|------------|

9

Accessories

| | |
|--|------------|
| Accessories for laboratory furniture system | 202 |
|--|------------|

10

General

| | |
|--|------------|
| General | 204 |
| Colours..... | 206 |
| Laboratory planning..... | 208 |
| Awards..... | 210 |
| Installation interfaces mechanical and electrical services..... | 212 |



1 Fume hoods and extraction devices

Energy efficiency, maximum ergonomics and a larger internal workspace make working with our new fume hoods even safer and more convenient.

Design together with an enlarged product range characterise the fume hoods of our **SCALA** laboratory range.

Combined with grid lengths up to 2400 mm of our fume hoods, we offer the most comprehensive product range available in the market. Almost all fume hoods are also available with Securflow technology.



| | | | |
|--|-----------|--|------------|
| Bench-mounted fume hoods with rear wall installation..... | 12 | Special fume hoods..... | 42 |
| Bench-mounted fume hood..... | 12 | Secuflow EN7 fume hood for high thermal loads..... | 42 |
| Low ceiling bench-mounted fume hood..... | 15 | EN7 fume hood for high thermal loads in connection with acid digestions (special application fume hood)..... | 45 |
| Secuflow bench-mounted fume hood..... | 18 | Fume hood for handling of perchloric acid..... | 48 |
| Secuflow low ceiling bench-mounted fume hood | 21 | Fume hood for handling of hydrofluoric acid..... | 51 |
| Bench-mounted fume hoods with side installation..... | 24 | Radio-isotope fume hood..... | 54 |
| Secuflow bench-mounted fume hood with side installation..... | 24 | Filter fume hood..... | 58 |
| Secuflow low ceiling bench-mounted fume hood with side installation..... | 27 | Safety weighing cabinet..... | 62 |
| Secuflow bench-mounted fume hood with side installation for work performed while seated..... | 30 | AKKURAT..... | 62 |
| Fume hood with side installation, made of steel SI 3 steel..... | 33 | Mobile fume hoods..... | 65 |
| Walk-in fume hoods..... | 36 | AeroEm..... | 65 |
| Walk-in fume hood with side installation | 36 | MobilAir..... | 68 |
| Low level fume hoods..... | 39 | Housings..... | 71 |
| Low level fume hood with side installation..... | 39 | Permanent enclosure..... | 71 |
| | | Local extraction devices..... | 73 |
| | | Underbench unit extraction system..... | 73 |
| | | Extraction arm..... | 74 |
| | | Snorkel hood..... | 75 |
| | | Extractor hood..... | 76 |
| | | Control and monitoring..... | 190 |



1 Fume hoods and extraction devices

All laboratory work during which gases, fumes, particles or liquids are handled in dangerous quantities and concentrations must be performed in fume hoods.

All our fume hoods ensure maximum safety, excellent ergonomics and maximum economy.

Reduced energy consumption – increased profitability

The fluid mechanics have been further optimised which means considerably reduced energy consumption of our fume hoods while maintaining the high safety level. Our bench-mounted fume hoods with Secuflow technology which are tested in accordance with EN 14175, e.g., require 270 m³/h/lfm.

As an important part of the overall laboratory ventilation scheme, our fume hoods can be perfectly integrated into the building ventilation concept.

The fact that our Secuflow fume hood technology also reduces the investment and operating costs for the ventilation system is another commercial advantage that is made possible by the integrated supportive flow technology. You will find further information on this topic in our Secuflow brochure.

Improved ergonomics with the inclined operating panel

The operating panel is inclined towards the user for easier handling and operation of all fittings and functions.



Fume hoods and extraction devices

1

Safety through the intake airflow profile on the front edge of the fume hood worktop

It prevents turbulence that could carry pollutant emissions.

Air flowing into the fume hood is guided via the airfoil-like profile geometry (with low turbulence) over the worktop to the rear panel low level extraction which ensures the safe removal of heavy gases, e.g. solvent fumes, directly above the worktop.

For more safety

Maximum user safety is provided by our toothed belt sash mounting along with significantly reduced maintenance effort. The stainless steel reinforced toothed belts prove maximum resistance during endurance tests with more than 200,000 load cycles. The shape of the sash frame offers maximum protection from splashes and splinters.

Anti-slip device for additional protection

In the unlikely case that both sash mountings fail, the sash is stopped in fractions of a second.

Largest possible access area

The slender side posts of our fume hoods offer an increased nominal width of the internal workspace and due to their special shape ensure that there is little turbulence in the intake air.

Larger capacity of the internal workspac

The internal workspace is 10 % higher thus increasing the entire internal workspace. Useful when working with tall and wide items of experimental equipment.



1 Fume hoods and extraction devices

Clear view of all processes in the workspace

The high level glazed panel enables tall experimental equipment and processes to be clearly seen.

Scaffold points

Scaffold rods with diameters of 12 and 13 mm can be firmly secured.

All functions at a glance

The Soft Touch control element integrated in the fume hood side post provides information on the operational state of the fume hood at eye level.

Sash handle with air guiding function

Air is pushed into the workspace when the sash is opened and pollutant emissions due to the opening sash are prevented. The balanced and free-moving sash mechanism including the release for the sash stop can be operated with one hand.

The automatic sash

The sash is closed automatically if there is nobody working on the fume hood. The photo-electric barrier stops the closing process if there are objects protruding from inside the workspace.

Various fume hood widths available

Our bench-mounted fume hoods are available in widths up to 2100 mm and the side-installed fume hoods up to 1800 mm - of course also in Secuflow technology.

Lighting for the internal workspace

LED lamps uniformly illuminate the fume hood interior - conveniently switchable from the side post.



The best for equipment and variability

Along with the convenient basic equipment, our fume hoods provide a wide range of variable equipment options. Depending on the application, the worktop is made of stoneware, epoxy resin, polypropylene or stainless steel. Our fume hoods are mounted on a steel support frame. You can install plinth mounted, mobile or solvent cabinets under the fume hood.

Service modules that can be equipped as desired

The replaceable service modules are integrated in the rear and side panels of our fume hoods and ensure the mechanical and electrical services supply. The integrated sink module for water offers more freedom when using the internal workspace.

Our certified test laboratory for fume hood measurements

We established our new test laboratory for fume hoods when the EN 14175 was published. The latest technical equipment and the certification by TÜV NORD CERT GmbH guarantee optimum measurement results with respect to accuracy and reproducibility.

We test fume hoods in accordance with EN 14175. We can also carry out measurements in accordance with ASHRAE 110.

With our ISO 9001 certification and the GS mark for our entire product range, we have closed the circle in relation to fume hood tests and had our test laboratory tested and certified by TÜV NORD CERT GmbH according to the German law on equipment safety (Gerätesicherheitsgesetz).


Bench-mounted fume hoods

Bench-mounted fume hood

Intended use

- Protective device for the user, tested in acc. with EN 14175
- Extraction of fumes, aerosols and dust from the internal workspace to prevent dangerous amounts of pollutants from escaping into the laboratory
- Reduced risk of the formation of a high concentration of hazardous substances / hazardous explosive atmosphere in the internal workspace
- Protection from splashes of hazardous substances
- Protection from flying particles, bodies or parts escaping from the internal workspace
- General fume hoods constructed in acc. with EN 14175 are normally not suited for use with radioactive substances or microorganisms
- Not suitable for openly breaking down chemicals
- Service outlets in the rear panel of the internal workspace
- Control units located horizontally on the service rail of the support unit

Design



- 1 Sash with handle and horizontal sashes
- 2 Worktop
- 3 FAZ or AC control panel
- 4 Upper sash window
- 5 Removable fascia panel
- 6 Exhaust hood
- 7 Baffle with service modules
- 8 Glass pane in the side wall
- 9 Material lock
- 10 Bench frame with push-in underbench units with support and service panels


Bench-mounted fume hoods

Bench-mounted fume hood

1

Fume hoods and extraction devices

Dimensional drawing



Technical data

| Dimensions | 1200 | 1500 | 1800 | 2100 |
|---------------------------------------|------|------|------|------|
| Width [mm] | 1200 | 1500 | 1800 | 2100 |
| Depth [mm] | 900 | | | |
| Height [mm] | 2700 | | | |
| Clear width, internal workspace [mm] | 1150 | 1450 | 1750 | 2050 |
| Clear height, internal workspace [mm] | 1550 | | | |
| Working height [mm] | 900 | | | |

| Weight | 1200 | 1500 | 1800 | 2100 |
|---------------------------|-------------|-------------|-------------|-------------|
| Without installation [kg] | Approx. 250 | Approx. 300 | Approx. 350 | Approx. 400 |

Bench-mounted fume hoods

Bench-mounted fume hood

| Design characteristics | 1200 | 1500 | 1800 | 2100 |
|---|---|------|---------------------|------|
| Supporting construction | H-frame with push-in underbench units | | | |
| Sash | 2 horizontal sashes | | 3 horizontal sashes | |
| Side panel of the fume hood | Glass pane on the left and/or right as an option; not with stoneware internal lining Material lock on the left and/or right as an option | | | |
| Number of devices for scaffold points, ø 12 to 13 mm | 9 | | 12 | |
| Max. load per scaffold point with scaffold rod length 300 mm [kg] | 5 | | | |
| Service modules | 2 | | 3 | |

| Electrics | |
|--------------------|---|
| Electrical supply | External sockets in service panels Internal sockets in service modules |
| Fuse box | Optional |
| Sash controller SC | Optional |

| Sanitary technology | |
|---------------------|--|
| Sanitary supply | Service modules with take-off valves for vacuum, gases and/or waters and integrated sink (PP) as an option |

| Ventilation technology | 1200 | 1500 | 1800 | 2100 |
|---|---|------|------|------|
| Minimum air exchange rate [m³/h] ¹⁾ | 480 | 600 | 720 | 840 |
| Function display | FAZ | | | |
| Airflow damper, constant | Airflow-Controller AC | | | |
| Airflow damper, variable | Airflow-Controller AC | | | |
| Detector of sash position | Only variable with Airflow-Controller AC | | | |
| Connection height [mm] for FAZ with extract air hood ø 250 mm | 2570 | | | |
| Connection height [mm] for FAZ with extract air hood ø 315 mm ²⁾ | 2570 | | | |
| Connection height [mm] for AC with extract air hood ø 250 mm | 2940 | | | |
| Connection height [mm] for AC with extract air hood ø 315 mm ²⁾ | 2910 | | | |
| Underbench unit extraction system | As an option, depending on requirements and regulations | | | |

¹⁾ All air volume specifications refer to an opening height of the sash window of 500 mm (test opening in acc. with EN 14175-3) and the maximum tracer gas values recommended by German Standard (BG RCI). Shown rates correspond to a face velocity of 0.24 m/s. For other design face velocities, please contact your Waldner sales representative.

²⁾ In order to minimise noise and pressure losses, for air volumes >1000 m³/h Waldner recommends using the extract air hood with a connection diameter of 315 mm.

A maximum admission pressure of 600 Pa should not be exceeded in the case of fume hoods with airflow dampers.

The indicated minimum air exchange rates were determined under specified test conditions in acc. with EN 14175-3. These minimum air exchange rates must be adapted when dimensioning the ventilation system

If on-site extract air monitoring systems or airflow dampers are used, the required air volumes may be different. The operating limitations must be agreed upon with Waldner.

| Material/surface | |
|------------------|--|
| Worktop | Stoneware Polypropylene Epoxy Stainless steel |
| Internal lining | Melamine resin facing Solid grade laminate Stoneware |


Bench-mounted fume hoods

Low ceiling bench-mounted fume hood

Intended use

- Protective device for the user, tested in acc. with EN 14175
- Extraction of fumes, aerosols and dust from the internal workspace to prevent dangerous amounts of pollutants from escaping into the laboratory
- Reduced risk of the formation of a high concentration of hazardous substances / hazardous explosive atmosphere in the internal workspace
- Protection from splashes of hazardous substances
- Protection from flying particles, bodies or parts escaping from the internal workspace
- General fume hoods constructed in acc. with EN 14175 are normally not suited for use with radioactive substances or microorganisms
- Not suitable for openly breaking down chemicals
- Service outlets in the rear panel of the internal workspace
- Control units located horizontally on the service rail of the support unit
- Suitable for rooms with low ceiling height

Design




- 1 Two-piece sash with handle and horizontal sashes
- 2 Worktop
- 3 FAZ or AC control panel
- 4 Removable fascia panel
- 5 Exhaust hood
- 6 Baffle with service modules
- 7 Glass pane in the side wall
- 8 Material lock
- 9 Bench frame with push-in underbench units with support and service panels

Bench-mounted fume hoods

Low ceiling bench-mounted fume hood

Dimensional drawing



Technical data

| Dimensions | 1200 | 1500 | 1800 | 2100 |
|---------------------------------------|------|------|------|------|
| Width [mm] | 1200 | 1500 | 1800 | 2100 |
| Depth [mm] | 900 | | | |
| Height [mm] | 2400 | | | |
| Clear width, internal workspace [mm] | 1150 | 1450 | 1750 | 2050 |
| Clear height, internal workspace [mm] | 1250 | | | |
| Working height [mm] | 900 | | | |

| Weight | 1200 | 1500 | 1800 | 2100 |
|---------------------------|-------------|-------------|-------------|-------------|
| Without installation [kg] | Approx. 220 | Approx. 260 | Approx. 300 | Approx. 350 |

| Design characteristics | 1200 | 1500 | 1800 | 2100 |
|---|---|------|---------------------|------|
| Supporting construction | H-frame with push-in underbench units | | | |
| Two-piece sash | 2 horizontal sashes | | 3 horizontal sashes | |
| Side panel of the fume hood | Glass pane on the left and/or right as an option; not with stoneware internal lining Material lock on the left and/or right as an option | | | |
| Max. number of devices for scaffold points, ø 12 to 13 mm | 9 | | 12 | |
| Max. load per scaffold point with scaffold rod length 300 mm [kg] | 5 | | | |
| Service modules | 2 | | 3 | |

Bench-mounted fume hoods

Low ceiling bench-mounted fume hood

| Electrics | |
|--------------------|---|
| Electrical supply | External sockets in service panels Internal sockets in service modules |
| Fuse box | Optional |
| Sash controller SC | Optional |

| Sanitary technology | |
|---------------------|--|
| Sanitary supply | Service modules with take-off valves for vacuum, gases and/or waters and integrated sink (PP) as an option |

| Ventilation technology | 1200 | 1500 | 1800 | 2100 |
|---|---|------|------|------|
| Minimum air exchange rate [m³/h] ¹⁾ | 480 | 600 | 720 | 840 |
| Function display | FAZ | | | |
| Airflow damper, constant | Airflow-Controller AC | | | |
| Airflow damper, variable | Airflow-Controller AC | | | |
| Detector of sash position | Only variable with Airflow-Controller AC | | | |
| Connection height [mm] for FAZ with extract air hood Ø 250 mm | 2270 | | | |
| Connection height [mm] for FAZ with extract air hood Ø 315 mm ²⁾ | 2270 | | | |
| Connection height [mm] for AC with extract air hood Ø 250 mm | 2640 | | | |
| Connection height [mm] for AC with extract air hood Ø 315 mm ²⁾ | 2610 | | | |
| Underbench unit extraction system | As an option, depending on requirements and regulations | | | |

¹⁾ All air volume specifications refer to an opening height of the sash window of 500 mm (test opening in acc. with EN 14175-3) and the maximum tracer gas values recommended by German Standard (BG RCI). Shown rates correspond to a face velocity of 0.24 m/s. For other design face velocities, please contact your Waldner sales representative.

²⁾ In order to minimise noise and pressure losses, for air volumes >1000 m³/h Waldner recommends using the extract air hood with a connection diameter of 315 mm.

A maximum admission pressure of 600 Pa should not be exceeded in the case of fume hoods with airflow dampers.

The indicated minimum air exchange rates were determined under specified test conditions in acc. with EN 14175-3. These minimum air exchange rates must be adapted when dimensioning the ventilation system.

If on-site extract air monitoring systems or airflow dampers are used, the required air volumes may be different. The operating limitations must be agreed upon with Waldner.

| Material/surface | |
|------------------|--|
| Worktop | Stoneware Polypropylene Stainless steel Epoxy |
| Internal lining | Melamine resin facing Solid grade laminate Stoneware |


Bench-mounted fume hoods

Secuflow bench-mounted fume hood

Intended use

- Protective device for the user, tested in acc. with EN 14175
- Extraction of fumes, aerosols and dust from the internal workspace to prevent dangerous amounts of pollutants from escaping into the laboratory
- Reduced risk of the formation of a high concentration of hazardous substances / hazardous explosive atmosphere in the internal workspace
- Protection from splashes of hazardous substances
- Protection from flying particles, bodies or parts escaping from the internal workspace
- General fume hoods constructed in acc. with EN 14175 are normally not suited for use with radioactive substances or microorganisms
- Not suitable for openly breaking down chemicals
- Active supportive flow technology (Secuflow technology) reduces the energy consumption while regulations and standards are observed
- Service outlets in the rear panel of the internal workspace
- Control units located horizontally on the service rail of the support unit

Design




- 1 Sash with handle and horizontal sashes
- 2 Worktop
- 3 FAZ or AC control panel
- 4 Upper sash window
- 5 Removable fascia panel
- 6 Exhaust hood
- 7 Baffle with service modules
- 8 Glass pane in the side wall
- 9 Material lock
- 10 Bench frame with push-in underbench units with support and service panels

Bench-mounted fume hoods

Secuflow bench-mounted fume hood

Dimensional drawing



Technical data

| Dimensions | 1200 | 1500 | 1800 | 2100 |
|---------------------------------------|------|------|------|------|
| Width [mm] | 1200 | 1500 | 1800 | 2100 |
| Depth [mm] | 900 | | | |
| Height [mm] | 2700 | | | |
| Clear width, internal workspace [mm] | 1150 | 1450 | 1750 | 2050 |
| Clear height, internal workspace [mm] | 1550 | | | |
| Working height [mm] | 900 | | | |

| Weight | 1200 | 1500 | 1800 | 2100 |
|---------------------------|-------------|-------------|-------------|-------------|
| Without installation [kg] | Approx. 250 | Approx. 300 | Approx. 350 | Approx. 400 |

Bench-mounted fume hoods

Secuflow bench-mounted fume hood

| Design characteristics | 1200 | 1500 | 1800 | 2100 |
|---|---|------|---------------------|------|
| Supporting construction | H-frame with push-in underbench units | | | |
| Sash | 2 horizontal sashes | | 3 horizontal sashes | |
| Side panel of the fume hood | Glass pane on the left and/or right as an option; not with stoneware internal lining Material lock on the left and/or right as an option | | | |
| Max. number of devices for scaffold points, ø 12 mm to 13 mm | 9 | | 12 | |
| Max. load per scaffold point with scaffold rod length 300 mm [kg] | 5 | | | |
| Service modules | 2 | | 3 | |

| Electrics | |
|--------------------|---|
| Electrical supply | External sockets in service panels Internal sockets in service modules |
| Fuse box | Optional |
| Sash controller SC | Optional |

| Sanitary technology | |
|---------------------|--|
| Sanitary supply | Service modules with take-off valves for vacuum, gases and/or waters and integrated sink (PP) as an option |

| Ventilation technology | 1200 | 1500 | 1800 | 2100 |
|---|---|------|------|------|
| Minimum air exchange rate [m ³ /h] ¹⁾ | 330 | 410 | 490 | 570 |
| Function display | FAZ | | | |
| Airflow damper, constant | Airflow-Controller AC | | | |
| Airflow damper, variable | Airflow-Controller AC | | | |
| Detector of sash position | Only variable with Airflow-Controller AC | | | |
| Connection height [mm] for FAZ with extract air hood ø 250 mm | 2570 | | | |
| Connection height [mm] for FAZ with extract air hood ø 315 mm ²⁾ | 2570 | | | |
| Connection height [mm] for AC with extract air hood ø 250 mm | 2940 | | | |
| Connection height [mm] for AC with extract air hood ø 315 mm ²⁾ | 2910 | | | |
| Underbench unit extraction system | As an option, depending on requirements and regulations | | | |

¹⁾ All air volume specifications refer to an opening height of the sash window of 500 mm (test opening in acc. with EN 14175-3) and the maximum tracer gas values recommended by German Standard (BG RCI). Shown rates correspond to a face velocity of 0.17 m/s (+/- 0.01 m/s). For other design face velocities, please contact your Waldner sales representative.

²⁾ In order to minimise noise and pressure losses, for air volumes >1000 m³/h Waldner recommends using the extract air hood with a connection diameter of 315 mm. A maximum admission pressure of 600 Pa should not be exceeded in the case of fume hoods with airflow dampers. The indicated minimum air exchange rates were determined under specified test conditions in acc. with EN 14175-3. These minimum air exchange rates must be adapted when dimensioning the ventilation system. If on-site extract air monitoring systems or airflow dampers are used, the required air volumes may be different. The operating limitations must be agreed upon with Waldner.

| Material/surface | |
|------------------|--|
| Worktop | Stoneware Polypropylene Stainless steel Epoxy |
| Internal lining | Melamine resin facing Solid grade laminate Stoneware |


Bench-mounted fume hoods

Secuflow low ceiling bench-mounted fume hood

Intended use

- Protective device for the user, tested in acc. with EN 14175
- Extraction of fumes, aerosols and dust from the internal workspace to prevent dangerous amounts of pollutants from escaping into the laboratory
- Reduced risk of the formation of a high concentration of hazardous substances / hazardous explosive atmosphere in the internal workspace
- Protection from splashes of hazardous substances
- Protection from flying particles, bodies or parts escaping from the internal workspace
- General fume hoods constructed in acc. with EN 14175 are normally not suited for use with radioactive substances or microorganisms
- Not suitable for openly breaking down chemicals
- Active supportive flow technology (Secuflow technology) reduces the energy consumption while regulations and standards are observed
- Service outlets in the rear panel of the internal workspace
- Control units located horizontally on the service rail of the support unit
- Suitable for rooms with low ceiling height

Design




- 1 Two-piece sash with handle and horizontal sashes
- 2 Worktop
- 3 FAZ or AC control panel
- 4 Removable fascia panel
- 5 Exhaust hood
- 6 Baffle with service panel
- 7 Glass pane in the side wall
- 8 Material lock
- 9 Bench frame with push-in underbench units with support and service panels

Bench-mounted fume hoods

Secuflow low ceiling bench-mounted fume hood

Dimensional drawing



Technical data

| Dimensions | 1200 | 1500 | 1800 | 2100 |
|---------------------------------------|------|------|------|------|
| Width [mm] | 1200 | 1500 | 1800 | 2100 |
| Depth [mm] | 900 | | | |
| Height [mm] | 2400 | | | |
| Clear width, internal workspace [mm] | 1150 | 1450 | 1750 | 2050 |
| Clear height, internal workspace [mm] | 1250 | | | |
| Working height [mm] | 900 | | | |

| Weight | 1200 | 1500 | 1800 | 2100 |
|---------------------------|-------------|-------------|-------------|-------------|
| Without installation [kg] | Approx. 220 | Approx. 260 | Approx. 300 | Approx. 350 |

| Design characteristics | 1200 | 1500 | 1800 | 2100 |
|---|---|------|---------------------|------|
| Supporting construction | H-frame with push-in underbench units | | | |
| Two-piece sash | 2 horizontal sashes | | 3 horizontal sashes | |
| Side panel of the fume hood | Glass pane on the left and/or right as an option; not with stoneware internal lining Material lock on the left and/or right as an option | | | |
| Max. number of devices for scaffold points, ø 12 to 13 mm | 9 | | 12 | |
| Max. load per scaffold point with scaffold rod length 300 mm [kg] | 5 | | | |
| Service modules | 2 | | 3 | |

Bench-mounted fume hoods

Secuflow low ceiling bench-mounted fume hood

| Electrics | |
|--------------------|---|
| Electrical supply | External sockets in service panels Internal sockets in service modules |
| Fuse box | Optional |
| Sash controller SC | Optional |

| Sanitary technology | |
|---------------------|--|
| Sanitary supply | Service modules with take-off valves for vacuum, gases and/or waters and integrated sink (PP) as an option |

| Ventilation technology | 1200 | 1500 | 1800 | 2100 |
|---|---|------|------|------|
| Minimum air exchange rate [m ³ /h] ¹⁾ | 330 | 410 | 490 | 570 |
| Function display | FAZ | | | |
| Airflow damper, constant | Airflow-Controller AC | | | |
| Airflow damper, variable | Airflow-Controller AC | | | |
| Detector of sash position | Only variable with Airflow-Controller AC | | | |
| Connection height [mm] for FAZ with extract air hood Ø 250 mm | 2270 | | | |
| Connection height [mm] for FAZ with extract air hood Ø 315 mm ²⁾ | 2270 | | | |
| Connection height [mm] for AC with extract air hood Ø 250 mm | 2640 | | | |
| Connection height [mm] for AC with extract air hood Ø 315 mm ²⁾ | 2610 | | | |
| Underbench unit extraction system | As an option, depending on requirements and regulations | | | |

¹⁾ All air volume specifications refer to an opening height of the sash window of 500 mm (test opening in acc. with EN 14175-3) and the maximum tracer gas values recommended by German Standard (BG RCI). Shown rates correspond to a face velocity of 0.17 m/s (+/- 0.01 m/s). For other design face velocities, please contact your Waldner sales representative.

²⁾ In order to minimise noise and pressure losses, for air volumes >1000 m³/h Waldner recommends using the extract air hood with a connection diameter of 315 mm.

A maximum admission pressure of 600 Pa should not be exceeded in the case of fume hoods with airflow dampers.

The indicated minimum air exchange rates were determined under specified test conditions in acc. with EN 14175-3. These minimum air exchange rates must be adapted when dimensioning the ventilation system.

If on-site extract air monitoring systems or airflow dampers are used, the required air volumes may be different. The operating limitations must be agreed upon with Waldner.

| Material/surface | |
|------------------|--|
| Worktop | Stoneware Polypropylene Epoxy Stainless steel |
| Internal lining | Melamine resin facing Solid grade laminate Stoneware |


Bench-mounted fume hood with side installation

Secuflow bench-mounted fume hood with side installation

Intended use

- Protective device for the user, tested in acc. with EN 14175
- Extraction of fumes, aerosols and dust from the internal workspace to prevent dangerous amounts of pollutants from escaping into the laboratory
- Reduced risk of the formation of a high concentration of hazardous substances / hazardous explosive atmosphere in the internal workspace
- Protection from splashes of hazardous substances
- Protection from flying particles, bodies or parts escaping from the internal workspace
- General fume hoods constructed in acc. with EN 14175 are normally not suited for use with radioactive substances or microorganisms
- Not suitable for openly breaking down chemicals
- Active supportive flow technology (Secuflow technology) reduces the energy consumption while regulations and standards are observed
- Service outlets in the service modules of the side panels of the internal workspace
- Control units located vertically on the side service panels

Design




- 1 Sash with handle and horizontal sashes
- 2 Worktop
- 3 Service panel
- 4 FAZ or AC control panel
- 5 Service modules in the side panel of the fume hood
- 6 Upper sash window
- 7 Removable fascia panel
- 8 Exhaust hood
- 9 Baffle with scaffold points
- 10 Bench frame with push-in underbench units as an option

Bench-mounted fume hood with side installation

Secuflow bench-mounted fume hood with side installation

Dimensional drawing



Technical data

| Dimensions | 1200 | 1500 | 1800 | 2100 | 2400 |
|---------------------------------------|------|------|------|------|------|
| Width [mm] | 1200 | 1500 | 1800 | 2100 | 2400 |
| Depth [mm] | 900 | | | | |
| Height [mm] | 2700 | | | | |
| Clear width, internal workspace [mm] | 950 | 1250 | 1550 | 1850 | 2150 |
| Clear height, internal workspace [mm] | 1550 | | | | |
| Working height [mm] | 900 | | | | |

| Weight | 1200 | 1500 | 1800 | 2100 | 2400 |
|---------------------------|-------------|-------------|-------------|-------------|-------------|
| Without installation [kg] | Approx. 320 | Approx. 390 | Approx. 450 | Approx. 510 | Approx. 570 |

Bench-mounted fume hood with side installation

Secuflow bench-mounted fume hood with side installation

| Design characteristics | 1200 | 1500 | 1800 | 2100 | 2400 |
|---|---|------|---------------------|------|------|
| Supporting construction | H-frame with push-in underbench units | | | | |
| Sash | 2 horizontal sashes | | 3 horizontal sashes | | |
| Side panel of the fume hood | Glass pane on the left and/or right as an option; not if service modules are installed in the side panel of the fume hood Material lock on the left and/or right as an option | | | | |
| Max. number of devices for scaffold points, ø 12 to 13 mm | 9 | 12 | | 15 | |
| Max. load per scaffold point with scaffold rod length 300 mm [kg] | 5 | | | | |
| Service modules | Service modules in the left and/or right side panel of the fume hood, depending on requirement | | | | |

| Electrics | |
|--------------------|---|
| Electrical supply | External sockets in service panels Internal sockets in service modules |
| Fuse box | Optional |
| Sash controller SC | Optional |

| Sanitary technology | |
|---------------------|--|
| Sanitary supply | Service modules with take-off valves for vacuum, gases and/or waters and integrated sink (PP) as an option |

| Ventilation technology | 1200 | 1500 | 1800 | 2100 | 2400 |
|---|---|------|------|------|------|
| Minimum air exchange rate [m ³ /h] ¹⁾ | 330 | 410 | 490 | 570 | 650 |
| Function display | FAZ | | | | |
| Airflow damper, constant | Airflow-Controller AC | | | | |
| Airflow damper, variable | Airflow-Controller AC | | | | |
| Detector of sash position | Only variable with Airflow-Controller AC | | | | |
| Connection height [mm] for FAZ with extract air hood ø 250 mm | 2570 | | | | |
| Connection height [mm] for FAZ with extract air hood ø 315 mm ²⁾ | 2570 | | | | |
| Connection height [mm] for AC with extract air hood ø 250 mm | 2940 | | | | |
| Connection height [mm] for AC with extract air hood ø 315 mm ²⁾ | 2910 | | | | |
| Underbench unit extraction system | As an option, depending on requirements and regulations | | | | |

¹⁾ All air volume specifications refer to an opening height of the sash window of 500 mm (test opening in acc. with EN 14175-3) and the maximum tracer gas values recommended by German Standard (BG RCI). Shown rates correspond to a face velocity of 0.17 m/s (+/- 0.01 m/s). For other design face velocities, please contact your Waldner sales representative.

²⁾ In order to minimise noise and pressure losses, for air volumes >1000 m³/h Waldner recommends using the extract air hood with a connection diameter of 315 mm. A maximum admission pressure of 600 Pa should not be exceeded in the case of fume hoods with airflow dampers. The indicated minimum air exchange rates were determined under specified test conditions in acc. with EN 14175-3. These minimum air exchange rates must be adapted when dimensioning the ventilation system. If on-site extract air monitoring systems or airflow dampers are used, the required air volumes may be different. The operating limitations must be agreed upon with Waldner.

| Material/surface | |
|------------------|--|
| Worktop | Stoneware Polypropylene Epoxy Stainless steel |
| Internal lining | Solid grade laminate Stainless steel Melamine resin facing |


Bench-mounted fume hood with side installation

Secuflow low ceiling bench-mounted fume hood with side installation

Intended use

- Protective device for the user, tested in acc. with EN 14175
- Extraction of fumes, aerosols and dust from the internal workspace to prevent dangerous amounts of pollutants from escaping into the laboratory
- Reduced risk of the formation of a high concentration of hazardous substances / hazardous explosive atmosphere in the internal workspace
- Protection from splashes of hazardous substances
- Protection from flying particles, bodies or parts escaping from the internal workspace
- General fume hoods constructed in acc. with EN 14175 are normally not suited for use with radioactive substances or microorganisms
- Not suitable for openly breaking down chemicals
- Active supportive flow technology (Secuflow technology) reduces the energy consumption while regulations and standards are observed
- Service outlets in the service modules of the side panels of the internal workspace
- Control units located vertically on the side service panels
- Suitable for rooms with low ceiling height


Design



- 1 Two-piece sash with handle and horizontal sashes
- 2 Worktop
- 3 Service panel
- 4 Service module in the side panel of the fume hood
- 5 FAZ or AC control panel
- 6 Removable fascia panel
- 7 Exhaust hood
- 8 Baffle with scaffold points
- 9 Bench frame with push-in underbench units as an option

Bench-mounted fume hood with side installation Secuflow low ceiling bench-mounted fume hood with side installation

Dimensional drawing



Technical data

| Dimensions | 1200 | 1500 | 1800 |
|---------------------------------------|------|------|------|
| Width [mm] | 1200 | 1500 | 1800 |
| Depth [mm] | 900 | | |
| Height [mm] | 2400 | | |
| Clear width, internal workspace [mm] | 950 | 1250 | 1550 |
| Clear height, internal workspace [mm] | 1250 | | |
| Working height [mm] | 900 | | |

| Weight | 1200 | 1500 | 1800 |
|---------------------------|-------------|-------------|-------------|
| Without installation [kg] | Approx. 220 | Approx. 260 | Approx. 300 |

| Design characteristics | 1200 | 1500 | 1800 |
|---|--|------|---------------------|
| Supporting construction | H-frame with push-in underbench units | | |
| Two-piece sash | 2 horizontal sashes | | 3 horizontal sashes |
| Side panel of the fume hood | Glass pane on the left and/or right as an option; not if service modules are installed in the side panel of the fume hood, not with stoneware internal lining Material lock on the left and/or right as an option | | |
| Max. number of devices for scaffold points, \varnothing 12 to 13 mm | 6 | 9 | |
| Max. load per scaffold point with scaffold rod length 300 mm [kg] | 5 | | |
| Service modules | Service modules in the left and/or right side panel of the fume hood, depending on requirement | | |

Bench-mounted fume hood with side installation

Secuflow low ceiling bench-mounted fume hood with side installation

| Electrics | |
|--------------------|---|
| Electrical supply | External sockets in service panels Internal sockets in service modules |
| Fuse box | Optional |
| Sash controller SC | Optional |

| Sanitary technology | |
|---------------------|--|
| Sanitary supply | Service modules with take-off valves for vacuum, gases and/or waters and integrated sink (PP) as an option |

| Ventilation technology | 1200 | 1500 | 1800 |
|---|---|------|------|
| Minimum air exchange rate [m ³ /h] ¹⁾ | 330 | 410 | 490 |
| Function display | FAZ | | |
| Airflow damper, constant | Airflow-Controller AC | | |
| Airflow damper, variable | Airflow-Controller AC | | |
| Detector of sash position | Only variable with Airflow-Controller AC | | |
| Connection height [mm] for FAZ with extract air hood Ø 250 mm | 2270 | | |
| Connection height [mm] for FAZ with extract air hood Ø 315 mm ²⁾ | 2270 | | |
| Connection height [mm] for AC with extract air hood Ø 250 mm | 2640 | | |
| Connection height [mm] for AC with extract air hood Ø 315 mm ²⁾ | 2610 | | |
| Underbench unit extraction system | As an option, depending on requirements and regulations | | |

¹⁾ All air volume specifications refer to an opening height of the sash window of 500 mm (test opening in acc. with EN 14175-3) and the maximum tracer gas values recommended by German Standard (BG RCI). Shown rates correspond to a face velocity of 0.17 m/s (+/- 0.01 m/s). For other design face velocities, please contact your Waldner sales representative.

²⁾ In order to minimise noise and pressure losses, for air volumes >1000 m³/h Waldner recommends using the extract air hood with a connection diameter of 315 mm.

A maximum admission pressure of 600 Pa should not be exceeded in the case of fume hoods with airflow dampers.

The indicated minimum air exchange rates were determined under specified test conditions in acc. with EN 14175-3. These minimum air exchange rates must be adapted when dimensioning the ventilation system.

If on-site extract air monitoring systems or airflow dampers are used, the required air volumes may be different. The operating limitations must be agreed upon with Waldner.

| Material | |
|-----------------|--|
| Worktop | Stoneware Polypropylene Epoxy Stainless steel |
| Internal lining | Solid grade laminate Stainless steel Melamine resin facing |


Bench-mounted fume hood with side installation

Secuflow bench-mounted fume hood with side installation for work performed while seated

Intended use

- Protective device for the user, tested in acc. with EN 14175
- Extraction of fumes, aerosols and dust from the internal workspace to prevent dangerous amounts of pollutants from escaping into the laboratory
- Reduced risk of the formation of a high concentration of hazardous substances / hazardous explosive atmosphere in the internal workspace
- Protection from splashes of hazardous substances
- Protection from flying particles, bodies or parts escaping from the internal workspace
- General fume hoods constructed in acc. with EN 14175 are normally not suited for use with radioactive substances or microorganisms
- Not suitable for openly breaking down chemicals
- Suitable for work performed while seated
- Active supportive flow technology (Secuflow technology) reduces the energy consumption while regulations and standards are observed
- Service outlets in the service modules of the side panels of the internal workspace
- Control units located vertically on the side service panels

Design




- 1 Sash with handle and horizontal sashes
- 2 Worktop
- 3 Service panel
- 4 FAZ or AC control panel
- 5 Service module in the side panel of the fume hood
- 6 Upper sash window
- 7 Removable fascia panel
- 8 Exhaust hood
- 9 Baffle with scaffold points
- 10 Bench frame with push-in underbench units as an option

Bench-mounted fume hood with side installation

Secuflow bench-mounted fume hood with side installation for work performed while seated

Dimensional drawing



Technical data

| Dimensions | |
|---------------------------------------|------|
| Width [mm] | 1500 |
| Depth [mm] | 900 |
| Height [mm] | 2550 |
| Clear width, internal workspace [mm] | 1250 |
| Clear height, internal workspace [mm] | 1550 |
| Working height [mm] | 750 |

| Weight | |
|---------------------------|-------------|
| Without installation [kg] | Approx. 390 |

Bench-mounted fume hood with side installation

Secuflow bench-mounted fume hood with side installation for work performed while seated

| Design characteristics | |
|---|---|
| Supporting construction | H-frame |
| Sash | 2 horizontal sashes |
| Side panel of the fume hood | Glass pane on the left and/or right as an option; not if service modules are installed in the side panel of the fume hood Material lock on the left and/or right as an option |
| Max. number of devices for scaffold points, \varnothing 12 to 13 mm | 12 |
| Max. load per scaffold point with scaffold rod length 300 mm [kg] | 5 |
| Service modules | Service modules in the left and/or right side panel of the fume hood, depending on requirement |

| Electrics | |
|--------------------|---|
| Electrical supply | External sockets in service panels Internal sockets in service modules |
| Fuse box | Optional |
| Sash controller SC | Optional |

| Sanitary technology | |
|---------------------|--|
| Sanitary supply | Service modules with take-off valves for vacuum, gases and/or waters and integrated sink (PP) as an option |

| Ventilation technology | |
|---|---|
| Minimum air exchange rate [m^3/h] ¹⁾ | 410 |
| Function display | FAZ |
| Airflow damper, constant | Airflow-Controller AC |
| Airflow damper, variable | Airflow-Controller AC |
| Detector of sash position | Only variable with Airflow-Controller AC |
| Connection height [mm] for FAZ with extract air hood \varnothing 250 mm | 2420 |
| Connection height [mm] for FAZ with extract air hood \varnothing 315 mm ²⁾ | 2420 |
| Connection height [mm] for AC with extract air hood \varnothing 250 mm | 2790 |
| Connection height [mm] for AC with extract air hood \varnothing 315 mm ²⁾ | 2760 |
| Underbench unit extraction system | As an option, depending on requirements and regulations |

¹⁾ All air volume specifications refer to an opening height of the sash window of 500 mm (test opening in acc. with EN 14175-3) and the maximum tracer gas values recommended by German Standard (BG Chemie). Shown rates correspond to a face velocity of 0.17 m/s (+/- 0.01 m/s). For other design face velocities, please contact your Waldner sales representative.

²⁾ In order to minimise noise and pressure losses, for air volumes >1000 m^3/h Waldner recommends using the extract air hood with a connection diameter of 315 mm.

A maximum admission pressure of 600 Pa should not be exceeded in the case of fume hoods with airflow dampers.

The indicated minimum air exchange rates were determined under specified test conditions in acc. with EN 14175-3. These minimum air exchange rates must be adapted when dimensioning the ventilation system.

If on-site extract air monitoring systems or airflow dampers are used, the required air volumes may be different. The operating limitations must be agreed upon with Waldner.

| Material/surface | |
|------------------|--|
| Worktop | Stoneware Polypropylene Epoxy Stainless steel |
| Internal lining | Solid grade laminate Melamine resin facing |

Bench-mounted fume hoods with side installation


Fume hood with side installation, made of steel

SI 3 steel

Intended use

- Protective equipment for users, tested in accordance with EN 14175 and ASHRAE 110.
- Extraction of fumes, aerosols and dust from the internal workspace to prevent dangerous concentrations of pollutants from escaping into the laboratory
- Reduced risk of the formation of a high concentration of hazardous substances / hazardous explosive atmosphere in the internal workspace
- Protection from splashes of hazardous substances
- Protection from flying particles, bodies or parts from the internal workspace
- General purpose fume hoods constructed in compliance with EN 14175 and ASHRAE 110 are normally not suitable for use with radioactive substances or micro-organisms
- Not suitable for openly breaking down chemicals
- Service outlets in the side panel of the internal workspace
- Control units located externally on the service panels

Design




- 1 Sash with sash handle
- 2 Worktop
- 3 Service panel
- 4 Side panel in fume hood wall
- 5 Function display control panel
- 6 Extract air spigot
- 7 Removable fascia panel
- 8 Upper sash window
- 9 Baffle with scaffold points
- 10 Bench frame

Bench-mounted fume hoods with side installation

Fume hood with side installation, made of steel

SI 3 steel

Dimensional drawing



Technical data

| Dimensions | 1200 | 1500 | 1800 | 2100 | 2400 |
|---------------------------------------|------|------|------|------|------|
| Width [mm] | 1200 | 1500 | 1800 | 2100 | 2400 |
| Depth [mm] | 900 | | | | |
| Height [mm] | 2400 | | | | |
| Clear width, internal workspace [mm] | 940 | 1240 | 1540 | 1840 | 2140 |
| Clear height, internal workspace [mm] | 1345 | | | | |
| Working height [mm] | 900 | | | | |

| Weight | 1200 | 1500 | 1800 | 2100 | 2400 |
|---------------------------|-------------|-------------|-------------|-------------|-------------|
| Without installation [kg] | Approx. 220 | Approx. 290 | Approx. 350 | Approx. 410 | Approx. 470 |

Bench-mounted fume hoods with side installation

Fume hood with side installation, made of steel

SI 3 steel

| Design characteristics | 1200 | 1500 | 1800 | 2100 | 2400 |
|--|--|------|---------------------|------|------|
| Supporting construction | H-frame with push-in underbench units | | | | |
| Sash | One-piece | | | | |
| Sash, divided (optional) | 2 horizontal sashes | | 3 horizontal sashes | | |
| Side panel of the fume hood | Without glazing and without equipment pass through hatch | | | | |
| Number of units for scaffold points, ø 12 to 13 mm | 6 | 6 | 6 | 8 | 10 |

| Electrics | |
|--------------------|--|
| Electrical supply | Sockets only external in service panel |
| Fuse box | Optional |
| Sash controller SC | Optional |

| Sanitary technology | |
|---------------------|---|
| Sanitary supply | Optional: Take-off valves for vacuum, gases and/or water and integrated sink (PP) in side panel |

| Ventilation technology | 1200 | 1500 | 1800 | 2100 | 2400 |
|---|--|------|------|------|------|
| EN 14175 minimum air exchange rate [m ³ /h] ¹⁾ | 380 | 460 | 500 | 650 | 750 |
| ASHRAE with 0.3 m/s / 60 fpm [m ³ /h] ²⁾ | 470 | 620 | 770 | 910 | 1060 |
| ASHRAE with 0.5 m/s / 100 fpm [m ³ /h] ³⁾ | 780 | 1030 | 1300 | 1520 | 1770 |
| Function display | FAZ / External control | | | | |
| Airflow damper, variable | Airflow-Controller AC | | | | |
| Connection height [mm] with function display with extract air spigot Ø 315 mm | 2420 | | | | |
| Connection height [mm] with AC with extract air spigot Ø 315 mm | 3262 | | | | |
| Floor/underbench extraction system | Optional depending on requirements and regulations | | | | |

¹⁾ Air volume specifications refer to an opening height of the sash window of 500 mm (test opening in line with EN 14175-3) and the maximum tracer gas values recommended by BG RCI.

²⁾ Air volume specifications refer to the prototype test in line with ASHRAE 110 with a face velocity of 60 fpm (0.3 m/s).

³⁾ Air volume specifications refer to the prototype test in line with ASHRAE 110 with a face velocity of 100 fpm (0.5 m/s).

The indicated minimum air exchange rates were determined under specified test conditions in compliance with EN 14175-3 and ASHRAE 110. Adapt these minimum air exchange rates when sizing the ventilation system.

The required air volumes may be different if on-site extract air monitoring systems or airflow dampers are used. Agree the operating limitations with Waldner.

| Material/surface | |
|------------------|---|
| Worktop | Epoxy, polypropylene, stainless steel |
| Internal lining | Polyresin, solid grade laminate, polypropylene, stainless steel |


Walk-in fume hoods

Walk-in fume hood with side installation

Intended use

- Protective device for the user, tested in acc. with EN 14175
- Extraction of fumes, aerosols and dust from the internal workspace to prevent dangerous amounts of pollutants from escaping into the laboratory
- Reduced risk of the formation of a high concentration of hazardous substances / hazardous explosive atmosphere in the internal workspace
- Protection from splashes of hazardous substances
- Protection from flying particles, bodies or parts escaping from the internal workspace
- General fume hoods constructed in acc. with EN 14175 are normally not suited for use with radioactive substances or microorganisms
- Not suitable for openly breaking down chemicals
- Suitable for barrier-free entering of the internal workspace
- Service outlets in the service modules of the side panels of the internal workspace
- Control units located vertically on the side service panels
- Suitable for high experimental set-ups

Design




- 1 Sash with handle and horizontal sashes
- 2 Service panel
- 3 Service module in the side panel of the fume hood
- 4 FAZ or AC control panel
- 5 Upper sash window
- 6 Removable fascia panel
- 7 Exhaust hood
- 8 Baffle with scaffold points

Walk-in fume hoods

Walk-in fume hood with side installation

Dimensional drawing



Technical data

| Dimensions | 1200 | 1500 | 1800 | 2100 | 2400 |
|---------------------------------------|------|------|------|------|------|
| Width [mm] | 1200 | 1500 | 1800 | 2100 | 2400 |
| Depth [mm] | 900 | | | | |
| Height [mm] | 2700 | | | | |
| Clear width, internal workspace [mm] | 950 | 1250 | 1550 | 1850 | 2150 |
| Clear height, internal workspace [mm] | 2450 | | | | |

| Weight | 1200 | 1500 | 1800 | 2100 | 2400 |
|---------------------------|-------------|-------------|-------------|-------------|-------------|
| Without installation [kg] | Approx. 320 | Approx. 390 | Approx. 450 | Approx. 510 | Approx. 570 |

Walk-in fume hoods

Walk-in fume hood with side installation

| Design characteristics | 1200 | 1500 | 1800 | 2100 | 2400 |
|---|---|------|---|------|------|
| Two-piece sash | 2 horizontal sashes at the top and bottom | | 3 horizontal sashes at the top and bottom | | |
| Side of fume hood | Glass pane on the left and/or right as an option; not if service modules are installed in the side panel of the fume hood Material lock on the left and/or right as an option | | | | |
| Number of devices for scaffold points, ø 12 to 13 mm | 9 | | 12 | | 15 |
| Max. load per scaffold point with scaffold rod length 300 mm [kg] | 5 | | | | |
| Service modules | In the left and/or right side panel of the fume hood, depending on requirement | | | | |

| Electrics | |
|--------------------|---|
| Electrical supply | External sockets in service panels Internal sockets in service modules |
| Fuse box | Optional |
| Sash controller SC | Optional |

| Sanitary technology | |
|---------------------|--|
| Sanitary supply | Service modules with take-off valves for vacuum, gases and/or waters and integrated sink (PP) as an option |

| Ventilation technology | 1200 | 1500 | 1800 | 2100 | 2400 |
|---|--|------|------|------|------|
| Minimum air exchange rate [m ³ /h] ¹⁾ | 480 | 600 | 720 | 840 | 960 |
| Function display | FAZ | | | | |
| Airflow damper, constant | Airflow-Controller AC | | | | |
| Airflow damper, variable | Airflow-Controller AC | | | | |
| Detector of sash position | Only variable with Airflow-Controller AC | | | | |
| Connection height [mm] for FAZ with extract air hood Ø 250 mm | 2570 | | | | |
| Connection height [mm] for FAZ with extract air hood Ø 315 mm ²⁾ | 2570 | | | | |
| Connection height [mm] for AC with extract air hood Ø 250 mm | 2940 | | | | |
| Connection height [mm] for AC with extract air hood Ø 315 mm ²⁾ | 2910 | | | | |

¹⁾ All air volume specifications refer to an opening height of the sash window of 500 mm (test opening in acc. with EN 14175-3) and the maximum tracer gas values recommended by German Standard (BG RCI). Shown rates correspond to a face velocity of 0.24 m/s. For other design face velocities, please contact your Waldner sales representative.

²⁾ In order to minimise noise and pressure losses, for air volumes >1000 m³/h Waldner recommends using the extract air hood with a connection diameter of 315 mm.

A maximum admission pressure of 600 Pa should not be exceeded in the case of fume hoods with airflow dampers.

The indicated minimum air exchange rates were determined under specified test conditions in acc. with EN 14175-3. These minimum air exchange rates must be adapted when dimensioning the ventilation system.

If on-site extract air monitoring systems or airflow dampers are used, the required air volumes may be different. The operating limitations must be agreed upon with Waldner.

| Material | |
|-----------------|---|
| Internal lining | Solid grade laminate Melamine resin facing |


Low level fume hoods

Low level fume hood with side installation

Intended use

- Protective device for the user, tested in acc. with EN 14175
- Extraction of fumes, aerosols and dust from the internal workspace to prevent dangerous amounts of pollutants from escaping into the laboratory
- Reduced risk of the formation of a high concentration of hazardous substances / hazardous explosive atmosphere in the internal workspace
- Protection from splashes of hazardous substances
- Protection from flying particles, bodies or parts escaping from the internal workspace
- General fume hoods constructed in acc. with EN 14175 are normally not suited for use with radioactive substances or microorganisms
- Not suitable for openly breaking down chemicals
- Suitable for experimental set-ups on an add-on table
- Service outlets in the service modules of the side panels of the internal workspace
- Control units located vertically on the side service panels

Design




- 1 Sash with handle and horizontal sashes
- 2 Service panel
- 3 Service module in the side panel of the fume hood
- 4 FAZ or AC control panel
- 5 Upper sash window
- 6 Removable fascia panel
- 7 Exhaust hood
- 8 Baffle with scaffold points
- 9 Add-on table

Low level fume hoods

Low level fume hood with side installation

Dimensional drawing



Technical data

| Dimensions | 1200 | 1500 | 1800 | 2100 | 2400 |
|---------------------------------------|-----------|------------|------------|------------|------------|
| Width [mm] | 1200 | 1500 | 1800 | 2100 | 2400 |
| Depth [mm] | 900 | | | | |
| Height [mm] | 2700 | | | | |
| Clear width, internal workspace [mm] | 950 | 1250 | 1550 | 1850 | 2150 |
| Clear height, internal workspace [mm] | 1950 | | | | |
| Add-on table with H-frame [mm] | 900 x 575 | 1200 x 575 | 1500 x 575 | 1800 x 575 | 2100 x 575 |
| Working height [mm] | 500 | | | | |

| Weight | 1200 | 1500 | 1800 | 2100 | 2400 |
|---------------------------|-------------|-------------|-------------|-------------|-------------|
| Without installation [kg] | Approx. 320 | Approx. 390 | Approx. 450 | Approx. 510 | Approx. 570 |

Low level fume hoods

Low level fume hood with side installation

| Design characteristics | 1200 | 1500 | 1800 | 2100 | 2400 |
|---|---|------|---|------|------|
| Work surface | Add-on table H-frame with surrounding increased edge | | | | |
| Two-piece sash | 2 horizontal sashes at the top and bottom | | 3 horizontal sashes at the top and bottom | | |
| Side of fume hood | Glass pane on the left and/or right as an option; not if service modules are installed in the side panel of the fume hood Material lock on the left and/or right as an option | | | | |
| Number of devices for scaffold points, ø 12 to 13 mm | 9 | | 12 | | 15 |
| Max. load per scaffold point with scaffold rod length 300 mm [kg] | 5 | | | | |
| Service modules | Service modules in the left and/or right side panel of the fume hood, depending on requirement | | | | |

| Electrics | |
|--------------------|---|
| Electrical supply | External sockets in service panels Internal sockets in service modules |
| Fuse box | Optional |
| Sash controller SC | Optional |

| Sanitary technology | |
|---------------------|--|
| Sanitary supply | Service modules with take-off valves for vacuum, gases and/or waters and integrated sink (PP) as an option |

| Ventilation technology | 1200 | 1500 | 1800 | 2100 | 2400 |
|---|---|------|------|------|------|
| Minimum air exchange rate [m³/h] ¹⁾ | 480 | 600 | 720 | 840 | 960 |
| Function display | FAZ | | | | |
| Airflow damper, constant | Airflow-Controller AC | | | | |
| Airflow damper, variable | Airflow-Controller AC | | | | |
| Detector of sash position | Only variable with Airflow-Controller AC | | | | |
| Connection height [mm] for FAZ with extract air hood ø 250 mm | 2570 | | | | |
| Connection height [mm] for FAZ with extract air hood ø 315 mm ²⁾ | 2570 | | | | |
| Connection height [mm] for AC with extract air hood ø 250 mm | 2940 | | | | |
| Connection height [mm] for AC with extract air hood ø 315 mm ²⁾ | 2910 | | | | |
| Underbench unit extraction system | As an option, depending on requirements and regulations | | | | |

¹⁾ All air volume specifications refer to an opening height of the sash window of 500 mm (test opening in acc. with EN 14175-3) and the maximum tracer gas values recommended by German Standard (BG RCI). Shown rates correspond to a face velocity of 0.24 m/s. For other design face velocities, please contact your Waldner sales representative.

²⁾ In order to minimise noise and pressure losses, for air volumes >1000 m³/h Waldner recommends using the extract air hood with a connection diameter of 315 mm.

A maximum admission pressure of 600 Pa should not be exceeded in the case of fume hoods with airflow dampers.

The indicated minimum air exchange rates were determined under specified test conditions in acc. with EN 14175-3. These minimum air exchange rates must be adapted when dimensioning the ventilation system.

If on-site extract air monitoring systems or airflow dampers are used, the required air volumes may be different. The operating limitations must be agreed upon with Waldner.

| Material | |
|---|--|
| Worktop H-frame with surrounding increased edge | Polypropylene Epoxy Stoneware Stainless steel |
| Internal lining | Solid grade laminate Melamine resin facing |


Special fume hoods

Secuflow EN7 fume hood for high thermal loads

Intended use

- For working with high thermal loads in the internal workspace (heat sources of 4 KW per metre of clear width in the fume hood)
- Protective device for the user, tested in accordance with DIN EN 14175-7
- Extraction of fumes, aerosols and dust from the internal workspace to prevent dangerous amounts of pollutants from escaping into the laboratory
- Reduced risk of the formation of a high concentration of hazardous substances / hazardous explosive atmosphere in the internal workspace
- Protection from splashes of hazardous substances
- Protection from flying particles and compounds from the internal workspace
- Fume hoods which are built in accordance with EN 14175-7 are not permitted for working with radioactive substances and microorganisms
- Not suitable for openly breaking down chemicals
- Active supportive flow technology (Secuflow technology) reduces the energy consumption while regulations and standards are observed
- Service outlets for sanitary supply in the rear panel of the internal workspace
- Control units located on the exterior of the support

Design




- 1 Sash with handle bar and horizontal sashes
- 2 Worktop
- 3 FAZ or AC control panel
- 4 Removable fascia panel
- 5 Exhaust hood
- 6 Bench frame with push-in underbench units with support and service panels

Special fume hoods

Secuflow EN7 fume hood for high thermal loads

Dimensional drawing



Technical data

| Dimensions | 1200 | 1500 | 1800 |
|--|------|------|------|
| Width [mm] | 1200 | 1500 | 1800 |
| Depth [mm] | 900 | | |
| Height [mm] | 2700 | | |
| Clear internal width of internal space [mm] | 1150 | 1450 | 1750 |
| Clear internal height of internal space [mm] | 1550 | | |
| Working height [mm] | 900 | | |

| Weight | 1200 | 1500 | 1800 |
|---------------------------|-------------|-------------|-------------|
| Without installation [kg] | Approx. 250 | Approx. 300 | Approx. 350 |

Special fume hoods

Secuflow EN7 fume hood for high thermal loads

| Design characteristics | 1200 | 1500 | 1800 |
|--|---------------------------------------|------|---------------------|
| Supporting construction | H-frame with push-in underbench units | | |
| Sash | 2 horizontal sashes | | 3 horizontal sashes |
| Side panel of the fume hood | full | | |
| Max. number of devices for scaffold points, dia. 12 mm up to 13 mm | 9 | | 12 |
| Max. load per scaffold point with scaffold rod length 300 mm [kg] | 5 | | |
| Service modules | 2 | | 3 |

| Electrics | |
|--------------------|------------------------------------|
| Electrical supply | External sockets in service panels |
| Fuse box | Optional |
| Sash controller SC | Optional |

| Sanitary technology | |
|---------------------|--|
| Sanitary supply | Service modules with take-off valves for vacuum, gases and/or waters and integrated sink (PP) as an option |

| Ventilation technology | 1200 | 1500 | 1800 |
|--|---|---------|---------|
| Air flow range without/with thermal load [m ³ /h] ¹⁾ | 450/700 | 450/750 | 540/900 |
| Function display with temperature monitoring | FAZ | | |
| Airflow damper, constant with temperature monitoring | Airflow-Controller AC | | |
| Airflow damper, variable with temperature control | Airflow-Controller AC | | |
| Detector of sash position | Only variable with Airflow-Controller AC | | |
| Connection height [mm] for FAZ with extract air hood dia. 250 mm | 2570 | | |
| Connection height [mm] for FAZ with extract air hood dia. 315 mm ²⁾ | 2570 | | |
| Connection height [mm] for AC with extract air hood dia. 250 mm | 2940 | | |
| Connection height [mm] for AC with extract air hood dia. 315 mm ²⁾ | 2910 | | |
| Underbench unit extraction system | As an option, depending on requirements and regulations | | |

¹⁾ All air volume specifications refer to an opening height of the sash window of 500 mm (Test opening in accordance with EN14175-3) and the recommended maximum trace gas values of BG RCI.

²⁾ In order to minimise noise and pressure losses, for air volumes >1000 m³/h Waldner recommends the extract air hood with a connection diameter of 315 mm.

A maximum admission pressure of 600 Pa should not be exceeded with fume hoods with airflow dampers.

The indicated minimum air exchange rates were determined under specified test conditions in acc. with EN 14175-3. These minimum air exchange rates must be adapted when dimensioning the ventilation system.

If on-site extract air monitoring systems or airflow dampers are used, the required air volumes may be different. The operating limitations must be agreed upon with Waldner.

| Material/surface | |
|------------------|--|
| Worktop | Stoneware Polypropylene Stainless steel Epoxy |
| Internal lining | Melamine resin facing Solid (grade) laminate Stoneware |


Special fume hoods

EN7 fume hood for high thermal loads in connection with acid digestions (special application fume hood)

Intended use

- Protective device for the user, tested in accordance with DIN EN 14175-7
- Suitable for open, thermal processes of breaking down chemicals with aggressive media such as, e.g., sulphuric acid, hydrochloric acid or aqua regia
- The construction of the fume hood and the materials of the inner lining of the internal workspace determine which aggressive media the device can be used for
- Extraction of fumes and aerosols from the internal workspace to prevent dangerous amounts of pollutants from escaping into the laboratory
- Reduced risk of the formation of a high concentration of hazardous substances / hazardous explosive atmosphere in the internal workspace
- Protection from splashes of hazardous substances in the internal workspace
- Protection from flying particles, bodies or parts from the internal workspace
- Fume hoods which are built in accordance with EN 14175-7 are not permitted for working with radioactive substances and microorganisms
- For working with high thermal loads in connection with acid digestions in the internal workspace (heat sources of 4 KW per metre of clear width in the fume hood)

Design




- 1 Sash with handle bar
- 2 Worktop
- 3 FAZ or AC control panel
- 4 Removable fascia panel
- 5 Extract air spigot integrated in fume-scrubber (as an option)
- 6 Baffle
- 7 Bench frame with push-in underbench units with support and service panels

Special fume hoods

EN7 fume hood for high thermal loads in connection with acid digestions (special application fume hood)

Dimensional drawing



Technical data

| Dimensions | 1200 | 1500 | 1800 |
|--|------|------|------|
| Width [mm] | 1200 | 1500 | 1800 |
| Depth [mm] | 900 | | |
| Height [mm] | 2700 | | |
| Clear internal width of internal space [mm] | 1150 | 1450 | 1750 |
| Clear internal height of internal space [mm] | 1060 | | |
| Working height [mm] | 900 | | |

| Weight | 1200 | 1500 | 1800 |
|--|--------------|-------------|-------------|
| Without installations and fume-scrubber [kg] | Approx. 250 | Approx. 300 | Approx. 350 |
| Fume-scrubber without filling [kg] | 55 (MAALS 1) | | |
| | 63 (MAALS 2) | | |

Special fume hoods

EN7 fume hood for high thermal loads in connection with acid digestions (special application fume hood)

| Design characteristics | |
|-------------------------------------|---------------------------------------|
| Supporting construction | H-frame with push-in underbench units |
| Extract manifold | Standard |
| Fume-scrubber | Optional |
| Neutralisation unit underbench unit | Optional |

| Electrics | |
|--------------------|------------------------------------|
| Electrical supply | External sockets in service panels |
| Fuse box | Optional |
| Sash controller SC | Optional |

| Sanitary technology | |
|---------------------|--|
| Sanitary supply | With take-off valves for vacuum, gases and/or waters and drip cup integrated in the worktop as an option |

| Ventilation technology | 1200 | 1500 | 1800 |
|--|---|---------|---------|
| Minimum air exchange rate [m ³ /h] ¹⁾ | 650 | 800 | 950 |
| Pressure loss in extract manifold with FAZ/ AC [Pa] | 45/120 | 50/120 | 85/150 |
| Pressure loss in fume hood with fume-scrubber [Pa] | 670/740 | 780/840 | 860/920 |
| Müller fume-scrubber type | MAALS 1 | | MAALS 2 |
| Function display with temperature monitoring | FAZ | | |
| Airflow damper, constant with temperature monitoring | Airflow-Controller AC | | |
| Connection height [mm] for FAZ and AC with extract air spigot dia. 250 mm with fume-scrubber | 3040 | | |
| Connection height [mm] for FAZ with extract manifold dia. 250 mm (with no fume-scrubber) | 2380 | | |
| Connection height [mm] for AC with extract manifold dia. 250 mm (with no fume-scrubber) | 2760 | | |
| Underbench unit extraction system | As an option, depending on requirements and regulations | | |

¹⁾ All air volume specifications refer to an opening height of the sash window of 500 mm (Test opening in accordance with EN 14175-3) and the recommended maximum trace gas values of BG RCI.

A maximum admission pressure of 600 Pa should not be exceeded with fume hoods with airflow dampers.

The indicated minimum air exchange rates were determined under specified test conditions in acc. with EN 14175-3. These minimum air exchange rates must be adapted when dimensioning the ventilation system.

If on-site extract air monitoring systems or airflow dampers are used, the required air volumes may be different. The operating limitations must be agreed upon with Waldner.

| Material/surface | |
|-----------------------------------|---|
| Internal lining including worktop | Stoneware (when sulphuric acid, hydrochloric acid, aqua regia are used) |


Special fume hoods

Fume hood for handling of perchloric acid

Intended use

- Protective device for the user, tested in accordance with DIN EN 14175
- Suitable for open, thermal processes of breaking down chemicals with aggressive media, particularly for perchloric acid
- The construction of the fume hood and the materials of the inner lining of the internal workspace determine which aggressive media the device can be used for
- Extraction of fumes and aerosols from the internal workspace to prevent dangerous amounts of pollutants from escaping into the laboratory
- Reduced risk of the formation of a high concentration of hazardous substances / hazardous explosive atmosphere in the internal workspace
- Protection from splashes of hazardous substances in the internal workspace
- Protection from flying particles, bodies or parts from the internal workspace
- Fume hoods which are built in accordance with EN 14175-7 are not permitted for working with radioactive substances and microorganisms
- For working with high thermal loads in connection with acid digestions in the internal workspace (heat sources of 4 KW per metre of clear width in the fume hood).

Design




- 1 Sash with handle bar
- 2 Worktop
- 3 FAZ or AC control panel
- 4 Removable fascia panel
- 5 Extract air spigot integrated in fume-scrubber (as an option)
- 6 Baffle
- 7 Bench frame with push-in underbench units with support and service panels

Special fume hoods

Fume hood for handling of perchloric acid

Dimensional drawing



Technical data

| Dimensions | 1200 | 1500 | 1800 |
|--|------|------|------|
| Width [mm] | 1200 | 1500 | 1800 |
| Depth [mm] | 900 | | |
| Height [mm] | 2700 | | |
| Clear internal width of internal space [mm] | 1150 | 1450 | 1750 |
| Clear internal height of internal space [mm] | 1060 | | |
| Working height [mm] | 900 | | |

| Weight | 1200 | 1500 | 1800 |
|--|--------------|-------------|--------------|
| Without installations and fume-scrubber [kg] | Approx. 250 | Approx. 300 | Approx. 350 |
| Fume-scrubber without filling [kg] | 55 (MAALS 1) | | 63 (MAALS 2) |

Special fume hoods

Fume hood for handling of perchloric acid

| Design characteristics | |
|-------------------------------------|--|
| Supporting construction | H-frame with push-in underbench units |
| Fume-scrubber | Optional |
| Extract manifold with sprinkler | Optional (only for fume hoods for handling of perchloric acid) |
| Neutralisation unit underbench unit | Optional |

| Electrics | |
|--------------------|------------------------------------|
| Electrical supply | External sockets in service panels |
| Fuse box | Optional |
| Sash controller SC | Optional |

| Sanitary technology | |
|---------------------|--|
| Sanitary supply | With take-off valves for vacuum, gases and/or waters and drip cup integrated in the worktop as an option |

| Ventilation technology | 1200 | 1500 | 1800 |
|---|---|---------|---------|
| Minimum air exchange rate [m ³ /h] ¹⁾ | 650 | 800 | 950 |
| Pressure loss in extract manifold with sprinkler for FAZ/AC [Pa] | 140/300 | 160/350 | 270/500 |
| Pressure loss in extract manifold with FAZ/AC [Pa] | 45/120 | 50/120 | 85/150 |
| Pressure loss in fume hood with fume-scrubber [Pa] | 670/740 | 780/840 | 860/920 |
| Müller fume-scrubber type | MAALS 1 | | MAALS 2 |
| Function display with temperature monitoring | FAZ | | |
| Airflow damper, constant with temperature monitoring | Airflow-Controller AC | | |
| Connection height [mm] for FAZ and AC with extract air spigot dia. 250 mm with fume-scrubber) | 3040 | | |
| Connection height [mm] with FAC/AC with extract manifold and sprinkler | 2430 / 2810 | | |
| Underbench unit extraction system | As an option, depending on requirements and regulations | | |

¹⁾ All air volume specifications refer to an opening height of the sash window of 500 mm (Test opening in accordance with EN14175-3) and the recommended maximum trace gas values of BG RCI.

A maximum admission pressure of 600 Pa should not be exceeded with fume hoods with airflow dampers.

The indicated minimum air exchange rates were determined under specified test conditions in acc. with EN 14175-3. These minimum air exchange rates must be adapted when dimensioning the ventilation system.

If on-site extract air monitoring systems or airflow dampers are used, the required air volumes may be different. The operating limitations must be agreed upon with Waldner.

| Material/surface | |
|-----------------------------------|---|
| Internal lining including worktop | Stoneware (when sulphuric acid, hydrochloric acid, aqua regia are used) |


Special fume hoods

Fume hood for handling of hydrofluoric acid

Intended use

- Protective device for the user, tested in accordance with DIN EN 14175-7
- Suitable for open, thermal processes of breaking down chemicals with aggressive media, particularly for hydrofluoric acid
- The construction of the fume hood and the materials of the inner lining of the internal workspace determine which aggressive media the device can be used for
- Extraction of fumes and aerosols from the internal workspace to prevent dangerous amounts of pollutants from escaping into the laboratory
- Reduced risk of the formation of a high concentration of hazardous substances / hazardous explosive atmosphere in the internal workspace
- Protection from splashes of hazardous substances in the internal workspace
- Protection from flying particles, bodies or parts from the internal workspace
- Fume hoods which are built in accordance with EN 14175-7 are not permitted for working with radioactive substances and microorganisms
- For working with high thermal loads in connection with acid digestions in the internal workspace (heat sources of 4 KW per metre of clear width in the fume hood)

Design




- 1 Sash with handle bar
- 2 Worktop
- 3 FAZ or AC control panel
- 4 Removable fascia panel
- 5 Extract air spigot integrated in fume-scrubber (as an option)
- 6 Baffle
- 7 Bench frame with push-in underbench units with support and service panels

Special fume hoods

Fume hood for handling of hydrofluoric acid

Dimensional drawing



Technical data

| Dimensions | 1200 | 1500 | 1800 |
|--|------|------|------|
| Width [mm] | 1200 | 1500 | 1800 |
| Depth [mm] | 900 | | |
| Height [mm] | 2700 | | |
| Clear internal width of internal space [mm] | 1150 | 1450 | 1750 |
| Clear internal height of internal space [mm] | 1060 | | |
| Working height [mm] | 900 | | |

| Weight | 1200 | 1500 | 1800 |
|--|--------------|-------------|--------------|
| Without installations and fume-scrubber [kg] | Approx. 250 | Approx. 300 | Approx. 350 |
| Fume-scrubber without filling [kg] | 55 (MAALS 1) | | 63 (MAALS 2) |

Special fume hoods

Fume hood for handling of hydrofluoric acid

| Design characteristics | |
|-------------------------------------|---------------------------------------|
| Supporting construction | H-frame with push-in underbench units |
| Extract manifold | Standard |
| Fume-scrubber | Optional |
| Neutralisation unit underbench unit | Optional |

| Electrics | |
|--------------------|------------------------------------|
| Electrical supply | External sockets in service panels |
| Fuse box | Optional |
| Sash controller SC | Optional |

| Sanitary technology | |
|---------------------|--|
| Sanitary supply | With take-off valves for vacuum, gases and/or waters and drip cup integrated in the worktop as an option |

| Ventilation technology | 1200 | 1500 | 1800 |
|--|---|---------|---------|
| Minimum air exchange rate [m ³ /h] ¹⁾ | 650 | 800 | 950 |
| Pressure loss in extract manifold with FAZ/ AC [Pa] | 45/120 | 50/120 | 85/150 |
| Pressure loss in fume hood with fume-scrubber [Pa] | 670/740 | 780/840 | 860/920 |
| Müller fume-scrubber type | MAALS 1 | | MAALS 2 |
| Function display with temperature monitoring | FAZ | | |
| Airflow damper, constant with temperature monitoring | Airflow-Controller AC | | |
| Connection height [mm] for FAZ and AC with extract air spigot dia. 250 mm with fume-scrubber | 3040 | | |
| Connection height [mm] for FAZ with extract manifold dia. 250 mm (with no fume-scrubber) | 2380 | | |
| Connection height [mm] for AC with extract manifold dia. 250 mm (with no fume-scrubber) | 2760 | | |
| Underbench unit extraction system | As an option, depending on requirements and regulations | | |

¹⁾ All air volume specifications refer to an opening height of the sash window of 500 mm (Test opening in accordance with EN 14175-3) and the recommended maximum trace gas values of BG RCI.

A maximum admission pressure of 600 Pa should not be exceeded with fume hoods with airflow dampers.

The indicated minimum air exchange rates were determined under specified test conditions in acc. with EN 14175-3. These minimum air exchange rates must be adapted when dimensioning the ventilation system.

If on-site extract air monitoring systems or airflow dampers are used, the required air volumes may be different. The operating limitations must be agreed upon with Waldner.

| Material/surface | |
|-----------------------------------|--|
| Internal lining including worktop | Polypropylene (when hydrofluoric acid is used) |


Special fume hoods

Radio-isotope fume hood

Intended use

- Protective device for the user, tested in accordance with DIN EN 14175-8
- Extraction during work with radioactive substances if increased requirements for radiation protection apply
- Protection from incorporation, contamination and external radiation exposure
- Extraction of fumes, aerosols and dust from the internal workspace to prevent dangerous amounts of pollutants from escaping into the laboratory
- Reduced risk of the formation of a high concentration of hazardous substances / hazardous explosive atmosphere in the internal workspace
- Protection from splashes of hazardous substances in the internal workspace
- Protection from flying particles, bodies or parts escaping from the internal workspace
- Fume hoods which are built in accordance with DIN EN 14175-8 are normally not permitted for use with microorganisms
- Not suitable for openly breaking down chemicals

Design




- 1 Differential pressure gauge
- 2 Sash with handle
- 3 Worktop
- 4 FAZ or AC control panel
- 5 Removable fascia panel
- 6 Extract air spigot integrated in filter housing
- 7 Baffle with scaffold points
- 8 Bench frame with push-in underbench units with support and service panels

Special fume hoods Radio-isotope fume hood

1

Fume hoods and extraction devices

Dimensional drawing



Technical data

| Dimensions | 1200 | 1500 |
|---|-----------------|------|
| Width [mm] | 1200 | 1500 |
| Depth [mm] | 900 | |
| Height [mm] | 2700 | |
| Clear width, internal workspace [mm] | 1150 | 1450 |
| Clear height, internal workspace [mm] | 1053 | |
| Working height [mm] | 900 | |
| Filter housing, width x depth x height [mm] | 820 x 775 x 674 | |

| Weight | 1200 | 1500 |
|--|-------------|-------------|
| Without installations and lead insert [kg] | Approx. 250 | Approx. 300 |
| Filter housing [kg] | 90 | |

Special fume hoods

Radio-isotope fume hood

| Design characteristics | |
|--|--|
| Supporting construction | H-frame with push-in underbench units |
| Sash | One-piece |
| Number of devices for scaffold points, ø 12 to 13 mm | 6 |
| Max. load per scaffold point with scaffold rod length 300 mm [kg] | 5 |
| Filter, fume hood roof | Standard equipment: Filter F7 / particle filter H13 |
| Filter, lateral cabinet (max. 3 filter housings) | Filter housing, top: Particulate filter Filter housing, centre: Active charcoal filter Filter housing, bottom: Filter and particle filter |
| Differential pressure gauges | Display of the degree of saturation of the filters (not for active charcoal filter) |
| Lead insert | Optional |
| Waste disposal system for radio-isotope residual material in the underbench unit | Canister for liquid radio-isotope residual material as an option Collapsible boxes for solid radio-isotope residual material as an option Level indicator and/or opening in the worktop as an option |

| Electrics | |
|--------------------|------------------------------------|
| Electrical supply | External sockets in service panels |
| Fuse box | Optional |
| Sash controller SC | Optional |

| Sanitary technology | |
|---------------------|--|
| Sanitary supply | With take-off valves for vacuum and gases as an option |

| Ventilation technology | 1200 | 1500 |
|--|---|--------|
| Minimum air exchange rate [m ³ /h] ¹⁾ | 480 | 650 |
| Pressure loss, filter [Pa] ²⁾ | 25/200 | 30/235 |
| Pressure loss, particle filter [Pa] ²⁾ | 50/300 | 60/350 |
| Pressure loss, active charcoal filter [Pa] ²⁾ | 25/25 | 30/30 |
| Pressure loss, particulate filter [Pa] ²⁾ | 30/250 | 35/290 |
| Function display | FAZ | |
| Airflow damper, constant | Airflow-Controller AC | |
| Airflow damper, variable | Airflow-Controller AC | |
| Connection height [mm] for FAZ and AC with filter housing Ø 250 mm | 3100 | |
| Underbench unit extraction system | As an option, depending on requirements and regulations | |

¹⁾ All air volume specifications refer to an opening height of the sash window of 500 mm (test opening in acc. with EN 14175-3) and the maximum tracer gas values recommended by German Standard (BG RCI).

²⁾ Pressure loss values refer to the states clean/contaminated.

A maximum admission pressure of 600 Pa should not be exceeded in the case of fume hoods with airflow dampers.

The indicated minimum air exchange rates were determined under specified test conditions in acc. with EN 14175-3. These minimum air exchange rates must be adapted when dimensioning the ventilation system.

If on-site extract air monitoring systems or airflow dampers are used, the required air volumes may be different. The operating limitations must be agreed upon with Waldner.

In the case of fume hoods with filters, the pressure loss of the integrated filter stages must be added to the pressure loss of the fume hood.

| Material/surface | |
|-----------------------------------|----------------------------------|
| Internal lining including worktop | Polypropylene Stainless steel |

Special fume hoods

Radio-isotope fume hood

1

Fume hoods and extraction devices

| Filter (filter in the filter cabinet or on the fume hood roof) | |
|---|--|
| Dimensions [mm] | 610 x 610 x 46 (+ 8 mm seal) |
| Pressure loss [Pa] at 1900 m ³ /h | 110 |
| Design characteristics | Filter element (fine particle filter); filter class EN 779: F7 Frame made of multilayer board with grip and type label on the 610-mm side; PU seal on the dust-laden air side |
| Use | Fine particle filter for particle adsorption, e.g.: Oil smoke and agglomerated soot, tobacco smoke, metal oxide smoke Average efficiency (Em) 80–90% |

| Particle filter (filter in the filter cabinet or on the fume hood roof) | |
|--|---|
| Dimensions [mm] | 610 x 610 x 292 (+ 7 mm seal) |
| Pressure loss [Pa] at 2435 m ³ /h | 250 |
| Design characteristics | Particle filter element type: Hepa H13; efficiency: MPPS Frame made of multilayer board with grip and type label on the 610-mm side; PU tight seat seal on the clean air side; filter medium flush on the clean air side |
| Use | Particle filter for the adsorption of particles up to H13; particle adsorption 99.95 %; transmittance 0.05% |

| Active charcoal filter (filter in the filter cabinet) | |
|--|---|
| Dimensions [mm] | 610 x 610 x 292 (+ 7 mm seal) |
| Pressure loss [Pa] at 600 m ³ /h | 9 |
| Design characteristics | Activated charcoal cell 7C for 16 x activated charcoal cartridges Frame galvanised sheet metal; 2 x U handle and type plate on the 610-mm side; PU tight seat seal on the clean air side |
| Use | Standard impregnation: for all common radioactive materials, radioactive iodine compounds, radioactive iodomethane, radioactive gases. (A filter with filters class F7 in acc. with EN 779 is recommended.) |

| Particulate filter (filter in the filter cabinet) | |
|--|--|
| Dimensions [mm] | 610 x 610 x 292 (+ 7 mm seal) |
| Pressure loss [Pa] at 1965 m ³ /h | 125 |
| Design characteristics | Particulate or Microtain filter element type: Hepa H11 in acc. with EN 1822 Frame made of multilayer board with grip and type label on the 610-mm side; PU tight seat seal on the clean air side; filter medium flush on the clean air side |
| Use | Particle filter for the adsorption of particles up to H11; particle adsorption 95 %; transmittance 5%; to be installed after active charcoal filters to bind the charcoal dust contamination from the charcoal filter. |


Special fume hoods

Filter fume hood

Intended use


- Before the extract air from the internal workspace is released into the environment, it is cleaned by a filter unit

Design



- 1 Differential pressure gauge
- 2 Sash with handle and horizontal sashes
- 3 Worktop
- 4 FAZ or AC control panel
- 5 Removable fascia panel
- 6 Extract air spigot
- 7 Baffle with scaffold points
- 8 Bench frame with push-in underbench units with support and service panels

Dimensional drawing



Technical data

| Dimensions | 1200 | 1500 | 1800 |
|---|-----------------|------|------|
| Width [mm] | 1200 | 1500 | 1800 |
| Depth [mm] | 900 | | |
| Height [mm] | 2700 | | |
| Clear width, internal workspace [mm] | 1150 | 1450 | 1750 |
| Clear height, internal workspace [mm] | 1060 | | |
| Working height [mm] | 900 | | |
| Filter housing, width x depth x height [mm] | 820 x 775 x 674 | | |

| Weight | 1200 | 1500 | 1800 |
|---|-------------|-------------|-------------|
| Filter fume hood without installations [kg] | Approx. 270 | Approx. 320 | Approx. 370 |
| Filter housing [kg] | 90 | | |

Special fume hoods

Filter fume hood

| Design characteristics | 1200 | 1500 | 1800 |
|--|--|------|---------------------|
| Supporting construction | H-frame with push-in underbench units | | |
| Sash | 2 horizontal sashes | | 3 horizontal sashes |
| Glass pane in the side wall | Possible on the left and/or right side of the fume hood; not with stoneware internal lining | | |
| Number of devices for scaffold points, ø 12 to 13 mm | 6 | | 8 |
| Max. load per scaffold point with scaffold rod length 300 mm [kg] | 5 | | |
| Material lock | Possible on the left and/or right side of the fume hood | | |
| Filter, fume hood roof | Standard equipment: Filter F7 / particle filter H13 | | |
| Differential pressure gauges | Display of the degree of saturation of the filters | | |

| Electrics | |
|--------------------|------------------------------------|
| Electrical supply | External sockets in service panels |
| Fuse box | Optional |
| Sash controller SC | Optional |

| Sanitary technology | |
|---------------------|--|
| Sanitary supply | With take-off valves for vacuum, gases and/or waters and drip cup integrated in the worktop as an option |

| Ventilation technology | 1200 | 1500 | 1800 |
|---|---|--------|---------|
| Minimum air exchange rate [m ³ /h] ¹⁾ | 480 | 600 | 720 |
| Pressure loss, filter [Pa] ²⁾ | 35/200 | 45/235 | 65/290 |
| Pressure loss, particle filter [Pa] ²⁾ | 70/300 | 95/365 | 130/430 |
| Function display | FAZ | | |
| Airflow damper, constant | Airflow-Controller AC | | |
| Airflow damper, variable | Airflow-Controller AC | | |
| Detector of sash position | Only variable with Airflow-Controller AC | | |
| Connection height [mm] for FAZ and AC with extract air spigot Ø 250 mm | 3100 | | |
| Underbench unit extraction system | As an option, depending on requirements and regulations | | |

¹⁾ All air volume specifications refer to an opening height of the sash window of 500 mm (test opening in acc. with EN 14175-3) and the maximum tracer gas values recommended by German Standard (BG RCI).

²⁾ Pressure loss values refer to the states clean/contaminated.

A maximum admission pressure of 600 Pa should not be exceeded in the case of fume hoods with airflow dampers.

The indicated minimum air exchange rates were determined under specified test conditions in acc. with EN 14175-3. These minimum air exchange rates must be adapted when dimensioning the ventilation system.

If on-site extract air monitoring systems or airflow dampers are used, the required air volumes may be different. The operating limitations must be agreed upon with Waldner.

In the case of fume hoods with filters, the pressure loss of the integrated filter stages must be added to the pressure loss of the fume hood.

| Material/surface | |
|------------------|--|
| Worktop | Stoneware Polypropylene Epoxy Stainless steel |
| Internal lining | Melamine resin facing Solid grade laminate |

Special fume hoods

Filter fume hood

1

| Filter | |
|--|--|
| Dimensions [mm] | 610 x 610 x 46 (+ 8 mm seal) |
| Pressure loss [Pa] at 1900 m ³ /h | 110 |
| Design characteristics | Filter element (fine particle filter); filter class EN 779: F7 Frame made of multilayer board with grip and type label on the 610-mm side; PU seal on the dust-laden air side |
| Use | Fine particle filter for particle adsorption, e.g.: Oil smoke and agglomerated soot, tobacco smoke, metal oxide smoke Average efficiency (Em) 80–90% |

| Particle filter | |
|--|---|
| Dimensions [mm] | 610 x 610 x 292 (+ 7 mm seal) |
| Pressure loss [Pa] at 2435 m ³ /h | 250 |
| Design characteristics | Particle filter element type: Hepa H13; efficiency: MPPS Frame made of multilayer board with grip and type label on the 610-mm side; PU tight seat seal on the clean air side; filter medium flush on the clean air side |
| Use | Particle filter for the adsorption of particles up to H13; particle adsorption 99.95 %; transmittance 0.05% |

Fume hoods and extraction devices


Safety weighing cabinet

AKKURAT

Intended use


- Protection against airborne particles or aerosols, which can be released when weighing powdered or dusty substances
- Accommodation of high-resolution laboratory balances from all popular manufacturers
- Fully vibration-free environment free from draughts, which micro-balances require
- Precise and safe working with highly active substances
- Completely vibration-decoupled and impact-dampened worktop
- A design focused on functionality and ergonomics meets the most exacting quality and safety requirements
- All functions are at your fingertips and can be selected using the integral control panel
- Air foil cill for optimum supply air flow, which doubles as an ergonomic armrest
- Disposal of any waste produced directly via the waste disposal system integrated in the worktop
- Arrangement of the filter and separate vibration-decoupled fan unit guarantees maximum legroom
- The fully welded steel support frame sub-structure, connected to the vibration-decoupled bearing of the solid ceramic worktop, guarantees fault-free operation, even when using highly sensitive micro-balances

Construction



- 1 Fan unit incl. second HEPA filter
- 2 Filter box with HEPA filter
- 3 Fume hood function display control panel
- 4 Double-walled rear wall with two cable through-puts
- 5 Plexiglass head unit and folding front screen
- 6 Hinge with integral lighting
- 7 Vibration-decoupled worktop
- 8 Ionisation bar
- 9 Waste disposal system
- 10 Welded steel support frame

Dimensional drawing



Technical data

| Dimensions | 1200 | 1500 |
|--|------|------|
| Width [mm] | 1200 | 1500 |
| Depth [mm] | 650 | |
| Height [mm] | 1450 | |
| Clear width of internal space [mm] | 1155 | 1455 |
| Effective depth of internal space [mm] | 575 | |
| Clear internal height of internal space [mm] | 530 | |
| Working height [mm] | 900 | |

| Weight | 1200 | 1500 |
|---------------------------|------|------|
| Without installation [kg] | 120 | 135 |

Safety weighing cabinet AKKURAT

| Design characteristics | 1200 | 1500 |
|----------------------------|---|------|
| Supporting construction | Steel support frame, fully welded with height-adjustable feet | |
| Desktop | Technical ceramic with raised edge on all sides | |
| Upper part | Fully plexiglass upper part Upwardly-hinged front screen LED lighting integrated into the hinge for the front screen Rear wall with integral exhaust system and 2 cable through-puts | |
| Rear wall | Double-walled, PP white, UV-resistant | |
| Ventilation | Fan in separate housing with frequency converter for continuously variable speed control | |
| Mean and second filter box | HEPA H14 particle filter as per DIN EN 1822 filter exchange with almost no contamination due to the bag exchange method | |
| Ionisation | The separately switched ionising system, which is fully integrated in the balance enclosure, neutralises the non-conductive surfaces in the interior workspace of the balance enclosure and diverts electrostatic charges from the samples. The capacitively working high-voltage electrodes in the inflow area of the front opening are designed to be touch-safe. | |

| Display / Operation | |
|---------------------|--|
| Control panel | Control panel with 5 illuminated keys integrated flush with the air foil cill, under which the flow passes, on the front edge of the operating opening |
| Displays | <ul style="list-style-type: none"> - Equipment On/Off - Lighting On/Off - Air flow monitoring – alarm triggered when the limit is transgressed - Front opening monitoring – alarm triggered when the front screen is opened - Filter change (time-dependent) - Colour temperature LED lighting – switches between daylight white and neutral white |
| Operation | <ul style="list-style-type: none"> - Equipment On/Off - Lighting On/Off - Colour temperature lighting - Air flow alarm acknowledgement - Front opening alarm acknowledgement - Filter change alarm acknowledgement |

| Handling | |
|--------------------|---|
| Working substances | Powdered or dusty substances (e.g active ingredients for drug production) |
| Equipment | Laboratory balances from all popular manufacturers |

| Ventilation technology | 1200 | 1500 |
|------------------------|--|---|
| Extract air monitoring | Air volume measurement using the measuring equipment of the fume hood function display | |
| Face velocity | 0,2 m/s (+/-10%) at 140 m ³ /h optional 0,3 m/s | 0,2 m/s (+/-10%) at 180 m ³ /h optional 0,3 m/s |
| Noise level | 54 dB | |
| Standard filtering | HEPA / ULPA filter H14 according to DIN EN 1822 | |
| Voltage | 230 V | |
| Power consumption | 250 W | |


| Options | 1200 | 1500 |
|---------------------------------------|--------|------|
| Ventilation connection Ø [mm] | DN 125 | |
| Air exchange rate [m ³ /h] | 170 | 210 |
| Pressure loss [PA] | 10 | 15 |

Intended use

- Can be used where required, with connections for the services supply, e. g. service wings
- Unrestricted view into the hood from all sides
- Service outlets in the internal workspace
- Control units located horizontally on the service rail of the support unit

Design


Front view



- 1 Trolley
- 2 Worktop with surrounding increased edge
- 3 Viewing window and baffle (safety glass)
- 4 Gas outlet
- 5 Water outlet with sink and waste water lifting unit
- 6 Openings for pipes and cables

Mobile fume hoods AeroEm

Rear view



- 1 Valve for water outlet
- 2 Handle with sash and horizontal sash
- 3 FAZ control panel incl. switch for interior sockets
- 4 Valve for gas outlet

Technical data

| Dimensions | |
|----------------------|------|
| Width [mm] | 1050 |
| Depth [mm] | 835 |
| Height [mm] | 1975 |
| Working height [mm] | 900 |
| Height, castors [mm] | 120 |

| Weight | |
|-------------|-----|
| Weight [kg] | 180 |

| Design characteristics | |
|-----------------------------|--|
| Sash | Two-piece, moves up and down with 2 horizontal sashes each |
| Glass pane in the side wall | All 4 sides of the fume hood |
| Lighting | Dazzle-free, can be switched from the outside |
| Roller shutter guiding | For pipes and cables on the left and right side of the fume hood |

| Electrics | |
|--|--|
| Electrical supply | 2 sockets in the internal workspace, can be switched individually from the outside |
| Total power of sockets [W] | 1000 |
| Connection voltage [V AC] | 230 |
| Voltage of waste water lifting unit [V] | 230 |
| Power of lighting [W] | 52 |
| Length, electrical connection cable [mm] | 2500 |

| Sanitary technology | |
|------------------------|--|
| Water connection | Optional |
| Waste water connection | Waste water quick release outlet as an option |
| Gas connection | Optional |
| Water fitting (tap) | Cold water WPC or WNC (EN) as an option, with drip cup, can be operated from the outside |
| Gas outlet | Optional |

| Ventilation technology | |
|---|----------------------------|
| Minimum air exchange rate [m ³ /h] ¹⁾ | 300 |
| Air-supply assistance fan | Can be switched on the FAZ |
| Function display | FAZ |
| 2 extract air spigots Ø [mm] | 90 |
| Length of extract air duct [mm] | 2500 |

¹⁾ All air volume specifications refer to an opening height of the sash window of 500 mm (test opening in acc. with EN 14175-3) and the maximum tracer gas values recommended by German Standard (BG RC1).

The indicated minimum air exchange rates were determined under specified test conditions in acc. with EN 14175-3.

These minimum air exchange rates must be adapted when dimensioning the ventilation system.

If on-site extract air monitoring systems or airflow dampers are used, the required air volumes may be different. The operating limitations must be agreed upon with Waldner.

| Material/surface | |
|------------------|--|
| Worktop | Stoneware-composite worktop with raised Polypropylene edge |

Mobile fume hoods


MobilAir

Intended use

- Can be used where required (only in air-circulating mode)
- Control units located externally
- Not suitable for openly breaking down chemicals
- Not suitable as a replacement for bench-mounted fume hoods in acc. with EN 14175


Design

Air-circulating mode



- 1 Sash with handle
- 2 FAZ control panel
- 3 Removable fascia panel
- 4 Filter housing with ventilator in air-circulating mode
- 5 Rear panel with air guiding profile
- 6 Material lock
- 7 Sockets

Extract air operation



- 1 Sash with handle
- 2 FAZ control panel
- 3 Removable fascia panel
- 4 Extract air spigot
- 5 Rear panel with air guiding profile
- 6 Material lock
- 7 Sockets

Technical data

| Dimensions | |
|--|-----------|
| Width [mm] | 900 |
| Depth [mm] | 600 |
| Height with sash closed/open [mm] | 1215/1620 |
| Access width [mm] | 730 |
| Clear width, internal workspace [mm] | 850 |
| Effective depth [mm] | 503 |
| Clear internal height up to lamp [mm] | 846 |
| Clear internal height up to ceiling [mm] | 935 |

| Weight | |
|---|------------|
| MobilAir for extract air operation [kg] | Approx. 70 |
| MobilAir for air-circulating mode incl. filter [kg] | Approx. 82 |

Mobile fume hoods

MobilAir

| Design characteristics | |
|------------------------|--|
| Air-circulating mode | With ventilator and filter |
| Extract air operation | Extract air spigot connected to on-site extract air system |
| Lighting | Dazzle-free, can be switched from the outside |
| Sash | Moves vertically |
| Material lock | Left and right solid side panel |

| Electrics | |
|----------------------------|--------------------|
| Electrical supply | 2 external sockets |
| Total power of sockets [W] | 1000 |
| Connection voltage [V AC] | 230 |
| Lighting [W] | 13 |
| Ventilator power [W] | 115 |

| Ventilation technology | |
|--|------------------|
| Minimum air exchange rate [m ³ /h] | 300 |
| Function display | FAZ as an option |
| Connection height [mm] Extract air spigot Ø 125 mm | 1137 |

| Material | |
|-------------------------|-----------|
| Side panel design, sash | Plexiglas |

| Filter type „A“ no.5, gas filter | |
|---|--|
| Dimensions [mm] | 610 x 305 x 150 (+ 8 mm seal) |
| Pressure loss [Pa] at 300 m ³ /h | 130 |
| Design characteristics | Gas filter cell with 5 layers of activated carbon mat, type „A“; MDF frame; with white-painted grid on both sides, with grip and type label on the 610-mm-side, PU seal on the dust-laden air side |
| Use | Separable substances: organic gases and vapours (e.g. solvents, petrol fumes, toluol, benzol, kerosine, odours, hydrocarbons with a molar mass (g/mol) of 30 and higher in cold form, non-boiling (VOC, high-boiling substances) |


| Filter type „BEP“, gas and particle filter | |
|---|--|
| Dimensions [mm] | 610 x 305 x 150 (+ 8 mm seal) |
| Pressure loss [Pa] at 300 m ³ /h | 240 |
| Design characteristics | Combination filter Hepa H13 with activated carbon mat and particle filter, type „BEP“; MDF frame, with white-painted grid on both sides, with grip and type label on the 610-mm-side, PU seal on the dust-laden air side |
| Use | Separable substances: inorganic gases and vapours (e.g. chlorine, hydrosulphides, sulphur dioxide, hydrogen chlorides, cold and heated). Molecules and particle separation 99.95 % MPPS |

| Filter type „P“, particle filter cell | |
|---|---|
| Dimensions [mm] | 610 x 305 x 150 (+ 8 mm seal) |
| Pressure loss [Pa] at 300 m ³ /h | 150 |
| Design characteristics | Particle filter, type „P“, Hepa H13, Midilar MDSA; MDF frame, with white-painted grid on both sides, with grip and type label on the 610-mm-side, fold height 45 mm, PU seal on the dust-laden air side, filter medium flush on the dust-laden air side |
| Use | Separable substances: Particle separation 99.95 % MPPS, Hepa H13 |

Intended use

- Extraction of thermal loads, gases, fumes, aerosols or dust escaping from the internal workspace of the housing
- Reduced sound emission
- Not suitable for openly breaking down chemicals
- Not suitable as a replacement for bench-mounted fume hoods in acc. with EN 14175

Design



- 1 Extract air spigot
- 2 Horizontal sash
- 3 Ventilation slots

Technical data

| Dimensions | 1200 | 1500 | 1800 | 2100 |
|--------------------------------------|------|--------------------------|------|------|
| Width [mm] | 1200 | 1500 | 1800 | 2100 |
| Depth [mm] | | 565 715 750 900 | | |
| Height [mm] | | 1450 | | |
| Height incl. extract air spigot [mm] | | 1550 | | |
| Height incl. extract manifold [mm] | | 1750 | | |

Housings

Permanent enclosure

| Design characteristics | 1200 | 1500 | 1800 | 2100 |
|---|---|---------------------|------|------|
| Construction | Shorter rear panel for using the services if combined with service spines | | | |
| Sash | 2 horizontal sashes | 3 horizontal sashes | | |
| Extract air operation | Connected to on-site extract air system Extract manifold as an option | | | |
| Material lock | Optional | | | |
| Lighting | Optional | | | |
| Shelf board, inside | Optional | | | |
| Ventilation technology | | | | |
| Function display | FAZ as an option | | | |
| Connection height [mm] for extract air spigot Ø 125 mm | 1550 | | | |
| Material | | | | |
| Side panel design, sash | Safety glass | | | |


Local extraction devices

Underbench unit extraction system

Intended use

- For the extraction of safety cabinets (underbench units) used for the storage of hazardous materials
- For the extraction of underbench units in service spines and fume hoods

Design



1 Extract air spigot

Technical data

| Ventilation technology | |
|--|----|
| Air exchange rate [m ³ /h] | 30 |
| Ventilation connection (ascending duct) Ø [mm] | 90 |

| Material | |
|------------------|-----|
| Ventilating pipe | PPS |


Local extraction devices

Extraction arm

Intended use

- For the extraction of a specific area
- For fixing to service wings, service spines or the wall

Design



Technical data

| Dimensions | 50 | 75 |
|----------------------------------|-----|----|
| Pipe system Ø [mm] ¹⁾ | 50 | 75 |
| Coupling hood Ø [mm] | 350 | |
| Extraction maximum [mm] | 50 | 75 |

¹⁾ Pipe system Ø 50 mm only for fastening to the service wing

| Ventilation technology | 50 | 75 |
|---|-----|-----|
| Minimum air exchange rate [m ³ /h] | 50 | 100 |
| Admission pressure [Pa] | 150 | |
| Admission pressure [Pa] with Waldner airflow damper | 200 | |

| Material | |
|----------------|--------------------|
| Pipe | Anodised aluminium |
| Hinged bracket | Polypropylene |
| Coupling hood | Polycarbonate |
| Suction tip | Anodised aluminium |

Local extraction devices


Snorkel hood

1

Intended use

- For the specific extraction of fumes
- Connection to extract air adapter in the service panel

Design



1 Snorkel hood

Technical data

| Dimensions | |
|---------------------------------------|------|
| Length of pipe system [mm] at Ø 40 mm | 1000 |
| Hood Ø [mm] | 120 |
| Suction tip [mm] | 50 |

| Ventilation technology | |
|----------------------------------|-----|
| Minimum air exchange rate [m³/h] | 5 |
| Admission pressure [Pa] | 200 |

| Material | |
|---------------|---------------|
| Pipe and hood | Polypropylene |


Local extraction devices

Extractor hood

Intended use

- For the extraction of a specific area
- For fixing to service spines and to the wall

Design



Technical data

| Dimensions | 1200 | 1500 |
|---------------------------|-----------|------|
| Width [mm] | 1200 | 1500 |
| Height x depth [mm] | 300 x 600 | |
| Extract air spigot Ø [mm] | 200 | |

| Ventilation technology | 1200 | 1500 |
|---|------|------|
| Minimum air exchange rate [m ³ /h] | 480 | 600 |
| Admission pressure [Pa] | 25 | 30 |
| Admission pressure [Pa] with Waldner airflow damper | 150 | |

| Material | |
|----------------|---------------|
| Extractor hood | Polypropylene |

Fume hoods and extraction devices



2 Service modules

Our **SCALA** range of laboratory furniture is defined by flexibility, mobility and ergonomic design to meet future requirements in the laboratory.

The supply of services plays a major part in a laboratory system.

Our service modules, including service spine, suspended service boom, service column and service wing, not only provide the services in the laboratory but also – more than ever – meet the ergonomic needs of the people working there. The service panels are inclined towards the user for easier accessibility of the fittings and control units.

Characterised by many useful details and a straightforward design, our service modules are fit to meet all requirements of laboratory design.

Our laboratory furniture system is made up of many fewer individual parts. Our service panels are fitted without joints, have even surfaces without edges, and the hidden accessory rail for supplementary functions is installed right where it is needed.

This simplifies cleaning and meets high hygienic requirements.



| | |
|------------------------------------|-----|
| Service duct element..... | 84 |
| Service spine..... | 86 |
| Service wing..... | 90 |
| Suspended service boom..... | 93 |
| Service column..... | 95 |
| Service distribution terminal..... | 98 |
| Service wall duct..... | 99 |
| Bench-mounted service duct..... | 100 |



2 Service modules

Space saving services installation

The services supply installations are integrated in the service duct to save space. The modular service panels are inclined towards the user for ergonomic access and handling. This, in turn, leads to a greater usable depth of the worktop.

The service spine

Our service spine gives the designer a basis for designing the laboratory environment and provides a large variety of options for different designs and rapid changes. The service spine is an autonomous unit and can be combined with freely selectable bench frames to form a wall bench or a double work bench.

The accessory rail for suitable accessories

The accessory rail below the service panel level is used for fitting useful accessories such as shelves, scaffold poles and towel rail. These „helpers“ can be moved over all grids and securely fastened.

Simple upgradability

The modular service panels without screws can be quickly replaced if necessary.

Supply pipes, for example for water or compressed air, can be rapidly expanded and fitted using a quick release coupling system without interrupting laboratory operation.

Configuration details of the service spine

The level above the service panels can be used as a shelf. The inserted glass shelves can easily be removed for cleaning. Above it, shelves can be fastened in the lateral pillars. The unit can always be expanded to the top by mounting overbench cabinets.



The service column

As a compact services supply, our service column enables the transparent design of the room.

The service column is equipped with removable panels and an accessory rail and can either be mounted directly to the laboratory ceiling or to the service ceiling.

The suspended service boom

The suspended service boom can be freely suspended from the laboratory ceiling which is useful for certain areas in the laboratory.

It is fitted with removable service panels and an accessory rail and can also be used for floor plans independent of the services. The suspended service boom can be height-adjusted when mounted to the ceiling. It is also possible to install the suspended service boom to the service ceiling.

The service wall duct

As an alternative to the service spine, the service wall duct can be mounted at different heights and directly to a wall, or connected after a service spine fitted against a wall. It is also equipped with panel technology and an accessory rail for variable configuration.



2 Service modules

The service wing system

Our service wing defines the term „freedom in the laboratory“ in a very special way: The new service wing is a major design element which integrates all services such as mechanical or electrical services, EDP, energy-saving lighting, extract air and the waste water disposal system, thus offering a high degree of flexibility.

The possibility of being able to plug in to the service wing for reliable supply and disposal connections practically everywhere means maximum freedom of movement and floor plan design in the laboratory.

The expansion stages of the service wing

The service wing has a modular design and offers four independent expansion stages for free combination. For every application. Using the removable service panels, fittings and connections can be placed as desired.

The accessory rail for useful accessories

The accessory rail accommodates useful accessories such as shelves, service distribution terminal and scaffold points. These can be moved over grids and securely fastened in every position.

Service wing for easy integration

Using the service wing simplifies the laboratory fitting out process and the coordination of different trades. One central feed point suffices.

Existing architectural features and building materials often require costly and time-consuming installations. Requiring minimum installation efforts, this is where the service wing is especially useful.



Energy-saving

The service wing is equipped with energy-saving lamps that illuminate the entire workspace and room and save up to 50 % power (with daylight-dependent control).

The service wing reaches the entire room

All areas of the laboratory are reached using T-elements and our wing segments of different lengths. For a large number of possible configurations. Thus it is possible to „dock“ anywhere, anytime.

All benches, racks, mobile sink units or mobile fume hoods can be used as required under the wing. For a flexible working environment.

Precise planning, pre-assembly in the factory

The service wing for your laboratory project is fully pre-assembled by our laboratory builders in accordance with the plans.

You save assembly time on-site and your service wing will be quickly installed and ready to use.

Uncomplicated modification and expansion

Since it is an individual system unit, the service wing can always be modified.


Expanding, upgrading and checking the system are possible with little effort.

Service duct element

Intended use


- Services supply at laboratory workstations
- Integration of all service outlets including sockets and multiple connectors for information technology
- Expansion and modification of the services supply through clip-in service panels
- Use in service spines, service wall ducts, suspended service booms, service columns and bench-mounted service ducts
- Tool-free installation of supplementary service duct add-on parts such as pegboard, monitor arm, pipette holder, paper towel dispenser, universal storage area, etc.

Design



- 1 Service panel with corner valves
- 2 Storage area
- 3 Service panel with sockets
- 4 Accessory rail for the tool-free installation of add-on parts

Service panel variants



- 1 Service panel with corner valves
- 2 Service panel with 8 sockets of the same type
- 3 Service panel with different types of sockets
- 4 Service panel with automatic circuit breakers

Technical data

| Dimensions | | | | | |
|--------------------------------------|-----------|-----|------|------|--------------------|
| Width [mm] | 600 | 900 | 1200 | 1500 | 1800 ¹⁾ |
| Depth [mm] without supporting system | 110 | | | | |
| Height [mm] | 252 | | | | |
| Front inclination [°] | 9 | | | | |
| Service panel, width x height [mm] | 300 x 200 | | | | |

¹⁾ The service duct can be extended as desired in grid lengths of 300 mm.

| Design characteristics | |
|--------------------------|--|
| Number of service panels | Depending on the width of the service duct Supply of electrics and information technology depending on the combination with other service modules |
| Service panel | Clip-in |
| Splash guard | Protection type IP 44 |

| Material | |
|--------------|---------------------------------|
| Storage area | Solid grade laminate shelf 5 mm |

| Electrics | |
|---|---------------------------|
| Electrical supply | Sockets in service panels |
| Fuse box | Optional |
| Max. number of sockets 230 V per service panel | 8 |
| Max. number of sockets 400 V per service panel | 2 |
| Max. number of automatic circuit breakers per service panel | 15 |

| Sanitary technology | |
|---|---|
| Sanitary supply | Service panel with take-off valves for vacuum, gases and/or waters Services supply depending on the combination with other service modules |
| Max. number of corner valves per service panel | 5 |
| Max. number of high purity gas valves per service panel | 3 to 5 depending on the type and function |


Service spine

Intended use

- For floor-mounted services supply of:
 - ▶ Wall benches
 - ▶ Double work benches
 - ▶ Laboratory equipment on mobile tables or underbench constructions
 - ▶ Floor-mounted laboratory equipment
- Design versions for genetical engineering areas
- Modular fastening of cell add-on parts to the multipurpose uprights, e.g. glass shelves and OSB board, overbench cabinets, scaffold points, etc.
- Tool-free installation of supplementary service duct add-on parts such as pegboard, monitor arm, pipette holder, paper towel dispenser, universal storage area, etc.

Design


Service spine for wall bench




- 1 Knee-hole cover panel
- 2 Sink module
- 3 Service panel with corner valves
- 4 Pillar for cell add-on parts
- 5 Service duct with service panel, glass shelf and accessory rail for add-on parts
- 6 Fascia panel of the service spine
- 7 Console
- 8 Multipurpose upright

Service spine


Service spine for wall bench with console and 2 glass shelves, working height 900 mm




Service spine for wall bench with C-frame, underbench units on castors and overbench cabinet, working height 750 mm



Service spine for wall bench with underbench units on plinth and media supply from above, working height 900 mm




Service spine for double bench with underbench units on plinth and media supply from above, working height




Service spine

Service spine for double bench with H-frame, underbench units on castors and overbench cabinet, working height 900 mm



Dimensional drawing

Service spine for wall bench/double bench



Technical data

| Dimensions | | | | | |
|---|-----------------------------|-----|------|------|------|
| Width [mm] | 600 | 900 | 1200 | 1500 | 1800 |
| Depth, service spine for wall bench [mm] (incl. wall bench) | 75 (750/900) | | | | |
| Depth, service spine for double bench [mm] (incl. double bench) | 92 (1500/1800) | | | | |
| Height [mm] | 1790 | | | | |
| Working height [mm] | 750 900 | | | | |
| Height, pillar extension [mm] for overbench cabinet, height 460 mm | 470 | | | | |
| Height, pillar extension [mm] for overbench cabinet, height 760 mm | 770 | | | | |
| Height, pillar extension [mm] up to ceiling height 3500 mm | Depending on ceiling height | | | | |
| Service panel, width x height [mm] | 300 x 200 | | | | |
| Glass shelf, width x depth [mm] | Width, service spine x 150 | | | | |
| Shelf of OSB board, width x depth [mm] | Width, service spine x 300 | | | | |

| Load bearing capacity | |
|--|----|
| Glass shelf [kg] | 20 |
| Shelf of OSB board [kg] | 30 |
| Max. load per scaffold point with scaffold rod length 300 mm [kg] | 5 |

| Design characteristics | |
|------------------------------------|--|
| Modular design | Wall bench can be equipped on one side, double bench can be equipped on two sides Multipurpose uprights can be extended with service duct, e.g. for overbench cabinets Worktop, cantilever and underbench unit can be replaced without dismantling the installations Grid-independent mounting of accessories |
| Scaffold points \varnothing [mm] | 12 to 13 |
| Number of service panels | Depending on the width of the service duct |

| Electrics | |
|-------------------|------------------------------|
| Electrical supply | Sockets in the service panel |
| Fuse box | Optional |


| Sanitary technology | |
|---------------------|---|
| Sanitary supply | Service panel with take-off valves for vacuum, gases and/or waters The supply pipes and cables are routed underneath the worktop or cantilever |

Service wing

Intended use


- Laboratory areas with technical devices for services
- Services supply and disposal via the ceiling for:
 - ▶ Laboratory benches and sinks below the service wing
 - ▶ Local extraction devices and AeroEm fume hood
 - ▶ Laboratory equipment on mobile tables or underbench constructions
 - ▶ Floor-mounted laboratory equipment
- Tool-free installation of supplementary service wing add-on parts

Design




Expansion stage 1

- Electrical trunking with service panels for the power supply




Expansion stage 2

- Electrical trunking with service panels for the power supply
- Wing edge designed as a lamp





Expansion stage 3

- Electrical trunking with service panels for the power supply
- Sanitary duct
- Ventilation duct
- Wing edge designed as a lamp

**Expansion stage 4**

- Electrical trunking with service panels for the power supply
- Sanitary duct
- Ventilation duct
- Wing edge designed as an accessory for the sanitary and ventilation routing

**Dimensional drawing****Service wing, expansion stage 3**

Service wing

Technical data

| Dimensions | | | | |
|---|------------|-----|------|------|
| Width [mm] | 600 | 900 | 1200 | 1500 |
| Depth [mm] with expansion stage 1 | 240 | | | |
| Depth [mm] with expansion stage 2 | 496 | | | |
| Depth [mm] with expansion stages 3 and 4 | 750 | | | |
| Height [mm] without dust cover for expansion stages 1 and 2 | 181 | | | |
| Height [mm] without dust cover for expansion stages 3 and 4 | 191 | | | |
| Outer dimensions of service panel [mm] | 300x220x29 | | | |

| Load bearing capacity | |
|-------------------------------|-----|
| Maximum permissible load [kg] | 120 |

| Design characteristics | |
|------------------------|---|
| Construction | Feeding, wing, T-element as an option Flexible bracing to prevent vibrations Can be equipped on both sides Dust protection through grid elements installed above |

| Electrics | |
|-------------------|---|
| Electrical supply | Electrical trunking with service panels for the power supply Connections for telephone, computer, monitor and loudspeaker as an option |
| Lighting | Lamps integrated in wing edges (direct and indirect lighting) as well as down light in the electrical trunking as an option |
| Fuse box | Optional |

| Sanitary technology | |
|---------------------|---|
| Sanitary supply | Service panels with take-off valves for vacuum, gases and/or waters Supply pipes and cables, ventilation duct guiding Local extraction system and/or extract air spigot for AeroEM as an option |


Suspended service boom

Intended use

- Services supply from the ceiling for:
 - ▶ Laboratory benches below the suspended service boom
 - ▶ Laboratory equipment on mobile tables or underbench constructions
 - ▶ Floor-mounted laboratory equipment
- Design versions for genetical engineering areas
- Modular fastening of boom add-on parts to the supporting construction, e.g. glass shelves and OSB board, scaffold points, etc.
- Tool-free installation of supplementary service duct add-on parts such as monitor arm, pipette holder, paper towel dispenser, universal storage area, etc.

Design

- 1 Service panel with sockets
- 2 Service panel with corner valves
- 3 Shelf of OSB board
- 4 Media supply duct
- 5 Ceiling stator
- 6 Functional element
- 7 Glass shelf
- 8 Pillar
- 9 Storage area (service duct)
- 10 Service duct element



Suspended service boom

Technical data

| Dimensions | | | | | |
|--|-------------------------------------|-----|------|------|------|
| Width [mm] | 600 | 900 | 1200 | 1500 | 1800 |
| Depth [mm] without pillars | 350 | | | | |
| Depth [mm] with pillars | 471 | | | | |
| Recommended min. height [mm] bottom edge of suspended service boom to upper edge of finished floor | 1750 | | | | |
| Height, supporting construction (max. up to ceiling height 4000 mm) | Depending on ceiling height | | | | |
| Service panel, width x height [mm] | 300 x 200 | | | | |
| Glass shelf, width x depth [mm] | Width, suspended service boom x 150 | | | | |
| Shelf of OSB board, width x depth [mm] | Width, suspended service boom x 300 | | | | |

| Load bearing capacity | |
|---|-----|
| Maximum permissible load [kg] | 120 |
| Additional max. load bearing capacity, suspended service boom [kg] per grid | 30 |
| Glass shelf [kg] | 20 |
| Shelf of OSB board [kg] | 30 |
| Max. load per scaffold point with scaffold rod length 300 mm [kg] | 5 |

| Design characteristics | |
|-------------------------------------|--|
| Construction | Functional elements to take up service ducts fastened to the ceiling and connected |
| Number of service panels (per side) | Depending on the width of the service duct |
| Scaffold points \varnothing [mm] | 12 to 13 |

| Material | |
|-----------------------------|---------------------------|
| Storage area (service duct) | Solid grade laminate 5 mm |

| Electrics | |
|-------------------|------------------------------|
| Electrical supply | Sockets in the service panel |
| Fuse box | Optional |


| Sanitary technology | |
|---------------------|---|
| Sanitary supply | Service panel with take-off valves for vacuum, gases and/or waters Supply pipes and cables in supply duct from above |

Intended use

- Services supply from the ceiling for:
 - ▶ Laboratory benches below the suspended service column
 - ▶ Laboratory equipment on mobile tables or underbench constructions
 - ▶ Floor-mounted laboratory equipment
- Version with one or two sides

Design

- A: SimplyMount, for fastening to bare ceilings or walls
 B: SimplyMove, for fastening to ceiling grids / beams
 C: SimplyMobile, for fastening to moveable base frames




- 1 Electrical cover for connecting terminals
- 2 Separate panels for corner valves, sockets and IT
- 3 Bare ceiling suspension with steel plate and profiled suspension
- 4 Service column on one/two sides of ceiling grid/beam with built-in anti-slip device and sliding element
- 5 Service ceiling, beam design by Waldner
- 6 Fastening for wall installation
- 7 Moveable service column

Service column


**SimplyMove fastening models
on 100x100 mm ceiling grid profile
on the room / double bench axis**

- 1 Ceiling grid / beam profile
- 2 Service column equipped on both sides on ceiling grid profile
- 3 Single-sided service column back to back, moveable against each other on the ceiling grid profile.



**SimplyMove fastening models
on the 50x100 mm ceiling grid profile
in front of partition wall / wall line:**

- 1 Service column equipped on one side on ceiling grid profile in the wall area
- 2 Service column equipped on one side fixed in front of the ceiling grid profile (fixed to the profile in the room with the partition wall model)



Technical data

| Dimensions | |
|---|---|
| Width [mm] | 130 / 190 |
| Depth [mm] single-sided | 100 |
| Depth [mm] two sides | 160 |
| Height [mm] | 910 + 1520 + 1820 |
| Height, supporting construction [mm] (max. up to ceiling height 4500 mm) | Adapted to ceiling height |
| Service panel, width x height [mm] | 100 / 160 x variable depending on equipment |


| Design characteristics | |
|-------------------------------|--|
| Construction | Vertical services carried, can be equipped on one or both sides Directly mounted to the aluminium support system, fitted to the bare slab or wall |
| Equipment | Corner valves, high purity gas valves, socket and IT sockets are possible depending on the requirement |

Service distribution terminal

Intended use

- Services supply for clamping to a laboratory workstation
- The station is supplied through a service module which is fastened to the ceiling, such as suspended service boom, service column, service wing, service ceiling or a floor-mounted service spine

Design



- 1 Connection to several service modules
- 2 Clamping system
- 3 Service distribution terminal with sanitary installation
- 4 Plinth element
- 5 Service distribution terminal with 4 sockets

Technical data

| Dimensions | |
|------------------------------------|-----------|
| Width [mm] | 158 |
| Depth [mm] | 118 |
| Height [mm] | 205 |
| Height incl. plinth element [mm] | 310 |
| Service panel, width x height [mm] | 150 x 200 |
| Clamping area [mm] | 10 – 100 |

| Design characteristics | |
|------------------------|--|
| Construction | Clamping system for worktop or other frames Services supply via service modules or service spines mounted to the ceiling Tension relief for pipes and cables between the service distribution terminal and service module unit through service beam and straps Cables and hoses are connected to the service module by means of plug-in couplings |


| Electrics | |
|-------------------|---|
| Electrical supply | Max. of 4 sockets 230 V per service panel |

| Sanitary technology | |
|---|---|
| Sanitary supply | Various take-off valves for vacuum, gases or compressed air |
| Max. number of corner valves per service panel | 2 |
| Max. number of high purity gas valves per service panel | 1 or 2 (depending on the type and function) |

Intended use

- Wall-mounted services supply for:
 - ▶ Laboratory benches under or in front of the service wall duct
 - ▶ Laboratory equipment on mobile tables or underbench constructions
 - ▶ Floor-mounted laboratory equipment
- Design versions for genetical engineering areas
- Tool-free installation of supplementary service duct add-on parts such as monitor arm, pipette holder, paper towel dispenser, universal storage area, etc.

Design



- 1 Service panel with sockets
- 2 Service panel with corner valves
- 3 Storage area, service wall duct
- 4 Media supply from above
- 5 Fastening profile/media supply
- 6 Service duct element

Technical data

| Dimensions | | | | | |
|------------------------------------|-----------|-----|------|------|--------------------|
| Width [mm] | 600 | 900 | 1200 | 1500 | 1800 ¹⁾ |
| Depth [mm] | 184 | | | | |
| Height [mm] | 252 | | | | |
| Service panel, width x height [mm] | 300 x 200 | | | | |

¹⁾ The service wall duct can be extended as desired in grid lengths of 300 mm.

| Load bearing capacity | |
|-----------------------|-----------------------|
| Storage area [kg] | 20 per installed grid |

| Design characteristics | |
|--------------------------|---|
| Construction | Service duct for wall mounting incl. solution for inside corner |
| Number of service panels | Depending on the width of the service duct |

| Material | |
|--------------|---------------------------------|
| Storage area | Solid grade laminate shelf 5 mm |

| Electrics | |
|-------------------|----------------------------|
| Electrical supply | Service panel with sockets |
| Fuse box | Optional |

| Sanitary technology | |
|---------------------|---|
| Sanitary supply | Service panel with take-off valves for vacuum, gases and/or waters Supply pipes in the fastening profile |


Bench-mounted service duct

Intended use

- Services supply of double work benches
- Design versions for genetical engineering areas
- Modular fastening of cell add-on parts to the multipurpose uprights, e.g. glass shelves and OSB board, overbench cabinets, scaffold points, etc.
- Tool-free installation of supplementary service duct add-on parts such as pegboard, monitor arm, pipette holder, paper towel dispenser, universal storage area, etc.
- Not suitable for double benches where separate work surfaces are required (see also BGI/GUV-I 850-0)

Design


Bench-mounted service duct with C-frame and suspended underbench unit




- 1 Worktop
- 2 Drip cup with water outlet
- 3 Service panel with sockets
- 4 Service panel with corner valves
- 5 Storage area, service duct
- 6 Service duct element
- 7 Multipurpose upright
- 8 Media supply duct
- 9 Suspended underbench unit

Bench-mounted service duct

Bench-mounted service duct with overbench cabinets, H-frame and underbench units on castors



Bench-mounted service duct with overbench cabinets, pillar extension and underbench units on plinth



Bench-mounted service duct

Technical data

| Dimensions | | | | | |
|---|---------------------------------|-----|------|------|------|
| Width [mm] | 600 | 900 | 1200 | 1500 | 1800 |
| Depth [mm] | 310 | | | | |
| Height [mm] | 1602 | | | | |
| Height, opening at working height 900 mm [mm] | 450 | | | | |
| Height, pillar extension [mm] (for overbench cabinet, height 460 mm) | 462 | | | | |
| Height, pillar extension [mm] (for overbench cabinet, height 760 mm) | 762 | | | | |
| Height, pillar extension [mm] (up to ceiling height 3500 mm) | Depending on ceiling height | | | | |
| Service panel, width x height [mm] | 300 x 200 | | | | |
| Glass shelf, width x depth [mm] | Width, bench-mounted unit x 150 | | | | |
| Shelf of OSB board, width x depth [mm] | Width, bench-mounted unit x 300 | | | | |

| Load bearing capacity | |
|--|----|
| Glass shelf [kg] | 20 |
| Shelf of OSB board [kg] | 30 |
| Max. load per scaffold point with scaffold rod length 300 mm [kg] | 5 |

| Design characteristics | |
|------------------------------------|--|
| Construction | Double-sided service duct as bench-mounted unit with opening above the worktop |
| Number of service panels | Depending on duct width |
| Scaffold points \varnothing [mm] | 12 to 13 |

| Material | |
|----------------------------|---------------------------------|
| Storage area, service duct | Solid grade laminate shelf 5 mm |

| Electrics | |
|-------------------|------------------------------|
| Electrical supply | Sockets in the service panel |
| Fuse box | Optional |

| Sanitary technology | |
|---------------------|--|
| Sanitary supply | Service panel with take-off valves for vacuum, gases and/or waters Supply pipes in the bench-mounted unit |



3 Laboratory benches and sinks

Laboratory benches are crucial in our **SCALA** laboratory furniture system..

The consequent separation of services supply and furniture creates flexibility in the laboratory.

All variants of our benches can be selected with various worktop materials for a large number of application possibilities everywhere in the laboratory.

High stability, straightforward design and perfect appearance characterise our laboratory benches.

Access to water must meet various requirements in the laboratory.

Large sink modules, integrated sinks, drip cups and sink modules in service modules or fume hoods are integrated in the laboratory as required by the specific situation.

Wherever mobility is required, our mobile units are used: under the service wing, for the suspended service boom, the service columns and the service ceiling – for fast moving in the laboratory.



| | | | |
|---|------------|---|------------|
| Laboratory benches | 110 | Special tables | 127 |
| Combinations of worktop and bench frame materials..... | 110 | Add-on table for low level fume hoods..... | 127 |
| Worktop material..... | 112 | Balance table | 128 |
| Bench with H-frame..... | 117 | Rack | 129 |
| Bench with C-frame..... | 118 | Heavy duty rack | 130 |
| Bench with self-supporting underbench units..... | 119 | Heavy duty mobile table | 131 |
| Mobile tables | 120 | Height-adjustable table | 132 |
| Laboratory sinks | 121 | | |
| Laboratory sink..... | 121 | | |
| End sink..... | 122 | | |
| Drip cup on service spine..... | 123 | | |
| Drip cup in worktop..... | 124 | | |
| Mobile sink..... | 125 | | |
| AquaEl..... | 126 | | |



3 Laboratory benches and sinks

Our benches offer a large number of possible applications.

Our new bench frames are made of precision rectangular tubes with reinforced cross-section. The bench frames can carry a load of 200 kg without any problem. Optimally protected against external effects through the entirely homogenous powder coating, our bench frames have a flawless appearance.

The same applies to the surfaces of our work-tops. You can choose from our wide range of materials according to your requirements.

Bench frames for different needs

With their constructional designs, C-frame and H-frame bench frames form the basis for our work benches depending on the requirement and application.

Different standard widths available

In order to be able to divide the workplaces in your laboratory to suit your needs, we offer a large number of frame widths.

Improved level compensation

Our new height-adjustable feet for H-frames offer up to 40 mm regulating distance. Easy access and adjustment - for steady positioning.

Easy cleaning

The height adjustment for the C-frame holds it approx. 30 mm above the floor. This makes cleaning the floor extremely easy.



H-frame

provides a high level of stability for add-on tables, mobile tables and analysis tables for working sitting or standing.

Underbench cabinets can be mobile or suspended and moved independent of modular size. Sitting niches are therefore possible anywhere.

C-frames

are extremely steady and can be loaded with 200 kg. They provide users with a large amount of knee and legroom with mobile and suspended underbench units.

Suspended underbench units that can be moved

Our new profile enables underbench units suspended in cantilever and C-frames to be moved across frames.

Movable knee-hole cover panels

For benches without underbench units we use movable and height-adjustable knee-hole cover panels. In this way, installations routed below the rear side of the table can be hidden.

Other useful helpers

Independent table frames can also be combined to create new situations – just as you need them. Our height-adjustable bench can be adjusted from 700 to 950 mm.

Our multi-talent: the rack

The rack is perfect for fitting items of equipment, AquaEl and others. The robust shelves are height-adjustable and the castors enable the fast changing of location.



3 Laboratory benches and sinks

There are no limits to the use of sinks and drip cups in the laboratory. Coordinated with our **SCALA** laboratory furniture system and manufactured from tried-and-tested materials, our sink units can ideally be integrated precisely where they are needed. Materials such as stoneware, poly-propylene, stainless steel and epoxy are extremely durable.

Sink modules and drip cups

Sink modules made of stoneware or polypropylene are integrated into the service spine above the bench. Drip cups are fitted directly in the worktop. They are made of stoneware, polypropylene, epoxy resin or stainless steel.

**Laboratory sink**

Sinks are permanently installed components of laboratory furniture and placed against the service spine or a wall. Sinks can be combined with various types of worktop materials in many versions.






Mobile sink and AquaEI

The mobile sink with castors supplements the variable laboratory below the service wing and service ceiling. The mobile sink is connected directly to the service wing or service ceiling system via flexible pipes. AquaEI is a ready to plug in compact system for the easy supply and disposal of water in service modules. A lifting unit disposes of the waste water through the respective system.

Laboratory benches

Combinations of materials/bench frames

Combinations of worktop and bench frame materials

| Material, worktop | Coated particle-board | Solid grade laminate | Solid grade laminate with EBC surface | Polypropylene |
|--|-----------------------|----------------------|---------------------------------------|-----------------|
|  H-frame | x | x | x | x |
|  C-frame | x | x | x | x |
|  Mobile table frame | x | x | x | x |
|  H-frame for low level fume hoods | - | - | - | x ¹⁾ |
|  Balance table | x | - | - | - |
|  Rack | x ²⁾ | - | - | - |






¹⁾ Material with surrounding increased edge

²⁾ Shelves white

Laboratory benches

Combinations of materials/bench frames


Combinations of worktop and bench frame materials


| Material, worktop | Epoxy | Stainless steel | Stoneware | Composite worktop | Glass |
|--|-----------------|-----------------|-----------------|-------------------|-------|
|  H-frame | x | x | x | x | x |
|  C-frame | x | x | x | x | x |
|  Mobile table frame | x | x | x | x | x |
|  H-frame for low level fume hoods | x ¹⁾ | x ¹⁾ | x ¹⁾ | - | - |
|  Balance table | - | - | - | - | - |
|  Rack | - | - | - | - | - |

¹⁾ Material with surrounding increased edge


Laboratory benches


Worktop material

| Melamine resin facing/postforming | |
|---|--|
| Critical substances | Acids in concentrations > 10 % |
| Damaging substances | Concentrated hydrochloric acids Nitric acid Heated sulphuric acid |
| Advantage | Flat |
| Limitations | Joints sensitive to moisture Medium chemical resistance |
| Use | Mobile table, add-on table, window benches Instrument benches and laboratory benches in the dry area Cannot be used in the moist or wet area |
| Weight [kg/m ²] | 19.6 |
| Overall thickness [mm] | 30 |
|  | Light grey NCS S 3005 R80B |

| Polypropylene | |
|---|---|
| Critical substances | Hydrocarbons Citric acid Oxalic acid Carbon tetrachloride Diesel oil |
| Damaging substances | Ozone Concentrated nitric acid Chloroform Petrol Benzol |
| Advantage | No joints Flat Light High chemical resistance to acids and many solvents Easy to dispose of Less breakage of glass |
| Limitations | Soft surface sensitive to scratches Sensitive to heat |
| Use | Areas with high resistance to chemicals Working with hydrofluoric acid Radio-isotope area Areas in which the lack of joints is important |
| Weight [kg/m ²] | 20.3 |
| Overall thickness [mm] | 30 |
| Increased edge [mm] | 7 |
|  | Light grey NCS S 3005 R80B |


Laboratory benches Worktop material


| Solid grade laminate | |
|---|---|
| Critical substances | Acids in concentrations > 10 % |
| Damaging substances | Concentrated hydrochloric acids Nitric acid Heated sulphuric acid |
| Advantage | Moisture-resistant Flat Easy to dispose of |
| Limitations | Reduced coating thickness |
| Use | Wet rooms Physical laboratories Tables with average load |
| Weight [kg/m ²] | 26.4 |
| Overall thickness [mm] | 19 |
| Increased edge (as an option) [mm] | 7 |
|  | Light grey NCS S 3005 R80B |

| Solid grade laminate with EBC surface | |
|---|---|
| Critical substances | Acids in concentrations > 10 % |
| Damaging substances | Concentrated hydrochloric acids Nitric acid Heated sulphuric acid |
| Advantage | Antibacterial Hardened to withstand electron beams High chemical resistance Moisture-resistant Flat Easy to dispose of |
| Limitations | Reduced coating thickness |
| Use | Chemical, microbiological, genetical-engineering laboratories |
| Weight [kg/m ²] | 26.4 |
| Overall thickness [mm] | 20 |
| Increased edge (as an option) [mm] | 7 |
|  | Grey Similar to NCS S 3000N |

Laboratory benches

Worktop material


| Epoxy | |
|---|---|
| Critical substances | Various solvents Diluted acids |
| Damaging substances | Hydrofluoric acid Concentrated warm mineral acids |
| Advantage | No joints Flat Solid panel High mechanical load capacity Easy to dispose of |
| Limitations | Surface sensitive to scratches Sensitive to concentrated acids |
| Use | Laboratory workstation of all type |
| Weight [kg/m ²] | 32 |
| Overall thickness [mm] | 19 |
| Increased edge [mm] | 7 |
|  | Platinum grey Similar to NCS S 4202-R |


| Stainless steel | |
|---|--|
| Critical substances | Cadmium Lactic acid Oxalic acid |
| Damaging substances | Compounds containing chlorine and bromine Formic acid Sulphuric acid |
| Advantage | No joints High resistance to solvents High temperature resistance |
| Limitations | Sensitive to halogens and their compounds |
| Use | For maximum loads in the area of decontamination and moisture resistance as well as solvent resistance Biology Microbiology Pharmacy Radio-isotope area Pathology |
| Weight [kg/m ²] | 27.5 |
| Overall thickness [mm] | 30 |
| Increased edge [mm] | 7 |
|  | |

Laboratory benches Worktop material

3


Laboratory benches and sinks

| Stoneware | |
|---|---|
| Critical substances | None |
| Damaging substances | Hydrofluoric acid |
| Advantage | Best chemical resistance Mechanically stable Easy to dispose of |
| Limitations | Evenness tolerances due to firing process Thermodynamic stress limited |
| Use | Areas subject to very high chemical and mechanical stress |
| Weight [kg/m ²] | 56 |
| Overall thickness [mm] | 26 |
| Increased edge [mm] | 7 |
|  | Light grey NCS S 3005 R80B |

| Composite worktop | |
|---|--|
| Critical substances | None |
| Damaging substances | Hydrofluoric acid |
| Advantage | Flat Lighter than stoneware Best chemical resistance Easy to dispose of |
| Limitations | Thermodynamic stress limited |
| Use | Areas with very high chemical stress |
| Weight [kg/m ²] | 40 |
| Overall thickness [mm] | 30 |
| Increased edge [mm] | 7 |
|  | White Similar to NCS S 0300-N |

Laboratory benches

Worktop material

| Glass | |
|---|---|
| Critical substances | None |
| Damaging substances | Hydrofluoric acid |
| Advantage | Flat High chemical resistance |
| Limitations | Sensitive to knocks at corners and edges |
| Use | Laboratory benches of all types subject to large amounts of chemicals |
| Weight [kg/m ²] | 38 |
| Overall thickness [mm] | 30 |
|  | Light green NCS S 2010 G10Y |


Laboratory benches

Bench with H-frame

Intended use

- Bench frame with worktop made of various materials as a work surface and storage area for laboratory work
- Supporting construction for analytical equipment and superstructures

Design



Technical data

| Dimensions | | | | | |
|---------------------|-------------------|-----|------|------|------|
| Width [mm] | 600 | 900 | 1200 | 1500 | 1800 |
| Depth [mm] | 600 750 900 | | | | |
| Working height [mm] | 750 900 | | | | |

| Load bearing capacity | |
|-----------------------|---|
| H-frame [kg] | 200 (for fixing to the wall or for fixing to a service spine) |

| Design characteristics | |
|------------------------|--|
| Construction | For suspended underbench units, cannot be moved for all kinds of frames For underbench units on castors |
| Height-adjustable feet | Individually adjustable |

| Material | |
|------------------------|------------------------------------|
| Bench frame | Steel profile 60/25/2 mm |
| Worktop | Depending on requirement |
| Height-adjustable feet | Plastic housing with steel spindle |


Laboratory benches

Bench with C-frame

Intended use

- Bench frame with worktop made of various materials as a work surface and storage area for laboratory work
- Supporting construction for analytical equipment and superstructures

Design



Technical data

| Dimensions | | | | | |
|---------------------|-----|-----|------|------|------|
| Width [mm] | 600 | 900 | 1200 | 1500 | 1800 |
| Depth [mm] | | | 600 | | |
| | | | 750 | | |
| | | | 900 | | |
| Working height [mm] | | | 750 | | |
| | | | 900 | | |

| Load bearing capacity | |
|-----------------------|-----|
| C-frame [kg] | 200 |

| Design characteristics | |
|------------------------|--|
| Construction | For suspended underbench units, can be moved for all kinds of frames For movable underbench units |
| Height-adjustable feet | Individually adjustable |

| Material | |
|------------------------|------------------------------------|
| Bench frame | Steel profile 70/25/3 mm |
| Worktop | Depending on requirement |
| Height-adjustable feet | Plastic housing with steel spindle |


Laboratory benches

Bench with self-supporting underbench units

Intended use

- Self-supporting underbench unit on plinth with worktop made of various materials as a work surface and storage area for laboratory work
- Supporting construction for analytical equipment and superstructures

Design



Technical data

| Dimensions | |
|-----------------------------|---------------------------|
| Overall width [mm] | Max. 2750 |
| Width, underbench unit [mm] | 450 600 900 1200 |
| Total depth [mm] | 750 900 |
| Working height [mm] | 750 900 |

| Material | |
|----------|------------------------------------|
| Worktop | Depending on width and requirement |

| Load bearing capacity | |
|---|-----|
| Bench with self-supporting underbench unit [kg] | 200 |


Mobile tables

Mobile table

Intended use

- Movable bench frame with worktop made of various materials as a work surface and storage area for laboratory work
- Movable supporting construction for analytical equipment and superstructures

Design



Technical data

| Dimensions | |
|---------------------|-------------------------|
| Width [mm] | 900 1200 1500 |
| Depth [mm] | 600 750 900 |
| Working height [mm] | 750 900 |

| Load bearing capacity | |
|----------------------------|-----|
| Mobile table [kg] | 150 |
| Per heavy load castor [kg] | 110 |

| Design characteristics | |
|------------------------|--|
| Heavy load castors | 4, of which 2 can be locked (castor and steering axle can be locked) |
| Shelf | Optional |
| Underbench unit | Optional |

| Material | |
|-------------|--------------------------|
| Bench frame | Steel profile 60/25/2 mm |
| Worktop | Depending on requirement |

Laboratory sinks


Laboratory sink

Intended use

- Water supply and disposal
- For cleaning operating equipment
- To take up large amounts of water
- Not suitable for the disposal of chemicals

Design

- 1 Underbench unit
2 Sink
3 Outlet



Technical data

| Material Worktop | Material Sinks | Sink dimensions Width x depth x height [mm] | Type of installation |
|---|-----------------|---|--|
| Stoneware | Stoneware | 380 x 380 x 250 530 x 380 x 250 680 x 380 x 250 | Sink installed flush with the worktop |
| Melamine resin facing, Solid grade laminate, Solid grade laminate with EBC surface | Polypropylene | 320 x 320 x 200 400 x 400 x 250 500 x 400 x 250 | Sink with surrounding increased edge installed in the worktop from above |
| Melamine resin facing, Solid grade laminate, Solid grade laminate with EBC surface | Stainless steel | 340 x 370 x 150 500 x 400 x 250 | Sink with surrounding increased edge installed in the worktop from above |
| Polypropylene | Polypropylene | 385 x 385 x 250 485 x 385 x 250 | Sink attached to the worktop from the bottom and welded |
| Stainless steel | Stainless steel | 400 x 400 x 250 500 x 400 x 250 | Sink welded in flush with the worktop |
| Composite worktop | Stoneware | 380 x 380 x 250 530 x 380 x 250 | Sink installed flush with the worktop |
| Epoxy | Epoxy | 406 x 305 x 203 406 x 406 x 190 457 x 380 x 279 | Sink installed flush with the worktop |

| Dimensions | | | | | |
|---------------------|---|-----|------|------|------|
| Width [mm] | 600 | 900 | 1200 | 1500 | 1800 |
| Depth [mm] | 600 ¹⁾ 675 ¹⁾ 705 ¹⁾ 750 825 855 900 | | | | |
| Working height [mm] | 900 | | | | |

¹⁾ Positioning of the outlets on the side of the sink, if required

| Sanitary technology | |
|------------------------|---|
| Water connection | Permanent connection |
| Waste water connection | Permanent connection with siphon |
| Water fitting (tap) | Bench-mounted service outlet as an option |
| Eye shower | Optional |

Laboratory sinks


End sink

Intended use

- Water supply and disposal
- For cleaning operating equipment
- To take up large amounts of water
- Not suitable for the disposal of chemicals

Design

- 1 Underbench unit
2 Sink
3 Outlet



Technical data

| Material Worktop | Material Sinks | Sink dimensions Width x depth x height [mm] | Type of installation |
|---|-----------------|---|--|
| Stoneware | Stoneware | 380 x 380 x 250 530 x 380 x 250 680 x 380 x 250 | Sink installed flush with the worktop |
| Melamine resin facing, Solid grade laminate, Solid grade laminate with EBC surface | Polypropylene | 320 x 320 x 200 400 x 400 x 250 500 x 400 x 250 | Sink with surrounding increased edge installed in the worktop from above |
| Melamine resin facing, Solid grade laminate, Solid grade laminate with EBC surface | Stainless steel | 340 x 370 x 150 500 x 400 x 250 | Sink with surrounding increased edge installed in the worktop from above |
| Polypropylene | Polypropylene | 385 x 385 x 250 485 x 385 x 250 | Sink attached to the worktop from the bottom and welded |
| Stainless steel | Stainless steel | 400 x 400 x 250 500 x 400 x 250 | Sink welded in flush with the worktop |
| Composite worktop | Stoneware | 380 x 380 x 250 530 x 380 x 250 | Sink installed flush with the worktop |
| Epoxy | Epoxy | 406 x 305 x 203 406 x 406 x 190 457 x 380 x 279 | Sink installed flush with the worktop |

| Dimensions | | |
|---------------------|------|------|
| Width [mm] | 1500 | 1800 |
| Depth [mm] | 740 | |
| Working height [mm] | 900 | |

| Sanitary technology | |
|------------------------|---|
| Water connection | Permanent connection |
| Waste water connection | Permanent connection with siphon |
| Water fitting (tap) | Bench-mounted service outlet as an option |
| Eye shower | Optional |


Laboratory sinks

Drip cup on service spine

Intended use

- Water supply and disposal
- For cleaning operating equipment
- Sink module underneath water fittings to take up small amounts of water
- Not suitable for the disposal of chemicals

Design



- 1 Service panel with corner valves
- 2 Sink module
- 3 Fascia panel for service spine

Technical data

| Dimensions | Polypropylene |
|--|-----------------------|
| Width [mm] | 294 |
| Depth [mm] | 132 |
| Height [mm] | 112 |
| Internal sink dimensions width x depth x height [mm] | Approx. 280 x 93 x 93 |

| Material | |
|-------------|---------------|
| Sink module | Polypropylene |

| Design characteristics | |
|------------------------|---|
| Construction | Attached to the fascia panel of the service spine |

| Sanitary technology | |
|------------------------|----------------------------------|
| Water connection | Permanent connection |
| Waste water connection | Permanent connection with siphon |
| Water fitting (tap) | Cell outlets as an option |


Laboratory sinks

Drip cup in worktop

Intended use

- Water supply and disposal
- For cleaning operating equipment
- Drip cup underneath water fittings to take up small amounts of water
- Not suitable for the disposal of chemicals

Design



- 1 Outlet
2 Drip cup
3 Worktop

Technical data

| Dimensions | |
|--|--|
| Width x depth [mm] | 295 x 145 |
| Height [mm] | Approx. 125 to 140 depending on material |
| Internal sink dimensions width x depth x height [mm] | Approx. 250 x 100 x 150 |

| Material, drip cup | Material, worktop |
|--------------------|---|
| Stoneware | Stoneware, composite worktop |
| Polypropylene | Polypropylene, melamine resin facing, solid grade laminate, solid grade laminate with EBC surface |
| Stainless steel | Stainless steel, melamine resin facing, solid grade laminate, solid grade laminate with EBC surface |
| Epoxy | Epoxy |

| Design characteristics | |
|------------------------|---|
| Construction | Installed in the worktop from the top or bottom |

| Sanitary technology | |
|------------------------|--|
| Water connection | Permanent connection |
| Waste water connection | Permanent connection with siphon |
| Water fitting (tap) | Bench-mounted service outlets as an option |

Laboratory sinks


Mobile sink

Intended use

- Mobile water and gas supply and disposal
- For cleaning operating equipment at any location
- Not suitable for the disposal of chemicals

Design

- 1 Connecting pipes
- 2 Fitting with two cold water outlet points
- 3 Worktop
- 4 Sink
- 5 Underbench unit on castors



Technical data

| Dimensions | |
|---|-----------------|
| Width [mm] | 605 |
| Depth [mm] | 600 |
| Height without outlet [mm] | 900 |
| Sink dimensions width x depth x height [mm] | 320 x 320 x 200 |
| Height, [mm] castors | 110 |
| Length, supply and drain pipes [mm] | 2500 |
| Length, connecting pipes [mm] | 2500 |

| Material | |
|----------|---|
| Worktop | Particle-board with melamine resin facing |
| Sink | Polypropylene |

| Load bearing capacity | |
|-----------------------|-----|
| Mobile sink [kg] | 150 |

| Design characteristics | |
|------------------------|---|
| Construction | <p>Mounted on underbench unit on castors with hinged door</p> <p>Sink installed in the worktop from above</p> <p>Pipes and cables routed out at the rear of the underbench unit</p> <p>Waste water lifting unit in the underbench unit</p> <p>Water supply is switched off in the case of a power failure</p> |

| Electrics | |
|------------------|-----|
| Power supply [V] | 230 |

| Sanitary technology | |
|------------------------|--|
| Water connection | Flexible with plug connector |
| Waste water connection | Flexible with plug connector |
| Gas connection | Flexible with plug connector as an option |
| Water fitting (tap) | Standard outlet |
| Gas outlet | Standard outlet combined with water fitting as an option |
| Mixer tap | Additional flexible water connection as an option |

Special tables


AquaEI

Intended use

- Mobile water and gas supply and disposal
- For cleaning operating equipment at the workplace at any mobile or stationary laboratory workstation
- Not suitable for the disposal of chemicals

Design

- 1 Connecting pipes
- 2 Outlet with water outlet point
- 3 Housing with pump
- 4 Sink



Technical data

| Dimensions | |
|--|-----------------|
| Width x depth x height (without outlet) [mm] | 317 x 585 x 268 |
| Sink, width x depth x height [mm] | 260 x 275 x 105 |
| Length, supply and drain pipes [mm] | 1500 |
| Length, connecting pipes [mm] | 1500 |

| Weight | |
|----------------------------|----|
| Weight without outlet [kg] | 14 |

| Material | |
|----------|---------------|
| Material | GFK varnished |

| Design characteristics | |
|------------------------|---|
| Construction | Compact system with flexible pipes and cables ready for connection Waste water lifting unit integrated in the housing Water supply is switched off in the case of a power failure |

| Electrics | |
|------------------|-----|
| Power supply [V] | 230 |

| Sanitary technology | |
|------------------------|--|
| Water connection | Flexible with plug connector |
| Waste water connection | Flexible with plug connector |
| Gas connection | Flexible with plug connector as an option |
| Water fitting (tap) | Standard outlet |
| Gas outlet | Standard outlet combined with water fitting as an option |
| Mixer tap | Additional flexible water connection as an option |


Special tables

Add-on table for low level fume hoods

Intended use

- For adding to low level fume hoods
- Bench frame with worktop made of various materials as a work surface and storage area for laboratory work
- Supporting construction for analytical equipment and superstructures

Design



Technical data

| Dimensions | | | | | |
|---------------------|-----|------|------|------|------|
| Width [mm] | 900 | 1200 | 1500 | 1800 | 2100 |
| Depth [mm] | 575 | | | | |
| Working height [mm] | 500 | | | | |

| Material | |
|------------------------|------------------------------------|
| Bench frame | Steel profile 60/25/2 mm |
| Worktop | Depending on requirement |
| Height-adjustable feet | Plastic housing with steel spindle |

| Load bearing capacity | |
|-----------------------|-----|
| H-frame [kg] | 200 |

| Design characteristics | |
|------------------------|----------------------------|
| Worktop | Surrounding increased edge |
| Height-adjustable feet | Individually adjustable |


Special tables

Balance table

Intended use

- For setting up analytical balances and other sensitive measuring equipment
- Bench frame with worktop and specially mounted, vibration-free plate

Design



- 1 Balance plate made of fine concrete
2 Table cover

Technical data

| Dimensions | |
|----------------------------------|------------|
| Width [mm] | 900 |
| Depth [mm] | 750 900 |
| Working height [mm] | 750 900 |
| Width x depth [mm] balance plate | 400 x 450 |

| Material | |
|-------------------------|---|
| Supporting construction | Steel profile |
| Worktop | Melamine Solid grade laminate Solid grade laminate with EBC surface |
| Balance plate | Fine concrete |


| Weight | |
|--------------------|-----|
| Total weight [kg] | 120 |
| Balance plate [kg] | 65 |

| Design characteristics | |
|------------------------|---|
| Construction | Specially mounted, heavy balance plate made of fine concrete Supporting construction with balance plate, vibration-decoupled |

Intended use

- Mobile flexible storage area
- Can be used with the 600 mm deep shelf as a mobile workplace for desk work
- Not suitable for storing hazardous substances

Design



- 1 Steel support frame with grid
- 2 Shelf, depth 450 mm
- 3 Shelf, depth 600 mm
- 4 Heavy load castors with brakes

Technical data

| Dimensions | |
|------------------------------------|------------|
| Width [mm] with shelf | 900 |
| Depth [mm] with shelf depth 450 mm | 600 |
| Height [mm] | 1790 |
| Depth, shelf [mm] | 450 600 |

| Material | |
|-------------------------|--------------------|
| Supporting construction | Steel profile |
| Shelf 22 mm | Shelf of OSB board |

| Load bearing capacity | |
|-----------------------|-----|
| Total [kg] | 150 |
| Shelf [kg] | 20 |

| Design characteristics | |
|--------------------------------|--|
| Heavy load castors | 4, of which 2 can be locked (castor and steering axle can be locked) |
| Shelf | Can be adjusted without tools with a pitch of 150 mm |
| Integrated distribution pillar | Optional |


Special tables

Heavy duty rack

Intended use

- Mobile multi-stage storage area
- With flexible work surfaces for free horizontal configuration
- Suitable for accommodating stackable and non-stackable measuring instruments / measuring instrument towers
- Suitable for heavy apparatus

Design



- 1 Lower shelf, fixed
- 2 Height-adjustable shelf, full width
- 3 Height-adjustable shelf, depth 590 mm
- 4 Steel support frame
- 5 Heavy load castors without brake
- 6 Heavy load castors with brake

Technical data

| Dimensions | | |
|---------------------------|-------------------------|-------------------------|
| Width [mm] | 1200 | 1800 |
| Depth [mm] | 770 | |
| Height [mm] | 1790 | |
| Shelf, width x depth [mm] | 400 x 590 1070 x 750 | 400 x 590 1670 x 750 |

| Material | |
|-------------------------|--------------------------|
| Supporting construction | Steel profile 70 x 40 mm |
| Shelf | Shelf of OSB board |

| Load bearing capacity | |
|-----------------------------|----------|
| Total | 500 [kg] |
| Shelf 400 x 590 [mm] | 30 [kg] |
| Shelf 1070 x 750 [mm] | 70 [kg] |
| Shelf 1670 x 750 [mm] | 70 [kg] |
| Lower shelf 1070 x 590 [mm] | 150 [kg] |
| Lower shelf 1670 x 590 [mm] | 150 [kg] |

| Design characteristics | |
|------------------------|--|
| Heavy load castors | 4, of which 2 can be locked (castor and steering axle can be locked) |
| Shelf | Can be adjusted with a grid of 75 mm |

Special tables

Heavy duty mobile table


Intended use

- Mobile heavy duty table with worktop and integrated sliding handle. Total load bearing capacity: 500 kg

Design

- 1 Worktop
- 2 Integrated sliding handle
- 3 Height-adjustable intermediate shelf
- 4 Lower shelf, fixed

3+4 optionally available as accessories,
not included in the basic version



Technical data

| Dimensions | | |
|--|----------------|----------------|
| Width [mm] | 1200 | 1500 |
| Depth [mm] | 750 / 900 | |
| Working height [mm] | 900 | |
| Worktop width x depth [mm] | 1070 x 750/900 | 1370 x 750/900 |
| Intermediate side shelf width x depth [mm] | 690 x 400 | 690 x 400 |
| Bottom shelf width x depth [mm] | 1070 x 690/840 | 1370 x 690/840 |

| Material | |
|-----------------------------------|--|
| Supporting construction | Steel profile 70 x 40 mm |
| Worktop | 30 mm particle-board, coated / 19 mm solid board |
| Intermediate side shelf and shelf | 30 mm particle-board, coated |

| Load bearing capacity | |
|-------------------------|--|
| Total | 500 [kg] |
| Worktop | 500 [kg] only if no intermediate shelves are installed |
| Intermediate side shelf | 30 [kg] |
| Bottom shelf | 150 [kg] |

Attention: The maximum load bearing capacity of, in total, 500 kg for the worktop, shelf and intermediate side shelf must not be exceeded with the additional shelves.

| Design characteristics | |
|-------------------------|---|
| Heavy load casters | 4x lockable (castor and steering axle can be locked), load bearing capacity 300 kg / caster |
| Intermediate side shelf | Can be adjusted with a grid of 75 mm |

| Options/accessories | |
|-------------------------|--|
| Bottom shelf | A shelf is placed on the securely welded supports. Load bearing capacity: 150 kg |
| Intermediate side shelf | A shelf (depth 40 mm) can be hung in the grid to the left and/or right. Load bearing capacity: 30 kg |


Special tables

Height-adjustable table

Intended use

- Bench frame with worktop made of various materials as a height-adjustable work surface and storage area for laboratory work
- Supporting construction for analytical equipment and superstructures

Design



Technical data

| Dimensions | | |
|---------------------|-----------|------|
| Width [mm] | 1200 | 1500 |
| Depth [mm] | 750 | 900 |
| Working height [mm] | 700 – 950 | |

| Material | |
|------------------------|------------------------------------|
| Bench frame | Steel profile 60/25/2 mm |
| Worktop | Depending on requirement |
| Height-adjustable feet | Plastic housing with steel spindle |

| Load bearing capacity | |
|-----------------------|-----|
| Total [kg] | 200 |

| Design characteristics | |
|------------------------|--------------------------------------|
| Working height | Can be adjusted with a grid of 25 mm |
| Bench frame | H-frame |



4 Storage cupboards

Our **SCALA** laboratory furniture system provides a vast selection of storage variants for fast access and safe storage.

All storage cupboards can be variably equipped and provide optimum space utilisation in all areas of the laboratory.

Designed with a high quality appearance and manufactured to Waldner's high quality requirements.

The melamine resin-coated surfaces are easy to care for and robust against the influences in the laboratory. The front edges on the carcass and the shelves are equipped with impact-resistant 2 mm polypropylene edges. Furthermore, the foil-coated plinth of our furniture is made of waterproof glued blockboard.

For laboratories with increased demands on cleanability, mechanical stress, moisture resistance or extreme dryness, we build our storage cupboards from solid grade laminate material. SGL compact panels do not require edge protection against moisture and are also suitable for extreme climate fluctuations. Solid grade laminate is also difficult to set on fire and has the property of delaying the spread of flames.

Wherever storage cabinets made of steel are required, we also offer our cabinet portfolio made of powder-coated sheet steel. A SEFA-certified version is optionally available. These storage spaces are tested and certified according to SEFA 8M and are suitable for the highest loads and loads.

With furniture made of melamine, solid grade laminate or powder-coated sheet steel, we offer the best selection for every need.



| | | | |
|-----------------------------------|------------|--|------------|
| Underbench units | 138 | Special cabinets | 159 |
| Underbench unit on plinth..... | 138 | Laboratory cabinet for storing | |
| Underbench unit on castors..... | 140 | acids and alkalis..... | 159 |
| Suspended underbench unit..... | 142 | Underbench safety unit for storing | |
| Push-in underbench unit | | acids and alkalis..... | 161 |
| for fume hoods..... | 144 | FWF 90 underbench safety unit for fume | |
| Underbench unit for sinks..... | 146 | hoods for storing flammable liquids..... | 163 |
| Overbench cabinets | 149 | FWF 90 safety cabinet for storing | |
| Laboratory cabinets | 151 | flammable liquids..... | 165 |
| Laboratory cabinet..... | 151 | G 90 gas cylinder cabinet..... | 167 |
| Emergency cabinet..... | 155 | | |
| Top-mounted cabinets | 156 | | |
| Pull-out cabinets | 157 | | |



4 Storage cupboards

Large number of variants

For maximum flexibility in the laboratory, we offer a large variety of cabinet and underbench unit variants. Push-in underbench units, either movable or on plinth, easily fit under C- and H-frames, or under fume hoods with their own supporting structure.

Suspended underbench units are integrated directly under the worktop as movable variants in C-frames.

Design and function go together

Seamless handles, optionally made from aluminium die-cast or stainless steel, are resistant to chemicals and easy to clean. Unique design features can be created in the laboratory by the use of oak trim fronts. Our overbench cabinets are fastened to the service spine or wall without a visible gap.

More mobility in the laboratory

Equipped with four smooth running swivelling castors – two of which can be locked – our movable underbench units can be simply pushed into the support frame of add-on tables or laboratory benches. The castor height is also harmonised and flush with the plinth height of our fixed cabinets.

More safety details

Due to the self-locking protection and change-pull-out catch of the drawers, our movable underbench units will not tilt over. Our top-mounted cabinets are fitted with a rail for safely securing ladders.

**More usable storage space**

With a depth of 550 mm for the underbench units and 500 mm drawer depth, the storage space is used to full capacity. The best solution offered in the market. We have also expanded the usable storage space of corner cabinets by implementing new fittings.

Optimal positioning

Due to four height-adjustable feet, our laboratory cabinets and underbench units on plinth can be set up straight and steady.

Fully extensible drawers with hidden roller rails

The double-wall steel frame with hidden roller rails is more robust, protected against soiling and thus runs a lot easier than single wall frames with open roll rails. Our standard fully-extensible drawers ensure a clear overview of their contents. All drawer units feature soft-closing mechanisms as standard.

Full extension drawers are available for steel cupboards in various options adapted to the respective markets. If required also with SEFA 8M certification.

Safety for problematic substances

Our safety cabinets for gases, acids, alkalis and flammable liquids meet the highest requirements on material properties and function. Of course the cabinets comply with the current standards.


Underbench units

Underbench unit on plinth


Intended use

- For storing equipment and chemicals in acc. with EN 16121 + EN 16122
- For working heights of 750 mm and 900 mm
- Not suitable for storing flammable liquids, gas cylinders and self-igniting or self-decomposing substances
- Not suitable for storing acids and alkalis

Design



Variants



Underbench units

Underbench unit on plinth

Technical data

| Dimensions | | | | |
|----------------------|---|-----|-----|------|
| Width [mm] | 450 | 600 | 900 | 1200 |
| Depth [mm] | 550 | | | |
| Overall height [mm] | 720 870 | | | |
| Height, drawers [mm] | 150 200 400 Combination possibilities see variants | | | |
| Height, plinth [mm] | 110 | | | |

| Load bearing capacity | |
|-----------------------|----|
| Per shelf/drawer [kg] | 30 |

| Design characteristics | |
|---------------------------|--|
| Construction | For working height 750 and 900 mm Hinged doors with 270° hinges Drawers, fully extensible Open at the top, rear panel can be removed Shelf, height-adjustable Without doors as a rack 4 height-adjustable feet |
| Combination possibilities | See variants |
| Handle | Handle bar SCALA U handle, stainless steel |
| Full-height drawers | Optional |
| Soft stop for drawer | Standard |
| Extract air spigot | Optional |
| Closing | Optional |


Underbench units

Underbench unit on castors





Intended use

- For storing equipment and chemicals flexibly in acc. with EN 16121 + EN 16122
- For working heights of 750 mm and 900 mm
- Not suitable for storing flammable liquids, gas cylinders and self-igniting or self-decomposing substances
- Not suitable for storing acids and alkalis

Design



Variants



Underbench units

Underbench unit on castors

Technical data

| Dimensions | | | | | | | |
|----------------------|---|-----|-----|-----|-----|------|------|
| Width [mm] | 450 | 545 | 600 | 845 | 900 | 1145 | 1200 |
| Depth [mm] | 550 | | | | | | |
| Overall height [mm] | 640 790 | | | | | | |
| Height, drawers [mm] | 150 200 350 Combination possibilities see variants | | | | | | |
| Height, castors [mm] | 110 | | | | | | |

| Load bearing capacity | |
|-----------------------|----|
| Per shelf/drawer [kg] | 30 |
| Per castor [kg] | 70 |

| Design characteristics | |
|------------------------------------|---|
| Construction | For working height 750 and 900 mm Hinged doors with 270° hinges Drawers, fully extensible and with change-pull-out catch Shelf, height-adjustable Without doors as a rack Covered at the top, rear panel permanently connected with the carcass 4 swivelling castors, front castors can be locked |
| Combination possibilities | See variants |
| Handle | Handle bar <i>SCALA</i> U handle, stainless steel |
| Soft stop for drawer | Standard |
| Closing | Optional |
| Drawers with change-pull-out catch | Standard |


Underbench units

Suspended underbench unit




Intended use

- For storing equipment and chemicals flexibly in acc. with EN 16121 + EN 16122
- For working heights of 750 mm and 900 mm
- Not suitable for storing flammable liquids, gas cylinders and self-igniting or self-decomposing substances
- Not suitable for storing acids and alkalis

Design



Variants



Underbench units

Suspended underbench unit

Technical data

| Dimensions | | | | | | | |
|----------------------|---|-----|-----|-----|-----|------|------|
| Width [mm] | 450 | 545 | 600 | 845 | 900 | 1145 | 1200 |
| Depth [mm] | 500 (depth of frame 572) 550 (depth of frame 672) | | | | | | |
| Height [mm] | 380 530 | | | | | | |
| Height, drawers [mm] | 150 200 350 Combination possibilities see variants | | | | | | |

| Load bearing capacity | |
|-----------------------|----|
| Per shelf/drawer [kg] | 30 |

| Design characteristics | |
|---------------------------|---|
| Construction | For working height 750 and 900 mm 2 fittings for attaching to the profile rail of the bench frame Hinged doors with 270° hinges Drawers, fully extensible Covered at the top, rear panel permanently connected with the carcass Shelf, height-adjustable For C-frame bench frame: Can be moved to the sides until it protrudes over the bench grid Hinged doors with 1 shelf at a height of 530 mm At a height of 530 mm without doors as a rack with 1 shelf |
| Combination possibilities | See variants |
| Handle | Handle bar SCALA U handle, stainless steel |
| Soft stop for drawer | Standard |
| Closing | Optional |


Underbench units

Push-in underbench unit for fume hoods


Intended use

- For storing equipment and chemicals in acc. with EN 16121 + EN 16122
- For fume hoods with rear panel installation and for fume hoods with side installation on a steel support frame
- Not permitted for storing flammable liquids, gas cylinders and self-igniting or self-decomposing substances
- Not permitted for storing acids and alkalis

Design



Variants



Underbench units

Push-in underbench unit for fume hoods

Technical data

| Dimensions | | | | | | |
|--|--|-----|-----|-----|------|------|
| Width [mm] | 545 | 600 | 845 | 900 | 1145 | 1200 |
| Depth [mm] | 550 | | | | | |
| Overall height [mm], push-in underbench unit for bench-mounted fume hoods with rear panel installation | 640 | | | | | |
| Overall height [mm], push-in underbench unit for bench-mounted fume hoods with side installation | 716 | | | | | |
| Height, plinth [mm] | 110 | | | | | |
| Height, drawers [mm] | 150 220 Combination possibilities see variants | | | | | |

| Load bearing capacity | |
|-----------------------|----|
| Per shelf [kg] | 30 |

| Design characteristics | |
|---------------------------|--|
| Construction | Hinged doors with 270° hinges Drawers, fully extensible Closed at the top, rear panel can be removed Shelf, height-adjustable 4 height-adjustable feet |
| Combination possibilities | See variants |
| Handle | Handle bar <i>SCALA</i> U handle, stainless steel |
| Extract air spigot | Optional |
| Closing | Optional |

Underbench units


Underbench unit for sinks

Intended use


- As an underbench unit for sinks for storing equipment and chemicals in acc. with EN 16121 + EN 16122
- Not suitable for storing flammable liquids, gas cylinders and self-igniting or self-decomposing substances
- Not suitable for storing acids and alkalis

Design

Sink with underbench unit for service spines or wall benches




End sink for double benches




Underbench units Underbench unit for sinks

Variants

Sink with underbench unit for service spines or wall benches



End sink for double benches



Underbench units

Underbench unit for sinks

Technical data

| Dimensions | | | | | |
|---------------------|-------------------|-------------------|--------------------|--------------------|--------------------|
| Width [mm] | 600 ¹⁾ | 900 ¹⁾ | 1200 ¹⁾ | 1420 ²⁾ | 1720 ²⁾ |
| Depth [mm] | 550 | | | 700 | |
| Overall height [mm] | 870 | | | | |
| Height, plinth [mm] | 110 | | | | |

¹⁾ For sinks on service spines or wall benches

²⁾ For end sinks


| Load bearing capacity | |
|-----------------------|----|
| Per shelf/drawer [kg] | 30 |

| Design characteristics | |
|------------------------|---|
| Construction | Hinged doors with 270° hinges 4 height-adjustable feet Inclined swivel door with waste bin 30 l Waste bin 2 x 15 l with full-height drawer Waste bin 2 x 35 l with full-height drawer Hinged door(s), full-height drawer Combination possibilities see variants |
| Handle | Handle bar SCALA U handle, stainless steel |
| Closing | Optional |




Intended use

- For storing equipment and chemicals in acc. with EN 16121 + EN 16122
- Not permitted for storing flammable liquids, gas cylinders and self-igniting or self-decomposing substances
- Not permitted for storing acids and alkalis

Design




Variants



Overbench cabinets

Overbench cabinet



Technical data

| Dimensions | | | | |
|-------------|------------|-----|-----|------|
| Width [mm] | 450 | 600 | 900 | 1200 |
| Depth [mm] | 350 | | | |
| Height [mm] | 460 760 | | | |

| Load bearing capacity | |
|-----------------------------------|----|
| Per shelf [kg] | 30 |
| Load bearing capacity, total [kg] | 60 |

| Design characteristics | |
|---------------------------|--|
| Construction | Height-adjustable fitting for fastening to the wall or to the service spine For a width of 1200 mm with central panel Shelf, height-adjustable |
| Combination possibilities | See variants |
| Handle | U handle SCALA U handle, stainless steel Glass sliding door with affixed plastic handle |
| Closing | Optional |


Laboratory cabinets

Laboratory cabinet

Intended use

- For storing equipment and chemicals in acc. with EN 16121 + EN 16122
- Not permitted for storing flammable liquids, gas cylinders and self-igniting or self-decomposing substances
- Not permitted for storing acids and alkalis

Design




Laboratory cabinets

Laboratory cabinet

Variants





Laboratory cabinets

Laboratory cabinet

Technical data

| Dimensions | | | | |
|---------------------|-----|-----|------|------|
| Width [mm] | 450 | 600 | 900 | 1200 |
| Depth [mm] | | | 350 | 550 |
| Overall height [mm] | | | 2090 | |
| Height, plinth [mm] | | | 110 | |


| Load bearing capacity | |
|-----------------------|----|
| Per shelf [kg] | 30 |

| Design characteristics | |
|---------------------------|--|
| Construction | Hinged doors with 270° hinges Shelves, height-adjustable Drawers, fully extensible 4 height-adjustable feet |
| Combination possibilities | See variants Drawers only with a depth of 550 mm |
| Handle | U handle SCALA U handle, stainless steel |
| Shelves, extendable | Optional (with a cabinet depth of 550 mm) |
| Drawers | Optional (with a cabinet depth of 550 mm) |
| Soft stop for drawer | Standard |
| Extract air spigot | Optional |
| Closing | Optional |

Intended use

- For storing protection and rescue materials (fire extinguisher, first aid case, etc.)
- Not permitted for storing flammable liquids, gas cylinders and self-igniting or self-decomposing substances
- Not permitted for storing acids and alkalis

Design



Technical data

| Dimensions | |
|---------------------|------------|
| Width [mm] | 600 |
| Depth [mm] | 350 550 |
| Overall height [mm] | 2090 |
| Height, plinth [mm] | 110 |

| Design characteristics | |
|------------------------|--|
| Construction | Hinged door with 270° hinges 4 shelves, height-adjustable 4 height-adjustable feet |
| Equipment | First aid case Fire extinguisher, 5 kg Sand boxes Shovel Fire blankets |


Top-mounted cabinets

Top-mounted cabinet

Intended use

- For storing equipment and chemicals in acc. with EN 16121 + EN 16122
- Only suitable as a permanently installed top part on the following Waldner cabinets: Laboratory cabinet, pull-out cabinet, emergency cabinet and acids and alkalis cabinet
- Not permitted for storing flammable liquids, gas cylinders and self-igniting or self-decomposing substances
- Not permitted for storing acids and alkalis

Design



Technical data

| Dimensions | | | | |
|-------------|-----|-----|-----|------|
| Width [mm] | 450 | 600 | 900 | 1200 |
| Depth [mm] | | | 350 | 550 |
| Height [mm] | | | 610 | 760 |

| Load bearing capacity | |
|-----------------------|----|
| Per shelf [kg] | 30 |

| Design characteristics | |
|------------------------|---|
| Construction | With rail for securing a ladder For laboratory cabinets with or without extract air spigot 1 shelf, height-adjustable Hinged doors |
| Handle | U handle SCALA U handle, stainless steel |
| Hook ladder | Optional |
| Closing | Optional |


Pull-out cabinets

Pull-out cabinet

Intended use

- For storing liquid or solid substances in suitable containers in acc. with EN 16121 + EN 16122
- Not permitted for storing flammable liquids, gas cylinders and self-igniting or self-decomposing substances
- Not permitted for storing acids and alkalis

Design



- 1 Wire basket with tray
2 Pull-out

Pull-out cabinets

Pull-out cabinet

Technical data

| Dimensions | |
|-----------------------------------|----------------|
| Width [mm] | 600 900 |
| Depth [mm] | 550 |
| Overall height [mm] | 2090 |
| Height, plinth [mm] | 110 |
| Tray, width x depth x height [mm] | 240 x 425 x 40 |

| Load bearing capacity | |
|-----------------------|-----|
| Per drawer [kg] | 120 |
| Per tray [kg] | 10 |

| Design characteristics | |
|--------------------------|--|
| Construction | 5 wire baskets with trays for each drawer, height-adjustable Fastened to the wall 4 height-adjustable feet Drawer doors with drawers accessible from both sides |
| Handle | U handle SCALA U handle, stainless steel |
| Soft stop for drawers | Optional |
| Compartment partitioning | Optional |
| Extract air spigot | Optional |
| Closing | Optional |

| Material | |
|----------|--------------|
| Tray | Polyethylene |


Special cabinets

Laboratory cabinet for storing acids and alkalis

Intended use

- For storing limited amounts of flammable acids and alkalis
- Not suitable for storing flammable liquids, gas cylinders and self-igniting or self-decomposing substances

Design



Special cabinets

Laboratory cabinet for storing acids and alkalis

Technical data

| Dimensions | |
|---------------------|------|
| Width [mm] | 600 |
| Depth [mm] | 550 |
| Overall height [mm] | 2090 |
| Height, plinth [mm] | 110 |

| Load bearing capacity | |
|-----------------------------------|----|
| Per shelf, height-adjustable [kg] | 30 |
| Per pull-out shelf [kg] | 20 |

| Design characteristics | |
|------------------------|--|
| Construction | Connection to the permanently active ventilation system 4 shelves, fixed or pull-out 4 height-adjustable feet Separate compartments for acids and alkalis Trays made of polypropylene Coated fittings Hinged doors |
| Handle | U handle SCALA U handle, stainless steel |

| Ventilation data | |
|---|------|
| Air exchange rate [m ³ /h] | 100 |
| Ventilation connection Ø [mm] | 90 |
| Connection height extract air spigot [mm] | 2176 |


Special cabinets

Underbench safety unit for storing acids and alkalis

Intended use


- Push-in or self-supporting underbench unit for bench-mounted fume hoods for storing limited amounts of acids and alkalis
- Not suitable for storing flammable liquids, gas cylinders and self-igniting or self-decomposing substances

Design




Variants

Push-in underbench units for fume hoods



Self-supporting underbench units for benches



Special cabinets

Underbench safety unit for storing acids and alkalis

Technical data

| Dimensions | |
|---|-----------------|
| Width [mm] | 600 900 |
| Width for push-in underbench units [mm] | 545/600/845/900 |
| Depth [mm] | 550/650 |
| Overall height [mm] self-supporting underbench units for benches | 720/870 |
| Overall height [mm], push-in underbench units for bench-mounted fume hoods with rear panel installation | 640 |
| Overall height [mm], push-in underbench units for bench-mounted fume hoods with side installation | 716 |
| Height, plinth [mm] | 110 |

| Load bearing capacity | |
|-----------------------|----|
| Extendable shelf [kg] | 20 |

| Design characteristics | |
|------------------------|---|
| Construction | Connection to the permanently active ventilation system 4 height-adjustable feet Coated fittings 2 extendable shelves with trays Hinged doors Combination possibilities see variants |
| Handle | Handle bar SCALA U handle, stainless steel |

| Ventilation data | |
|--|----|
| Air exchange rate [m ³ /h] | 30 |
| Ventilation connection to the ascending duct Ø [mm] | 90 |


Special cabinets

FWF 90 underbench safety unit for fume hoods for storing flammable liquids


Intended use

- Push-in underbench unit for bench-mounted fume hoods for storing limited amounts of flammable liquids
- Not suitable for storing gas cylinders and self-igniting or self-decomposing substances
- Not suitable for storing acids and alkalis

Design



Variants



Special cabinets

FWF 90 underbench safety unit for fume hoods for storing flammable liquids

Technical data

| Dimensions | | | | |
|---------------------|-----|-----|------|------|
| Width [mm] | 600 | 890 | 1100 | 1400 |
| Depth [mm] | 600 | | | |
| Overall height [mm] | 635 | | | |
| Height, plinth [mm] | 35 | | | |
| Max. weight [kg] | 130 | 170 | 220 | 290 |

| Load bearing capacity | |
|-----------------------|----|
| Rigid shelf [kg] | 30 |
| Drawers [kg] | 25 |

| Design characteristics | |
|---------------------------|---|
| Construction | Connection to the permanently active ventilation system Connection to the earth wire with potential equalisation With closing Tray with perforated plate insert Self-closing through current-independent thermal activation in the case of fire Hinged doors Drawer |
| Combination possibilities | See variants |
| Handle | U handle, stainless steel |
| Additional tray pull-out | Optional for drawers |
| Regulations and standards | EN 14470-1 TRGS 510 |

| Ventilation data | |
|--|----|
| Air exchange rate [m³/h] | 30 |
| Ventilation connection to the ascending duct Ø [mm] | 90 |

| Material | |
|------------------------|---|
| Underbench unit | Powder-coated stainless steel on the outside Colour: Pure white RAL 9010 |
| Ventilation connection | PPS |


Special cabinets

FWF 90 safety cabinet for storing flammable liquids

Intended use


- For storing limited amounts of flammable liquids
- Not suitable for storing gas cylinders and self-igniting or self-decomposing substances
- Not suitable for storing acids and alkalis

Design



Storage cupboards

Variants



Special cabinets

FWF 90 safety cabinet for storing flammable liquids

Technical data

| Dimensions | | | |
|---------------------|--------------|-----|------|
| Width [mm] | 600 | 900 | 1200 |
| Depth [mm] | Approx. 600 | | |
| Overall height [mm] | Approx. 2000 | | |
| Height, plinth [mm] | Approx. 80 | | |
| Max. weight [kg] | 290 | 360 | 470 |

| Load bearing capacity | |
|-----------------------|----------------------|
| Basin bed [kg] | Depending on version |

| Design characteristics | |
|-----------------------------------|--|
| Construction | Connection to the permanently active ventilation system Connection to the earth wire with potential equalisation Self-closing through current-independent thermal activation in the case of fire 3 basin beds, height-adjustable Tray with perforated plate insert With closing 4 height-adjustable feet Hinged doors |
| Combination possibilities | See variants |
| Other versions and configurations | On request |
| Regulations and standards | EN 14470-1 TRGS 510 |

| Ventilation data | |
|---------------------------------------|----|
| Air exchange rate [m ³ /h] | 30 |
| Ventilation connection Ø [mm] | 75 |

| Material | |
|------------------------|---|
| Laboratory cabinet | Powder-coated stainless steel on the outside Colour: Pure white RAL 9010 |
| Ventilation connection | Galvanised steel |


Special cabinets

G 90 gas cylinder cabinet


Intended use

- For storing gas cylinders in buildings
- Not suitable for storing flammable liquids and self-igniting or self-decomposing substances
- Not suitable for storing acids and alkalis

Design



Variants



Special cabinets

G 90 gas cylinder cabinet

Technical data

| Dimensions | 600 | 900 | 1200 | 1400 |
|----------------------|--------------|-----|------|------|
| Width [mm] | 600 | 900 | 1200 | 1400 |
| Depth [mm] | Approx. 600 | | | |
| Overall height [mm] | Approx. 2000 | | | |
| Max. net weight [kg] | 390 | 530 | 660 | 740 |

| Design characteristics | 600 | 900 | 1200 | 1400 |
|---|--|-----|------|------|
| Construction | Connection to the permanently active ventilation system Mounting rail to take up gas reduction units Roll-in ramp for gas cylinders With closing 4 height-adjustable feet Feed-throughs for pipes and cables in the cabinet ceiling Hinged door(s) | | | |
| Max. number of 50 l gas cylinders for cabinet width | 1 | 3 | 4 | 4 |
| Other versions and configurations | On request | | | |
| Regulations and standards | EN 14470-2 | | | |

| Ventilation data | 600 | 900 | 1200 | 1400 |
|---|-----|-----|------|------|
| Air exchange rate [m ³ /h] for cabinet width | 60 | 90 | 120 | 140 |
| Ventilation connection Ø [mm] | 75 | | | |

| Material | |
|------------------------|---|
| Laboratory cabinet | Powder-coated stainless steel on the outside Colour: Pure white RAL 9010 |
| Ventilation connection | Galvanised steel |



5 Supply and disposal

For the disposal of liquid and solid substances, we offer our TÜV-certified systems for use in corresponding underbench units.

As a standard feature, our underbench units for waste disposal are equipped with safety trays to accommodate suitable containers. For more container replacement convenience.

Acids, alkalis and flammable liquids can be disposed of directly into the containers through screw-mounted safety funnels, or from the internal workspace through the funnels in the worktop.

Mechanical or electronic level indicators and suitable ventilation systems make these systems complete.

Our latest underbench units for the disposal of solid substances are supplied with two robust waste bins with a capacity of 35 l in a fully extensible drawer, or as a tilting door variant with a waste bin that holds 30 l.

Supply system for flammable liquids

For the cyclic and continuous supply with flammable liquids, suitable safety cabinets are used that are connected to a permanent exhaust air system.

Our cabinets are in accordance with the relevant standards and regulations.

With the safety pistol-grip nozzle with flexible stainless steel supply pipe, flammable liquids can be safely drawn.



| | |
|--|------------|
| Supply system for flammable liquids | 172 |
| Waste disposal system for acids and alkalis | 175 |
| Waste disposal system for flammable liquids | 178 |
| Waste disposal system for solid matter and domestic waste | 181 |
| Waste disposal system for radio-isotope residual material | 183 |


Supply system for flammable liquids

Intended use

- For safely storing and providing flammable liquids at the laboratory workstation in accordance with EN 14470-1 (type 90) and TRGS 510 (appendix L)
- For transferring flammable liquids from containers into small containers (max. 2 containers with 30 l each)
- Not permitted for supplying the following hazardous substances:
 - ▶ Acids and alkalis
 - ▶ Gas cylinders
 - ▶ Radioactive substances
 - ▶ Microorganisms

Design


Cyclic supply



- 1 Pistol-grip nozzle in the internal workspace
- 2 Outlet pipe
- 3 Top-mounted cabinet
- 4 3-way valve
- 5 Inert gas pipe
- 6 Pressure regulator
- 7 Container
- 8 Safety cabinet

Supply system for flammable liquids

Continuous supply with automatic container changeover



- 1 Pistol-grip nozzle in the internal workspace
- 2 Outlet pipe
- 3 Top-mounted cabinet
- 4 3-way valve
- 5 Inert gas pipe
- 6 Electric module of the monitoring system
- 7 Pressure regulator
- 8 Safety cabinet
- 9 Container 1
- 10 Container 2

Supply system for flammable liquids

Technical data

| Dimensions | |
|--------------------------------------|-------------|
| Width [mm] | Approx. 600 |
| Depth [mm] | Approx. 600 |
| Height [mm] with top-mounted cabinet | 2700 |
| Container 30 l, height [mm] | 440 |
| Container 30 l, Ø [mm] | 370 |

| Design characteristics | |
|--|---|
| Construction | Safety cupboard with: Connection to the ventilation system Connection to potential equalisation with earth wire Self-closing through current-independent thermal activation in the case of fire Shelves, height-adjustable Tray Hinged door |
| Number of containers 30 l | 1-2 |
| Cyclic supply | With different flammable liquids Separate pipes to 1-2 containers in the safety cabinet |
| Continuous supply | With automatic changeover to the second container Common pipe connected to no more than 2 containers in the safety cabinet Monitoring system: automatic changeover to the second container if container is empty |
| Pressure regulator, solvent tapping system | Defined pressure of 0.2 bar for transporting the flammable liquid Safety valve from 0.5 bar |
| Outlet, solvent tapping system | Solvent pistol flexibly mounted in the internal workspace Solvent pistol rigidly mounted in the internal workspace |

| Material | |
|--|--------------------------------|
| Safety cabinet | Stainless steel, powder-coated |
| Container | Stainless steel |
| Connection spigot, ventilation Ø 75 mm | Galvanised steel |

| Ventilation data | |
|---|----|
| Air exchange rate [m ³ /h] | 50 |
| Ventilation connection to the ascending duct Ø [mm] | 90 |


Waste disposal system for acids and alkalis

Intended use

- For safely storing the remnants of acids and alkalis at the laboratory workstation temporarily
- Not permitted for the disposal of the following hazardous substances:
 - ▶ Flammable liquids
 - ▶ Gas cylinders
 - ▶ Radioactive substances
 - ▶ Microorganisms


Design

Filling through funnel in the underbench unit



- 1 Funnel
- 2 Canisters
- 3 Tray
- 4 Underbench unit with full-height drawer

Filling through funnel in the internal workspace



- 1 Funnel on the worktop
- 2 Electric module with level indicator and control units
- 3 Canisters
- 4 Underbench unit with hinged door (without drawer)
- 5 Tray

Waste disposal system for acids and alkalis

Technical data

| Dimensions for underbench unit on plinth | |
|--|-----|
| Width [mm] | 600 |
| Depth [mm] | 550 |
| Height [mm] at working height 750 mm | 720 |
| Height [mm] at working height 900 mm | 870 |
| Max. height [mm] | 530 |
| Height, plinth [mm] | 110 |

| Dimensions for self-supporting/push-in underbench unit for bench-mounted fume hoods | |
|---|-----|
| Width [mm] | 600 |
| Depth [mm] | 550 |
| Height [mm] at working height 900 mm | 639 |
| Max. height [mm] | 425 |
| Height, plinth [mm] | 110 |

| Dimensions for self-supporting/push-in underbench unit for bench-mounted fume hoods with side installation | |
|--|-----|
| Width [mm] | 600 |
| Depth [mm] | 550 |
| Height [mm] at working height 900 mm | 716 |
| Max. height [mm] | 530 |
| Height, plinth [mm] | 110 |

| Dimensions, canister | |
|----------------------------------|---|
| 5 l width x depth x height [mm] | 160 x 185 x 230, connection thread S 55 |
| 10 l width x depth x height [mm] | 190 x 230 x 340, connection thread S 60 |
| 20 l width x depth x height [mm] | 260 x 285 x 390, connection thread S 60 |

| Design characteristics | |
|-----------------------------------|---|
| Construction | Extracted underbench unit with full-height drawer (max. 2 containers) or extracted underbench unit with hinged door and without drawer (max. 2 containers) Coated fittings Tray made of polypropylene |
| Funnel | Underbench unit with full-height drawer: Funnel, fastened to canister with screws Underbench unit with hinged door: Funnel on worktop with filling pipe between funnel and canister |
| Filling | Funnel fastened with screws on canister: optical check of the filling level when the canister is transparent Funnel on the worktop: Electronic level indicator, acoustic and visual indication when the maximum level is reached |
| Approval, canister 5l, 10 l, 20 l | UN 3H1/Y1,9 |
| Resistance | Based on consultation with Waldner |

Waste disposal system for acids and alkalis

| Funnel in the underbench unit | Canister 5 l | Canister 10 l | Canister 20 l | Canister 10 l and 20 l |
|---|--------------|---------------|---------------|------------------------|
| Underbench unit on plinth for service spine | – | 4 | 2 | 2 x 10 l and 1 x 20 l |
| Push-in underbench unit for service spine | – | 4 | – | – |
| Push-in underbench unit for bench-mounted fume hoods | – | 4 | – | – |
| Push-in underbench unit for bench-mounted fume hoods with side installation | – | 4 | 2 | 2 x 10 l and 1 x 20 l |

| Funnel in the internal workspace | Canister 5 l | Canister 10 l | Canister 20 l | Canister 10 l and 20 l |
|--|--------------|---------------|---------------|------------------------|
| Underbench unit on plinth for bench-mounted fume hoods | 2 | 2 | – | – |
| Underbench unit on plinth for bench-mounted fume hoods with side installation | 2 | 2 | 1 | 1 x 10 l and 1 x 20 l |
| Push-in underbench unit for bench-mounted fume hoods and fume hoods with side installation | 2 | 2 | – | – |

| Material | |
|--------------------------------|-------------------------------|
| Canisters | PP |
| Ventilation connection Ø 90 mm | PPS |
| Tray | PP |
| Components for installation | Electrically conductive PE-HD |

| Ventilation data | |
|---|----|
| Air exchange rate [m³/h] | 50 |
| Ventilation connection to the ascending duct Ø [mm] | 90 |

Waste disposal system for flammable liquids


Intended use

- For safely storing remnants of flammable liquids at the laboratory workstation temporarily in accordance with EN 14470-1 (type 90) and TRGS 510
- For waste disposal using screw-mounted funnels in the underbench safety unit or through funnels on the worktop in the internal workspace
- Disposal of HPLC equipment via universal connection
- Not permitted for the disposal of the following hazardous substances:
 - ▶ Acids and alkalis
 - ▶ Gas cylinders
 - ▶ Radioactive substances
 - ▶ Microorganisms

Design


Filling through funnel in the underbench unit

- 1 Safety cabinet with full-height drawer
- 2 Funnel
- 3 Earthing cable
- 4 Mechanical level indicator
- 5 Canisters
- 6 Tray



Filling through funnel in the internal workspace


- 1 Funnel on the worktop
- 2 Electric module with level indicator and control units
- 3 Canisters
- 4 Safety cabinet with hinged door
- 5 Tray



Waste disposal system for flammable liquids

Filled through mounting spigot
for capillaries in the internal workspace

- 1 Universal connection on worktop
- 2 Electrical panel with level indicator and control units
- 3 Canisters
- 4 Safety cabinet with hinged door
- 5 Tray



Waste disposal system for flammable liquids

Technical data

1. Filling via funnel in underbench unit
2. Filling through funnels in the internal workspace

| Dimensions | |
|---|-------------------|
| Underbench safety unit, width x depth [mm] | Approx. 595 x 600 |
| Underbench safety unit, overall height [mm] | Approx. 600 |
| Canister 5 l, width x depth x height [mm] | 160 x 185 x 230 |
| Canister 10 l, width x depth x height [mm] | 198 x 298 x 264 |

| Design characteristics | |
|-----------------------------------|---|
| Construction | With funnel in the underbench unit: Underbench safety unit with full-height drawer with max. 2 containers With funnel in the internal workspace: Underbench safety unit with hinged door with max. 2 containers Connection to the ventilation system Connection to potential equalisation with earth wire Funnel, grounded |
| Canisters | 2 canisters, 5 l (insulated) 2 containers 10 l, conductive |
| Funnel | Underbench safety unit with full-height drawer: Funnel, fastened to canister with screws Underbench safety unit with hinged door, transfer system: Funnel on the worktop is connected with the canister through one filling pipe per funnel |
| Approval, canister 5l, 10 l, 30 l | UN 3H1/Y1,6 |
| Filling, level indicator | Funnel in the underbench safety unit: mechanical level indicator integrated in 10 l canister Funnel in the internal workspace: Electric level indicator, acoustic and visual indication when the maximum level is reached Connection for liquid chromatographic instrument (HPLC) with spigot instead of funnels and electric level indicator, as an option Filling head is connected to extract air via gas suspension cord |
| Resistance | Based on consultation with Waldner |

| Material | |
|--------------------------------|--------------------------------|
| Underbench safety unit | Stainless steel, powder-coated |
| Canister 5 l | PP |
| Canister 10 l | Electrically conductive PE-HD |
| Ventilation connection Ø 90 mm | PPS |
| Components for installation | Electrically conductive PE-HD |
| Components for transfer system | Stainless steel |

| Ventilation data | |
|---|----|
| Air exchange rate [m ³ /h] | 50 |
| Ventilation connection to the ascending duct Ø [mm] | 90 |


Waste disposal system for solid matter and domestic waste

Intended use

- For the disposal of remnants of solid matter and garbage from laboratory work
- Not suitable for the permanent storage of solid matter and garbage
- Not permitted for the disposal of hazardous substances, especially:
 - ▶ Acids and alkalis
 - ▶ Flammable liquids
 - ▶ Gas cylinders
 - ▶ Radioactive substances
 - ▶ Microorganisms


Design

Waste bin with full-height drawer



- 1 Full-height drawer
- 2 Waste bin 2 x 15 l
- 3 Waste bin 2 x 35 l

Waste bin with tilting door



- 1 Tilting door
- 2 Waste bin 30 l
- 3 Waste bin 2 x 35 l

Waste disposal system for solid matter and domestic waste

Technical data

| Dimensions for underbench unit on plinth | | | | |
|--|-------------------------|---------------|-------------------------|---------------|
| Width x height [mm] | 450 x 870 | 600 x 870 | 450 x 720 | 600 x 720 |
| Depth [mm] | 550 | | | |
| Height, plinth [mm] | 110 | | | |
| Capacity with full-height drawer | 2 x 15 l or 2 x 35 l | 2 x 15 l – | 2 x 15 l or 2 x 35 l | 2 x 15 l – |
| Capacity with tilting door | 1 x 30 l | | | |

| Dimensions for underbench unit for sinks | | | |
|--|-----------|-----------|------------|
| Width x height [mm] | 600 x 870 | 900 x 870 | 1200 x 870 |
| Depth [mm] | 550 | | |
| Height, plinth [mm] | 110 | | |
| Capacity with full-height drawer | – | 2 x 15 l | |
| Capacity with tilting door | 1 x 30 l | 2 x 30 l | |

| Dimensions for self-supporting underbench unit for bench-mounted fume hoods | |
|---|-----------|
| Width x height [mm] | 600 x 820 |
| Depth [mm] | 550 |
| Height, plinth [mm] | 110 |
| Capacity with full-height drawer | 2 x 15 l |
| Capacity with tilting door | 1 x 30 l |

| Dimensions for push-in underbench unit for bench-mounted fume hoods | | |
|---|-----------|-----------|
| Width x height [mm] | 545 x 639 | 600 x 639 |
| Depth [mm] | 550 | |
| Height, plinth [mm] | 110 | |
| Capacity with full-height drawer | 2 x 15 l | |
| Capacity with tilting door | 1 x 30 l | |

| Design characteristics | |
|---------------------------------|--|
| Door | Full-height drawer Tilting door |
| Automatic foot-operated opening | Optionally for full-height drawers up to a width of 600 mm |
| Extract air spigot | Optional |

| Material | |
|------------------------|-----|
| Ventilation connection | PPS |

| Ventilation data | |
|---|----|
| Air exchange rate [m ³ /h] | 30 |
| Ventilation connection to the ascending duct Ø [mm] | 90 |


Waste disposal system for radio-isotope residual material

Intended use

- Waste canister at the workplace for the safe disposal of slightly radioactive material
- Not permitted for the disposal of the following hazardous substances:
 - ▶ Acids and alkalis
 - ▶ Flammable liquids
 - ▶ Gas cylinders
 - ▶ Microorganisms

Design

Filling through funnel in the underbench unit (funnel with mechanical level indicator)



- 1 Underbench unit on castors
- 2 10 l Canisters
- 3 Funnel

Technical data

| Dimensions of underbench units for radio-isotope residual material | | |
|--|---|-----|
| Width [mm] | 450 | 600 |
| Depth [mm] | 550 | |
| Overall height [mm] | 639 | |
| Height, castors [mm] | 110 | |
| Canister 10 l, width x depth x height [mm] | 190 x 230 x 340, connection thread S 60 | |
| Collapsible box, width x depth x height [mm] | 300 x 300 x 500 | |

| Design characteristics of underbench units for radio-isotope residual material | |
|--|--|
| Construction | Front side with lead shield on the inside With castors Max. 2 canisters of 10 l in tray made of polypropylene to take up slightly radioactive, liquid residual material Collapsible box to take up solid radio-isotope residual material as an option |



6 Education

HELLO, WE ARE WALDNER EDUCATION!

Waldner Education provides innovative all-in-one solutions for educational institutions around the globe.

Our product systems stand for high quality and ergonomic flexibility while promoting trend-setting didactics.

LEARNING SPACES FOR THE EDUCATION OF THE FUTURE

www.waldner-education.com





7 Services

We are the only manufacturer of laboratory equipment who offers you fume hoods and variable fume hood control all from one supplier. Make the most of our knowledge and expertise on fume hood and laboratory control.

We have installed a large number of projects of different sizes worldwide, all of which are operating to our customers' great satisfaction. This underlines the philosophy of our global technology company.


Furthermore, you as a customer will find it convenient and economical to have only one contact for all questions on the issue and also for maintenance.

Being a full-range supplier, we will plan and implement your project in no time – in the typical Waldner way. Being a market leader, we have the necessary capacity for your project – no matter how big. Please contact us. We will be glad to help you.



Services

7



7 Services

Large cost savings in every operating state

From an economic point of view, the laboratory furniture and the ventilation of the entire laboratory building are no longer separate entities today. Waldner's intelligent laboratory control significantly reduces the operating costs of the ventilation system and ensures maximum work safety.



Sophisticated technology for optimum operation

Our fume hoods are an important part of laboratory ventilation and can be integrated into the building ventilation concept in an ideal way. The measuring and control unit incorporated in our Airflow-Controller reliably detects the operating status of the fume hood at all times, and precisely and safely regulates the air flow rate within seconds for safe extract air operation.

If required, the user can increase or decrease the air exchange rate at the fume hood manually at any time.

Investing in our laboratory control will quickly pay for itself

A cost-benefit analysis clearly speaks for our laboratory control: Since the ventilation system is efficiently used while the energy supply is reduced, investing into this laboratory control system will pay off within one to two years. Considering continuously increasing energy prices, this is an important advantage.



Ventilation and control as an overall concept

Being a leading system partner, we will develop an overall concept for your laboratory – from the appropriate sizing of the central ventilation system and the ducts to the selection and use of the appropriate process measuring and control technology.



- A** Air flow damper for local AC4 Compact extraction
- B** Mechanical airflow damper
- C** Airflow damper extract air AC4 Compact
- D1** Airflow-Controller AC4 v Standard
- D2** Airflow-Controller AC4 v pipe controller
- E** Airflow damper Supply air AC4 Compact
- F** CANopen bus communication
- G** Airflow-Controller with activated master function for laboratory control
- H** The following types of communication with DDC/GLT are possible, among others: Analogue I/O, LON bus, Modbus, Profibus, BACnet, Ethernet
- I** Sash controller SC

Control and monitoring

Control

Control – Airflow-Controller (AC) for fume hoods according to DIN EN 14175-6

Airflow-Controller AC

The central control unit is a microprocessor-based electronic control unit and forms the heart of the Waldner control system.


Our Airflow-Controller calculates quickly and precisely the required air flow rate for safe operation from the sash position based on a defined process. The required air flow rate is then precisely adjusted for you within three seconds. The damper takes less than one second to make a 90° change of angle.

To achieve a precise result, the controller uses a characteristic curve that is based on the damper position and the effective pressure.

Apart from its pure extract air control system, the controller also monitors the optionally available Secuflow technology. If the air flow falls below the required extract air volume, the supportive flow technology automatically switches off and an acoustic and visual warning is issued.

A visual and acoustic alarm signal is also issued when the permissible sash opening area is exceeded in accordance with the specifications of EN 14175.

The control flap is used as standard with the extract air hood. Use the motorised damper as a pipe controller to conserve space if the room height is less than 3.30 m.



1 Display and operating device




2 Sensors for detector of sash position



3 Exhaust hood with actuator, measuring system and measurement acquisition



4 Central control unit AC



Control and monitoring Control

The fume hood and controller are an entity

The systems are precisely matched to each other, thus ensuring maximum reliability during laboratory operation.

The fume hood and variable air volume control are type-approved in accordance with EN 14175-6 as a complete safety system. Thus, the time-consuming and costly coordination of different trades becomes unnecessary and legal security and warranty are provided by one supplier, if need be.


Our measurement method and measuring system

Due to the variable measuring diaphragm coefficient and the special principle of operation of the measuring system, an airflow stroke of 1:15 can be realised. During night operation, the air volume at the fume hood can thus be reduced to 100 m³/h.

A measuring accuracy of +/- 5 % of the current actual value of the air exchange rate is also guaranteed. This is necessary to ensure that the directed airflow in the laboratory is maintained even if the air exchange rates are low.



EN 14175-6 type tested fume hood control in acc. with 5.4
Measuring in the outer measuring level




Schematic diagram of effective pressure development



Control panel AC

- Menu key
- Light On/Off
- Visual and acoustic alarm
- Flushing function
(increasing the air volume)
- Lowered operation
- Monitoring and control
on / off
- Service connection



Measuring system

Airflow-Controller

- 1 Control flap
- 2 Pressure measuring ducts
- 3 Venturi aperture
- 4 Exhaust hood
- 5 Electronics with pressure sensor and solenoid valve

Technical data

| Characteristics | |
|---|------------------------------|
| Air exchange rate range for diameter DN 250 | 100 - 1500 m ³ /h |
| Air exchange rate range for diameter DN 315 | 200 - 3000 m ³ /h |
| Measuring accuracy to the actual value | +/- 5 % |
| Nominal capacity | 35 VA |
| Motor run time for 0-90° | 1 second |
| Control time | less than 3 seconds |
| Permitted system pressure | 100 - 500 Pa |

| Inputs | |
|----------------|--|
| Voltage supply | 230 V |
| Digital input | 3 pieces (freely parametrisable) |
| Analogue input | 2 pieces (freely parametrisable) |
| Sash detector | 2 pieces (sash and horizontal sash detector) |
| Connection | RJ 45 |
| CAN bus | RJ 45 |

| Outputs | |
|--------------------------|----------------------------------|
| Digital output | 4 pieces (freely parametrisable) |
| Analogue output | 1 piece (freely parametrisable) |
| Control panel connection | RJ 12 |
| CAN bus | RJ 45 |
| Motor control | RJ 45 |

| Design | |
|-------------------------------|----------------------|
| Airflow damper and monitoring | Constant or variable |

Control and monitoring

Laboratory control

Master function for room control

Intelligent room control



Our master function forms a total extract air volume from the individual extract air volumes of all the extraction units in the laboratory. If the minimum air exchange is not achieved by the minimum air values of the fume hoods, then the module automatically raises the required minimum air exchange rate of the fume hoods or room extract air. If the minimum room air exchange is exceeded by opening a fume hood, then the increase of the minimum extract air of the fume hoods is reduced again to the minimum air of the fume hoods at most. If this is not sufficient to achieve the minimum room air exchange, the next step is to reduce the room extract air. Energy-efficient air through-ventilation can therefore always be guaranteed.

This module is also capable of providing temperature- and room pressure-dependent control.


In the event of the maximum simultaneous extract air volume in a room being exceeded, an alarm is emitted at the specific fume hoods contributing to the excess.

Various bus systems, including BACnet or Modbus, can make available data points, such as the set-points and actual values of the air flow damper,

motorised damper positions, error messages, operating statuses and fume hood sash positions. Remote diagnostics can also be implemented for even faster troubleshooting by means of a simple interface.



Two examples of laboratory control variants



Status window in the parametrisation software

Control and monitoring

Airflow damper for room supply air and extract air

AC4 Compact

Areas of application

- Room input air controller
- Room extract air controller
- Airflow measuring system/measuring panel (without control flap and actuator)
- Extractor hoods, fume hoods

AC4 Compact

AC4 Compact, the microprocessor-based electronic control unit, controls the air volume infinitely.

It quickly and precisely regulates the air flow rate according to the setpoint.

Performance criteria

- Control parameters are adaptively optimised on-line
- Standard tolerances are predictively corrected using a theoretical process model
- Control of the position of the motorised damper
- Floating time: 5 sec.
- Freely parametrisable on a PC basis
- Integrated pressure sensor 0-250 Pa (pressure-resistant up to 2500 Pa)
- Motorised damper housing: galvanised, stainless steel, PPs
- Motorised damper housing with and without insulation case

Connections

- 1 x analogue output
- 1 x analogue input
- 2 x digital input
- 1 x control unit input RJ12
- 2 x CAN interfaces
- 1 x motor output RJ45
- 1 x 24 V double terminal



AC4 Compact



Actuator



Galvanised controller housing with AC4 Compact and fast actuator



Galvanised controller housing with AC4 Compact and fast actuator - rectangular version

Control and monitoring

Airflow damper for room supply air and extract air

Technical data

Flow rate range for rectangular air flow damper, housing of galvanized steel

| Construction dimensions | | Installation length | Width without and with insulation case | Height without and with insulation case | Vmin | V 7 m/s | Vmax (10 m/s) |
|-------------------------|-------------|---------------------|--|---|---------------------|---------------------|---------------------|
| Width [mm] | Height [mm] | [mm] | [mm] | [mm] | [m ³ /h] | [m ³ /h] | [m ³ /h] |
| 200 | 100 | 135 | 282 | 182 | 72 | 504 | 720 |
| 300 | 100 | 135 | 382 | 182 | 108 | 756 | 1080 |
| 400 | 100 | 135 | 482 | 182 | 144 | 1008 | 1440 |
| 500 | 100 | 135 | 582 | 182 | 180 | 1260 | 1800 |
| 600 | 100 | 135 | 682 | 182 | 216 | 1512 | 2160 |
| 300 | 150 | 170 | 382 | 232 | 162 | 1134 | 1620 |
| 400 | 150 | 170 | 482 | 232 | 216 | 1512 | 2160 |
| 500 | 150 | 170 | 582 | 232 | 270 | 1890 | 2700 |
| 600 | 150 | 170 | 682 | 232 | 324 | 2268 | 3240 |
| 200 | 200 | 220 | 282 | 282 | 144 | 1008 | 1440 |
| 300 | 200 | 220 | 382 | 282 | 216 | 1512 | 2160 |
| 400 | 200 | 220 | 482 | 282 | 288 | 2016 | 2880 |
| 500 | 200 | 220 | 582 | 282 | 360 | 2520 | 3600 |
| 600 | 200 | 220 | 682 | 282 | 432 | 3024 | 4320 |
| 800 | 200 | 220 | 882 | 282 | 576 | 4032 | 5760 |
| 300 | 250 | 270 | 382 | 332 | 270 | 1890 | 2700 |
| 400 | 250 | 270 | 482 | 332 | 360 | 2520 | 3600 |
| 500 | 250 | 270 | 582 | 332 | 450 | 3150 | 4500 |
| 600 | 250 | 270 | 682 | 332 | 540 | 3780 | 5400 |
| 800 | 250 | 270 | 882 | 332 | 720 | 5040 | 7200 |
| 300 | 300 | 325 | 382 | 382 | 324 | 2268 | 3240 |
| 400 | 300 | 325 | 482 | 382 | 432 | 3024 | 4320 |
| 500 | 300 | 325 | 582 | 382 | 540 | 3780 | 5400 |
| 600 | 300 | 325 | 682 | 382 | 648 | 4536 | 6480 |
| 800 | 300 | 325 | 882 | 382 | 864 | 6048 | 8640 |
| 1000 | 300 | 325 | 1082 | 382 | 1080 | 7560 | 10800 |
| 400 | 400 | 430 | 482 | 482 | 576 | 4032 | 5760 |
| 500 | 400 | 430 | 582 | 482 | 720 | 5040 | 7200 |
| 600 | 400 | 430 | 682 | 482 | 864 | 6048 | 8640 |
| 800 | 400 | 430 | 882 | 482 | 1152 | 8064 | 11520 |
| 1000 | 400 | 430 | 1082 | 482 | 1440 | 10080 | 14400 |
| 1200 | 400 | 430 | 1282 | 482 | 1728 | 12096 | 17280 |

Control and monitoring

Airflow damper for room supply air and extract air

Technical data

Flow rate range for round air flow damper, housing of galvanized steel

| Nominal size | Installation length | External diameter without insulation case | External diameter with insulation case | Vmin | V 7 m/s | Vmax (10 m/s) |
|------------------|---------------------|---|--|---------------------|---------------------|---------------------|
| Nominal diam. DN | [mm] | [mm] | [mm] | [m ³ /h] | [m ³ /h] | [m ³ /h] |
| 100 | 195 | 99 | 199 | 27 | 190 | 272 |
| 125 | 195 | 124 | 224 | 43 | 300 | 428 |
| 160 | 215 | 159 | 259 | 71 | 494 | 706 |
| 200 | 215 | 199 | 299 | 111 | 776 | 1108 |
| 250 | 260 | 249 | 349 | 174 | 1217 | 1739 |
| 315 | 260 | 314 | 414 | 277 | 1939 | 2770 |
| 400 | 315 | 399 | 499 | 448 | 3135 | 4479 |

Flow rate range for round air flow damper, housing of stainless steel

| Nominal size | Installation length | External diameter without insulation case | External diameter with insulation case | Vmin | V 7 m/s | Vmax (10 m/s) |
|------------------|---------------------|---|--|---------------------|---------------------|---------------------|
| Nominal diam. DN | [mm] | [mm] | [mm] | [m ³ /h] | [m ³ /h] | [m ³ /h] |
| 100 | 195 | 99 | 199 | 27 | 190 | 272 |
| 125 | 225 | 124 | 224 | 43 | 300 | 428 |
| 160 | 260 | 159 | 259 | 71 | 494 | 706 |
| 200 | 300 | 199 | 299 | 111 | 776 | 1108 |
| 250 | 375 | 249 | 349 | 174 | 1217 | 1739 |
| 315 | 470 | 314 | 414 | 277 | 1939 | 2770 |
| 400 | 555 | 399 | 499 | 448 | 3135 | 4479 |
| 500 | 800 | 564 | 599 | 701 | 4908 | 7012 |
| 630 | 800 | 704 | 729 | 1115 | 7806 | 11151 |

Flow rate range for round air flow damper, housing of polypropylene

| Nominal size | Installation length | External diameter without insulation case | External diameter with insulation case | Vmin | V 7 m/s | Vmax (10 m/s) |
|------------------|---------------------|---|--|---------------------|---------------------|---------------------|
| Nominal diam. DN | [mm] | [mm] | [mm] | [m ³ /h] | [m ³ /h] | [m ³ /h] |
| 110 | 300 | 110 | 198 | 31 | 214 | 306 |
| 125 | 325 | 125 | 219 | 40 | 280 | 400 |
| 160 | 360 | 160 | 259 | 67 | 470 | 671 |
| 200 | 400 | 200 | 298 | 106 | 745 | 1064 |
| 250 | 475 | 250 | 348 | 367 | 1169 | 1670 |
| 315 | 570 | 315 | 414 | 263 | 1841 | 2630 |
| 400 | 655 | 400 | 499 | 426 | 2980 | 4257 |
| 500 | 850 | 500 | 599 | 662 | 4636 | 6623 |
| 630 | 1045 | 630 | 729 | 1052 | 7365 | 10521 |

Control and monitoring

Airflow damper for room supply air and extract air

Technical data

Flow rate range for rectangular air flow damper, housing of polypropylene

| Construction dimensions | | Installation length | Width without insulation case | Height without insulation case | Width with insulation case | Height with insulation case | Vmin | V 7 m/s | Vmax (10m/s) |
|-------------------------|------|---------------------|-------------------------------|--------------------------------|----------------------------|-----------------------------|------|---------|--------------|
| [mm] | [mm] | | | | | | | | |
| 140 | 140 | 530 | 200 | 200 | 240 | 240 | 135 | 473 | 676 |
| 200 | 140 | 530 | 260 | 200 | 300 | 240 | 194 | 680 | 972 |
| 250 | 140 | 530 | 310 | 200 | 350 | 240 | 244 | 853 | 1218 |
| 160 | 160 | 530 | 220 | 220 | 260 | 260 | 177 | 621 | 887 |
| 280 | 160 | 530 | 340 | 220 | 380 | 260 | 313 | 1096 | 1566 |
| 180 | 180 | 580 | 240 | 240 | 280 | 280 | 226 | 789 | 1128 |
| 315 | 180 | 580 | 375 | 240 | 415 | 280 | 398 | 1392 | 1988 |
| 200 | 200 | 580 | 260 | 260 | 300 | 300 | 279 | 978 | 1397 |
| 355 | 200 | 580 | 415 | 260 | 455 | 300 | 495 | 1734 | 2477 |
| 630 | 200 | 580 | 690 | 260 | 730 | 300 | 883 | 3092 | 4417 |
| 224 | 224 | 580 | 284 | 284 | 324 | 324 | 348 | 1220 | 1742 |
| 400 | 224 | 580 | 460 | 284 | 500 | 324 | 627 | 2195 | 3136 |
| 250 | 250 | 580 | 310 | 310 | 350 | 350 | 436 | 1525 | 2179 |
| 280 | 280 | 580 | 340 | 340 | 380 | 380 | 548 | 1920 | 2742 |
| 400 | 280 | 580 | 460 | 340 | 500 | 380 | 787 | 2754 | 3935 |
| 315 | 315 | 620 | 375 | 375 | 415 | 415 | 696 | 2437 | 3482 |
| 355 | 355 | 620 | 415 | 415 | 455 | 455 | 887 | 3105 | 4435 |
| 400 | 400 | 620 | 460 | 460 | 500 | 500 | 1129 | 3952 | 5645 |
| 500 | 400 | 620 | 580 | 480 | 600 | 500 | 1408 | 4927 | 7039 |
| 630 | 400 | 620 | 710 | 480 | 730 | 500 | 1770 | 6196 | 8851 |
| 800 | 400 | 620 | 880 | 480 | 900 | 500 | 2252 | 7883 | 11262 |

Control and monitoring

Monitoring


Control – Function display (FAZ) for fume hoods according to EN 14175-2

EN 14175-2 requires continuous monitoring of the ventilation functions of fume hoods to warn laboratory personnel by visual and acoustic signals in the event of a fault. The visual signals cannot be deleted.

Our function display takes over this task. It continuously measures the extract air volume flow and alerts users visually and acoustically if the level falls below the required minimum air exchange rate. It is also capable of monitoring our optional Secuflow technology. As with the AC, the control unit for the function display is ergonomically incorporated in the profile. Statuses are signalled here by means of different coloured lights and a display, and can be acknowledged if necessary.

Airflow measurement FAZ


The air flow rate present is determined by a differential pressure signal that is measured at the extract air hood. This measuring method is independent of room pressure fluctuations and the opening of the sash. In night mode, a second air volume can be monitored.



Differential pressure measurement FAZ

1 Exhaust hood, available in two designs:
Diameter of 250 mm and diameter of 315 mm
2 Electronics with pressure sensor

Diameter of measuring tube of 250 mm
for fume scrubber and filter fume hoods



Technical data

| Monitoring | Function display (FAZ) |
|-------------------|--|
| Power supply | 230 V |
| Outputs | Alarm output Operating message Light switch |
| Inputs | On Off Acoustic alarm acknowledgement Night operation |
| Diameter [mm] | 250, 315 |
| System connection | Analogue I/O, Modbus |

Sash controller SC

The SC supports the opening and closing of the sash by means of a motor, resulting in three benefits. All that is needed is for the sash to be gently touched and the opening and closing process is automatically continued/completed. It can also be operated by a foot switch, providing hands-free control.

The SC is coupled with a motion detector to increase the safety and energy efficiency of the fume hood and monitor the area in front of the fume hood. If no one is detected in this area, the SC automatically closes the sash after a pre-set time. Objects in the sash area are reliably detected by a light barrier, and the sash stops to protect the experiment set-up. The SC also automatically puts into practice the specifications of Technical Regulation TRGS 526. According to this Technical Regulation, fume hoods not currently in operation need to be closed.

In combination with an Airflow Controller, the SC can also be connected to the DDC/GLT.

Component parts:

- 1) Processor-controlled central control unit
- 2) Motor drive (closes and opens the sash)
- 3) The photo-electric barrier integrated in the sash frame serves to detect obstacles in the path of the sash when the sash is automatically closed
- 4) The motion detector stops the sash when working in front of the sash



Technical data SC

| Closing device | Sash controller SC |
|------------------|--------------------|
| Power supply | 24 V DC |
| Nominal capacity | 48 VA |
| Digital inputs | 3 pieces |
| Analogue output | 1 piece |

DIMENSIONS

8

Dimensions

DIMENSIONS BY WALDNER THE ALL-IN-ONE SOLUTION FOR CONVERTIBLE ROOMS

DIMENSIONS provides you with the technical infrastructure you need to quickly convert rooms. This gives you the flexibility to use rooms one way today and another way tomorrow. It's fast, functional, cost-effective and aesthetic.


You can obtain everything from a single source: innovative technology, excellent product quality, occupational safety, great user-friendliness, ergonomics, sustainability, appealing aesthetics and a feel-good factor.



DIMENSIONS

www.waldner-dimensions.com

Dimensions





9 Accessories

We have designed useful accessories to individually equip certain areas of your working environment as required for our **SCALA** laboratory furniture system..

We will be pleased to show you our accessories that are perfectly adapted to our system.

Make your choice. The complete range of Waldner original accessories can be found in our special catalogue which is available on the Internet at www.waldner.de.

We will also be pleased to send you a printed copy.





10 General

Our innovative developments have made us the market leader in laboratory equipment.

Our products have set the standard for the laboratory workplace worldwide.

We know what our customers expect and we are constantly improving.


We reserve the right to make technical changes in the context of further development. Illustrations, drawings and text content are copyright protected. Re-printing, even of extracts, only with express approval of
WALDNER Laboreinrichtungen SE & Co. KG.



| | |
|--|------------|
| Colours..... | 206 |
| Laboratory planning..... | 208 |
| Awards..... | 210 |
| Installation interfaces mechanical and electrical services..... | 212 |



With respect to design and colour, we placed the emphasis on a balanced appearance with consistency in the application for optimal orientation in the surroundings in which the user spends many hours a day. As a result laboratories can be clearly and timelessly designed for pleasant working.



Pure white

RAL 9010

Similar to NCS S 0602 G91Y

- Storage cupboards
- Internal workspace

Standard for

- Metal parts service module
- Bench frames
- Fume hood fronts



Oak

- Optional as emphasis for storage cupboard fronts



Light grey

NCS S 3005 R80B

Similar to RAL 7040

- Worktops

Option for

- Metal parts service module
- Bench frames



Anthracite metallic effect
Similar to NCS S 5502 R
 Option for
 ■ Fume hood fronts



Glass
NCS S 1010 G10Y
 ■ Worktops back-
 varnished



Dark grey
NCS S 7502 B
 Similar to RAL 7015
 ■ Storage cupboard
 plinth




Stainless steel
 ■ Handles
 ■ Worktops
 ■ Sinks



Pictograms
CMYK 0/16/65/0

■ Emphasising all markings for hazardous goods and special storage units



Our services go way beyond the pure manufacture of laboratory furniture. Due to our many years of experience in the project business, we have acquired fundamental planning competence. We not only equip your laboratory, but on request we will also take over the planning and coordination of all related trades.

The start of planning

The layout planning defines with two-dimensional clarity the intended space utilisation, requirements and existing features, connections, area dimensions, interfaces and other information.

Clear idea using an additional dimension

The laboratory will become clearly conceivable for you through the 3D drawing. We will then refine the details together with you.

In the next stage of the presentation, your laboratory will be almost “accessible” in colour and with clear, differentiated depth in the rendered representation. You will be able to see your laboratory from all angles.

As a logical conclusion to our precise planning and design work, the laboratory will be installed in your building – of course with the usual Waldner quality and on time.




Awards

We have been further developing laboratory furniture for more than 70 years. Over this long period of time, we have had a significant impact on the laboratory workplace with our innovations.

As a result of our attention to detail during development and manufacturing, we have an impressive pool of experience in development, manufacture, planning, installation and service.

Numerous patents, brands, design patents and registered designs clearly demonstrate our innovative power. As the market leader, we will continue to do everything to impress our customers with new and innovative ideas.



Quality right down into the detail is defined not only by our claims about what we do.

We are the first German manufacturer of laboratory furniture to be certified to the quality standard ISO 9001.

ISO 9001 gives you the assurance that you will receive the highest quality products and professional support from the planning phase through to service. Of course, this aspect also covers procurement, development, the technical areas, production and installation.

In-house quality checks and regular training ensure exact observance of the high criteria in ISO 9001.

The products for the **SCALA** laboratory furniture system have been tested by TÜV Product Service GmbH based on all applicable standards and regulations in accordance with the German law on equipment safety and have the GS marking.

These test certificates are only awarded if the manufacturing process is continuously monitored. We have undertaken the obligation to monitor production in several ways: all materials, components and individual parts used in our factory are continuously tested, in some cases also in external test institutes.

Waldner Laboreinrichtungen are environmentally certified. Our active environmental management system meets the EN ISO 14001 guidelines. To us, all aspects matter: From the materials used to the energy efficiency in the production processes, we strive to ensure environmental safety. The renewable resource "wood", for example, is exclusively supplied by regional distributors, our powder coatings do not contain any solvents, the wood left over in the production process covers 85% of our heating requirements, all employees receive continuous training in environmental issues, and the EN ISO 14001 conformity is tested by TÜV Süd at regular intervals.







Installation interfaces mechanical and electrical services

- For water and technical gases, shut-off valves with 1/2" internal threads must be provided on-site according to DIN EN 10226-2
- For pure gases, shut-off valves with 10 mm clamp ring connections must be provided on-site
- For waste water connections, a 56 mm plug sleeve must be provided on-site.

- Electrical supply pipe in acc. with DIN VDE 0100-430
- Type of cable/pipe with on-site fusing upon agreement with Waldner

- Waldner will indicate the transfer points for the on-site trades for each project in the corresponding positional drawings



.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....



A series of 20 horizontal dotted lines spaced evenly down the page, providing a template for writing.



WALDNER Laboreinrichtungen SE & Co. KG
Haidoesch 1 · 88239 Wangen · Germany
Phone +49 7522 986-0 · Fax +49 7522 986-418
labor@waldner.de · www.waldner.de

WALDNER